
- ES - Pulsioxímetro de dedo
- EN - Fingertip Pulse Oximeter
- IT - Pulsossimetro da dito
- FR - Oxyètre de pouls au doigt
- PT - Oxímetro de pulso de dedo
- DE - Fingerpuls-Oximeter
- NL - Vingerpulsometer
- SWE - Fingerpulsoskizmer
- PL - Pulsoksymetr palcowy
- DK - Fingerpulsokszimeter

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## USER MANUAL EN

### FINGERTIP PULSE OXIMETER

This Fingertip Pulse Oximeter is a kind of innovated medical device with non-invasive and continuous features for artery SpO2 and PR detection. Being portable, it is able to measure SpO2 and PR values quickly and precisely.

**NOTICE**  
Notice to the user and/or patient: any serious incident which has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**INDICATIONS**  
The finger pulse oximeter can be used to measure human haemoglobin saturation and pulse rate through the finger, can be used in hospitals, families, schools and medical centres.

**GENERAL DESCRIPTION**  
Hemoglobin Saturation is the percentage between the capacity of Oxyhemoglobin (HbO2) that is compounded with oxygen and that of all combinable hemoglobin (HbO2) in blood. In other words, it is the saturation of Oxyhemoglobin in blood. It is a very important physiological parameter for Respiratory and Circulation Systems. Many respiratory diseases could reduce hemoglobin saturation in human blood. Moreover, factors such as Automatic Organic Regulation Malfunction caused by anesthesia, trauma resulted from major operation and some medical examination can also cause problems in oxygen supply, which might reduce human hemoglobin saturation. As a result, such symptoms as megrium, vomiting and asthma might appear to patients. Hence, it is very important to know hemoglobin saturation of patient timely in clinical medical aspects.  
The Fingertip pulse oximeter features in small volume, low power consumption, convenient operation and portability. It is only necessary that the patient puts one finger into the fingertip photoelectric sensor for diagnosis, and the display screen will directly show measured value of hemoglobin saturation. It has been proved in clinical experiments that it

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possesses rather high precision and repeatability.


**MEASUREMENT PRINCIPLE**  
The principle of the oximeter is as follows: An experience formula of data process is established by exerting Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (R Hb) and oxyhemoglobin (O2 Hb) in glow and near-infrared zones. Operation principle of the instrument is to combine Photoelectric Oxyhemoglobin Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that two lights with different wavelength (660nm glow and 940nm near infrared light) can be focused onto human nail through perspective clamp, finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on two groups of LEDs through process in electronic circuits and microprocessor.

**PRECAUTIONS FOR USE**  
1. Do not use the pulse oximeter together with MRI or CT equipment.  
2. Explosion hazard: Do not use the pulse oximeter in an explosive atmosphere.  
3. The pulse oximeter is intended only as an adjunct in patient assessment. Doctors should make diagnosis in conjunction with clinical manifestation and symptoms.  
4. Check the pulse oximeter sensor application site frequently to make sure that the circulation and skin integrity of patient are under good condition.  
5. Do not stretch the adhesive tape while applying the pulse oximeter sensor. This may cause inaccurate reading or skin blisters.  
6. Please read the manual carefully before your operation.  
7. The pulse oximeter has no SpO2 prompt; it is not for continuous monitoring.  
8. Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity,

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
circulatory status, and correct alignment at least every 2 hours.  
9. Inaccurate measurements may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid.  
10. Significant levels of dysfunctional hemoglobin (such as carboxyl-hemoglobin or methemoglobin) may cause inaccurate reading.  
11. Inaccurate readings may also be caused by indocyanine green or methylene blue may cause inaccurate reading.  
12. SpO2 measurements may be adversely affected in the presence of high ambient light. Please shield the sensor area (with a surgical towel or direct sunlight, for example) if it is necessary.  
13. Unexpected action may cause inaccurate reading.  
14. Medical signal with high frequency or interference caused by defibrillator may lead to inaccurate reading.  
15. Venous pulsations may cause inaccurate reading.  
16. It may cause inaccurate reading when the positions of sensor and blood pressure cuff are on the same arterial catheter or intravascular line.  
17. Hypotension, severe vasoconstriction, severe anemia, or hypothermia may cause inaccurate reading.  
18. It may cause inaccurate reading by giving use of cardiotoxic to patient after his cardiac arrest or when he is in quiver.  
19. Bright nail or painted nail may cause inaccurate SpO2 reading.  
Follow local ordinances and recycling instructions regarding to disposal or recycling of the device and device components.  
Scope of application / intended use  
The Fingertip pulse oximeter is suitable for children under 12 years old to detect blood oxygen saturation and pulse rate. it can be used in hospitals, families, schools and medical centres.

## OPERATION INSTRUCTIONS

1. Hold the lock button on the battery cover, push the battery cover down at the same time and the battery cover is opened.
2. Install two AAA batteries into battery cassette before closing its cover.
3. Plug one finger into the rubber hole of the oximeter (it is best to plug the finger thoroughly), before releasing the clamp with the nail upwards.
4. Press the button on the front panel.
5. Don't tremble with your finger when the oximeter is working. It is recommended to stay quiet while using your oximeter.
6. Press the button on the front panel, if you want to change the display direction.
7. Read relevant data from the display screen.
8. The instrument has a standby mode, where no signal will be shown.
9. Please replace the batteries with new ones when the OLED indicates that the batteries are low in power.
10. When plugging your finger into the oximeter, your nail surface must show upward.

**Brief Description of Front Panel FEATURES:**

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
## OPERATION INSTRUCTIONS

Description of the key function: In standby mode, the instrument can be put into operating mode by pressing the start button. If this key is pressed during operation, the display mode can be changed.  
NOTE: The illustration is for reference purposes only, specific in kind prevail.  
**Parameter setup:**  
Press the Power button again (less than 1 second) to scroll down the mode you want to set up. The signal "" will move each time you press the Power button. Press the Power button more than 3 seconds to start setting up.  
Move the signal "" on the "Alm Setup"; press the Power button more than 3 seconds to enter the second page "Settings"; press Figure 2. The operation way is the same as Figure 1.  
"H" means Higher, "L" means Lower.

"H" means Higher, "L" means Lower.		Settings Sounds Setup *	
Alm Setup *	off	SpO2 Alm Hi	100
Alm off	off	SpO2 Alm Lo	100
Resto OK	OK	PR Alm Hi	100
Restore OK	OK	PR Alm Lo	100
Exit	Exit	PR Alm Lo	100

**Detailed descriptions of product functions:**

1. Display Type: OLED display
2. SpO2 Measurement range: 70%~99%. Accuracy: ±2% on the stage of 80%~99%. ±3% on the stage of 70%~79%. Below 70% no requirement
- Resolution: 11%
3. PR: Measurement range: 30BPM~240BPM Accuracy: ±1BPM or 1% (the largest one)
4. Power: rechargeable lithium battery

- Power consumption: below 0.1mA  
 6. Automatic power-off: the product shuts off by itself when no finger is in the product x8 seconds  
 7. Dimension: 44mm x 28.3mm x 26.5mm  
 8. Operation Environment: Operation Temperature: 5°C ~ 40°C  
 Storage Temperature: -10°C ~ 40°C  
 Ambient Humidity: 15% ~ 80% on operation  
 10% ~ 80% in storage  
 Air Pressure: 86kPa ~ 106kPa  
 9. Declaration: EMC of this product complies with IEC60901-1-2 standard.  
 10. Measurement Performance in Low Perfusion Condition: required test equipment (BIO-TEK INDEX Pulse Oximeter tester) can measure the available pulse wave with the amplitude of 6% of the simulation pulse wave amplitude.  
 11. Interference Resistance Capacity against Ambient Light: Device works normally when BIO-TEK INDEX Pulse Oximeter tester exerts interfering signal test.

#### Classification

- Class of medical device: Class II
- Type of protection against electric shock: Internally powered equipment internal power supply
- Degree of protection against electrical shock: Type BF

#### Maintenance and Preservation

- Replace the batteries timely when the low voltage lamp blinks.
- Clean the surface of fingertip oximeter before it is used to diagnose patients.
- Remove the batteries inside if you will not use the oximeter for a long time.
- It is recommended to preserve the product in -10~40°C (14~104°F) and a humidity of 10%~80%.
- It is recommended to keep the product dry at any time. A wet ambience might affect its lifetime and may even damage the product. Please follow the law of the local government to deal with the disposal of used batteries.

#### Product Accessories

- Lanyard for hanging the device around the neck.
- Two AAA batteries.
- User manual.

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
#### Guidance and manufacturer's declaration-electromagnetic radiation-for other EQUIPEMENTS and SYSTEMS

The Pulse Oximeter is designed to be used in specified electromagnetic environment. Users of the Pulse Oximeter must use it in the following environments:		
Radation Test	Compliance	Electromagnetic environment-guidance
RF interference CISPR 11	Group 1	RF signal of Pulse Oximeter is simply created by its internal function. Therefore, its RF interference is very low and is not likely to cause any interference to nearby electronic equipment.
RF interference CISPR 11	Class B	The Pulse Oximeter applies to all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

#### Possible Problems and Resolutions

Problem	Possible reason	Solution
SPO2 or PR cannot be shown normally	1. Finger is not plugged correctly. 2. Patient's Oxyhemoglobin value is too low to be shown.	1. Beby by plugging the finger correctly. 2. Try those times. If it still doesn't work, contact your dealer.
S602 or PR is shown untimely	1. The finger might not be plugged deep enough. 2. Finger is trembling, or patient is in movement status.	1. Try again inserting your finger correctly. 2. Please remain at rest.
The pulse oximeter cannot be turned on	1. Inadequate power or power cover. 2. The pulse oximeter might be damaged.	1. Please charge the batteries timely. 2. Please contact with local customer service Centre.
Indication lamps are suddenly off	1. The product automatically shut down because it is detected in 8 seconds. 2. Inadequate power.	1. Normal. 2. Please charge the batteries timely.

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#### Symbols and Definitions

Symbol	Definition	Serial number	
	Protection and isolation of medical device		Serial number
	Generate collection		Date of manufacture
	See the Instruction		Manufacturer
	Warning		Warning of no SpO2
	IP degree		Temperature range
	Humidity range		Avoid sunlight
	Keep dry		This way up

Reserves the right to technical change appearance, our products are subject to change without prior notice, please forgive me!

#### HOW TO REPLACE BATTERIES:

- To remove the back cover compartment, push the white button and follow the direction of the printed arrows.
- Place two AAA batteries in the battery compartment. Match the plus (+) and minus (-) signs in the compartment. If the polarities are not matched, damage may be caused to the oximeter.
- Slide the battery door cover horizontally along the arrow shown as the picture.

#### Note:

Please remove the batteries if the pulse oximeter will not be used for long periods of time. Please replace the battery when the power indicator starting flickering.

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
#### Installation of hanging rope

- Insert the thin end of the rope through the hole.
- Insert the coarser end of the rope through the thin end part and tighten it.

#### Warnings!

- Keep the oximeter away from young children. Small items such as the battery cover, battery, and lanyard are choking hazards.
- Be wary that the lanyard does not get twisted with the device's electrical wire.
- Please notice that the lanyard tied to the oximeter may cause strangulation due to excessive length.

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