MANUAL DE INSTRUCCIONES

PM10

Electrocargiógrafo portátil

Last revision date: May 25th 2022

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CMS2.782.350(CE,USA)ESS/1.6

Aviso

Notice to the user and/or patient: any serious incident that has occurred in relation to the product must be informed to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Prologue

Read the user manual carefully before using the product and strictly follow its operating rules. The user manual indicates operations to which users should pay close attention, which may lead to abnormalities, or may endanger the device or human body during use. Our company shall not be liable for the safety, reliability and performance of any abnormality or damage to the device and human body caused by not following this user manual to use, maintain and store, nor provide free service for any above situation.

We apologise for the content of the manual is subject to change according to product updates without prior notice.

The product is reusable as a medical instrument.

Warning:

Reliability depends on whether or not users follow the operation and maintenance of the user manual. Our company's website: <http://www.contecmed.com> is the only way to download the APP software and PC software as well as firmware updates, if the user downloads the software and firmware updates from other unauthorised channels, it will cause the risks associated with cybersecurity, our company will not be responsible for the consequences it may cause. All revisions and future updates of the device must be carried out by personnel

The schematic diagram and list of components can only be provided to the service station or service personnel designated by our company. The user must be aware of the battery life cycle, and dispose of the battery in accordance with local laws when the battery life cycle is exhausted. Do not expose the device to sunlight. Do not expose the device to a magnetic resonance (MR) environment. The device may present a risk of projectile injury due to the presence of ferromagnetic materials that may be attracted to the core of the MR magnet. Thermal injury and burns may occur because the metal components of the device may become hot during MR scanning. The device may generate artefacts in the MR image. The device may not function properly due to the strong magnetic and radiofrequency fields generated by the MRI scanner. This user manual contains proprietary information, which is protected by copyright. All rights are reserved. No part of the manual may be reproduced, adapted or translated without prior written permission. Our company assumes the following responsibilities To provide users with qualified products in accordance with the company's standard. To provide installation, debugging and training services in accordance with the contract. Provide one year warranty and product maintenance after the warranty period in accordance with the contract. Respond to user requests in a timely manner.

Chapter 1. Notice

1.1 Indications for use

The device is a handheld personal electrocardiograph, which can measure the electrical activities of the heart easily and conveniently. It is immediately available at any time to manually record transient cardiac events, it is suitable for home healthcare use, it can detect, display and store the ECG signal, and if possible, provide the average heart rate message after ECG measurement. The users can use it themselves to check their cardiac status. It is suitable for adult users, who suffer from cardiovascular diseases, or adult persons who are concerned about their heart working conditions during their daily life. This device is not intended to be used as a conventional diagnostic tool, but as a health tool that can provide the physician with the recorded data for reference.

1.2 Generic notice

Do not use the device in locations subject to high temperature or humidity. Use at a temperature not exceeding 5~40 °C and humidity not exceeding 25%~80% RH. Do not wash the device with water. The presetting time is within 30mins, in condition of: warm from minimum storage temperature until ready for use at a room temperature of 20 °C cool from maximum storage temperature until ready for use at a temperature of 20 °C. Do not use or store the device in the following environmental conditions near fires or open flames Places exposed to strong vibrations Places exposed to strong electromagnetic fields Do not disinfect the device in an autoclave or gas steriliser. In case of allergy or skin damage, do not use this device. The service life of the device is 3 years. Do not throw away the device and accessories

when they are no longer functional. If the device needs to be disposed of, it must comply with the requirements of local laws and regulations. The responsible organisation should contact local authorities to determine the appropriate method of disposal of potentially biohazardous parts and accessories. Please do not use multiple wireless devices connected to the product at the same time. This device has no contraindications. The parameters displayed by the ECG are to be interpreted by a professional professional. Please do not use the device for children weighing less than 10 kg. The device contains some tiny components that can be put in the mouth by children, which will cause choking or damage to organs (including the oesophagus and stomach), so please keep it out of the reach of children.

1.3 Measuring notice

If your skin is dry, wipe it with rubbing alcohol or electrical balm to strengthen the electrical capacity. It is best to sit comfortably, stand up and start measuring when your heart rate stabilises. When measuring, the finger and chest electrodes should touch your skin exactly, round and well.

1.4 Security Notice

Sampling should not be performed during battery charging (Sampling means acquisition of patient ECG data at designated anatomical sites). When the battery is charging, the device shall not record ECG. The following symbol shall be presented on the user interface to remind the charging status, the device cannot be operated in the battery charging status. Place the device in a shaded and cool environment when not in use when not in use for a long period of time, and electrify the device every three months. Do not use the device in an environment where flammable objects, such as an anaesthetic, are present. Other equipment connected to it must comply with national safety standards. That the conductive parts of the ELECTRODES and associated connectors for TYPE BF, including the NEUTRAL ELECTRODE, must not come into contact with other conductive parts, including earth.

1.5 EMC Notice

Please be aware of the effect of EMC when using the device, as it may be influenced by portable or mobile RF devices with high electromagnetic compatibility. This equipment must be installed and commissioned in accordance with the information provided in the accompanying documents. Wireless communications equipment may affect the equipment and must be kept at least a distance d from the equipment. The distance d is calculated by the manufacturer from the 800 MHz to 2,5 GHz column of Table 5 or Table 6 of IEC 60601-1-2:2007.

1.6 RF Instruction

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions:
A. This device may not cause harmful interference;
B. This device must accept any interference received, including interference that may use undesired operation.
Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Some electronic devices are susceptible to electromagnetic interference from this equipment if they are inadequately shielded. Please operate this equipment at a minimum distance of 20 cm or as far away as possible from the TV, radio and other automated office equipment to avoid interference. This device is a radio transmitter and receiver. It is designed and manufactured not to exceed the limits for exposure to radio frequency (RF) energy set by the Federal Communications Commission (FCC) of the United States Government. These limits are part of comprehensive guidelines and establish permitted levels of RF energy for the general population. The guidelines are based on standards developed by independent scientific organisations through regular and thorough evaluation of scientific studies. The standards include a substantial safety margin designed to ensure the safety of all people, regardless of age or health.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help, for assistance.

A minimum separation distance of at least 0.2 m must be ensured between this equipment and all persons to satisfy RF exposure compliance.

1.7 Quality of service and safety

The device ensures timely wireless data and information transfer, reliable, accurate and secure thanks to the following design. When you wish to establish a wireless connection to the portable ECG monitor, you must enter the correct communication instructions. Therefore, unauthorised access to ECG data is prevented.

1.8 Manufacturer's information

Manufacturer's name: Contec Medical Systems Co., Ltd.
Manufacturer's address: No.112 Qinhuang West Street, Economic & Technical Development Zone, Qinhuangdao, Hebei Province, PEOPLE'S REPUBLIC OF CHINA
Website: <http://www.contecmed.com>
Telephone: +86-335-8015430
Fax: +86-335-8015588
Technical support: +86-335-8015431
Email: cms@contecmed.com.cn
Website: <http://www.contecmed.com>

1.9 Intended operator

Non-professional person, with the following requirements:

Be able to read and understand the contents of the user manual;
Be able to distinguish the following anatomical sites: thorax, left/right palm, upper extremity and lower extremity.

Chapter 2. Introduction

The portable ECG monitor is designed for family and individual users. It is a good help for family members to prevent cardiovascular diseases. The device can record and display the user's ECG waveform and heart rate anytime and anywhere with easy operation.

2.1 Characteristics

Attractive shape, practical operation, comfortable to wear.
Real-time heart rate recording anytime, anywhere.
Built-in high capacity rechargeable lithium battery, continuously display 200 ECG waveforms after charging once.
Measurement of QRS and VEB intervals

2.2 Implementation

Occasion: family, medical clinic and hospital. The device cannot be used as a general electrocardiogram for clinical examination.
Subject: Persons under high pressure and long term workload, patients with heart disease, medium and patients with cardiac diseases, middle-aged and elderly people, people with elderly, persons with health problems.
Purpose: The device is only used for heart rate recording. The method of operation is simple and requires less requirements for operating personnel.

Chapter 3. Primary technical orders

3.1 Working environment

Temperature: +5 ~ +40 °C
Relative humidity: 25% ~80%.
Atmospheric pressure: 70 kPa ~106 kPa
Power supply: built-in rechargeable lithium battery, voltage: 3.7V
Transport and storage environment
Temperature: -40 ~+55 °C
Relative humidity: <=95%.
Atmospheric pressure: 50 kPa ~106 kPa


3.2 Basic parameters

Lead: Standard 3-wire, II III
Calibration voltage: 1 mV±5%.
Standard sensitivity: 10 mm/mV±5%.
Amplitude frequency characteristic: standard: 10 Hz; 1 Hz 20 Hz; (+0.4 dB, -3 dB)
Noise level: ≤30 μV
Input impedance: ≥50 MΩ
CMRR: ≥60 dB
Scan rate: 25 mm/s±5%.
Sampling rate: 250 points/s
HR measurement range: 30 bpm 300 bpm, error: ±1 bpm or 1%.
Battery voltage: DC 3.7 V
Type of protection against electric shock: Internal power supply device
Degree of protection against electric shock: Applied part type BF
Degree of watertightness: IP22
Display: 1.77" TFT-LCD colour display
Size: 100 mm (L)*45 mm (W)*15 mm (H)
Weight: approx. 60 g
Bluetooth Specifications:


Bluetooth Protocol: Bluetooth V4.0, Classic and Low Energy
Operating frequency: 2.4 GHz ISM band
Modulation: GFSK (Gaussian Frequency Shift Keying)
Transmit Power: <10 dBm, Class 1
Sensitivity: -92 dBm at BER<.001 with DHS
Transfer rate: 1.5 Mbps data throughput
Security features: Authentication and encryption
Support Services: Bluetooth SPP/IAP/HID/OBEX
Effective RF radiated power <10 dB
Operating range: 2400 MHz - 2485 MHz
Receive section bandwidth: 85 M
FCC ID: 2A80GPM10

Chapter 4. Instructions for use

4.1 The sketch map and the name of components

4.2 How to use it

Caution: Ensure that the electrode is in full contact with the skin.

4.3 Menu operations

Switching on:
Press and hold the on/off button for 3 seconds, you will hear a beep and see the display illuminated. The device will keep the level off when not measuring.
Starting the measurement:
After start-up, the device will enter the pre-sampling interface. Please use the correct measurement method according to the guide, the ECG waveform and heart rate will be displayed on the screen, as shown in Fig.

4.1. The method of heart rate calculation: the number of heartbeats without interference in the ECG fragment is set as N, then the formula of heart rate calculation is as follows: HR=60000/(Sum of R-R intervals over the number(N) of heartbeats/N)


Figure 4.1 Pre-sample interface

When the waveform stabilises, the device will start formal sampling automatically, the sampling time countdown in the lower right corner starts until the sampling is finished at once and the colour of the sampling time changes to red. See figure 4.2:


Figure 4.3 Case review interface

The device will enter the case review interface after completing the sampling. The case review interface displays the sampling start time and the heart rate, see figure 4.3.


Figure 4.2 Example of a formal interface

When the device enters the case review interface, it shall display the last sampled case. Click the button to review the information of other cases. The device can store a maximum of 100 cases. If the limit is reached, the newly stored case will cover the original case, the one that was stored first, piece by piece. The device shall automatically switch to the sampling interface to continue if the user holds the electrode at both ends again when the device is in the case review interface.

Notes on battery operation

The device can work continuously for more than 2 hours when the battery is fully charged. battery life up to 200 times.

Two charging methods:

Connect the device with a computer using a Micro USB cable, charging is completed after about 2 hours.
Use a Micro USB to connect the device with a power adapter (output current >500 mA, 5V), charging is completed after about 2 hours.
When the battery is charging, the device will not record ECG. The following symbol shall be presented on the user interface to remind the charging status, the device cannot be operated in the battery charging status.

Figure 4.4 Battery display

No.	Indicator	Maximum capacity Description
a		Capacity: 1/2
b		Capacity: 1/4
c		Capacity: 1/4
d		If the battery is low, it is recommended that it is recharged and the appliance will switch off automatically.
e		

Auto Power Off

The device shall automatically shut down if no operation is performed within 1 minute.

4.4 PC Sync Software Operation and Communication

The intended use of the PC software:
The PC management software is intended to be used as support software for portable ECG device. Functions include configuration of device parameters (language settings and acquisition types, etc.), downloading of ECG data from the portable ECG device, data management.
PC software:
Users can operate on the PC synchronous software as required, which includes sampling mode and time settings, case upload, case review, measurement, etc.
Software installation:
Run the installation software, and you can see a window as below, click "OK" button.


Figure 4.5

Click the "Next" button, and the dialog box shown in Fig.4.6 will appear. If you click "Browse...", you can set the installation path; otherwise, the default installation path is "CPN-PM10". Click the "Next" button again, and the dialog box shown in Fig. 4.7 will appear. Click "Browse...", you can set the appeal to the folder in the start menu, the default folder will be "PM10".


Figure 4.6

Figure 4.7

Click the "install" button, and the software will be installed at the designated location.

Chapter 4. Instructions for use

When the installation is complete, the dialog box shown in Fig. 4.9 will appear. Fig. 4.9.


Figure 4.8

Figure 4.9


Figure 4.10


Figure 4.11

Data communication: Turn on the device, insert it into the USB port, click "New Search", then the software starts searching for the device, as shown in Figure 4.13.


Figure 4.12

Figure 4.13

After the search, the device information will be displayed in the form of "List", including: case ID, time duration, check time, heart rate, which is shown in Figure 4.14. Click the "?" button, you can get operation help.


Figure 4.14

Download cases: double click on a selected case to download it, or select multiple cases, then click on "Download selected" to download these cases, or click on "Download all" to download all cases.


Figure 4.15

Delete a case: select a case or several cases, then click "Delete selected" to delete the selected case, or click "Delete all" to delete all cases. To avoid errors, before deleting, the system will ask the user, the system will delete the records after selecting "Yes".



Figure 4.16

Firmware update: click on "Update", select a firmware (naming rules: file name_version.bin, for example: xxx_USB_Vxxx.bin). If the firmware is to be upgraded from the low version to the high version, click "OK". Otherwise, first select "Force upgrade" and then click "OK". During the update, maintain the connection status between the computer and the device.

Setting device parameters: Languages and sampling time can be configured by the user. The interface for successful setting is shown as in the following picture:

The PC management software can run on the WIN8/WIN 10 operating system, the PC hardware requirements for running the WIN10/WIN 8 operating system are listed below.
Hardware specifications for WIN 10 operating system
 Processor: 1.0 GHz or faster
 Display: 800x600+ resolution
 Firmware: UEFI 2.3.1, secure boot supported
 Bootable memory: 2 GB(64-bit), 1 GB(32-bit)
 Hard disk space: not less than 16 GB (32-bit), not less than 20 GB (64-bit)
 Graphics card: Direct X9 tablet support
Hardware specification for WIN 8 operating system
 Windows 8 PC (32 bit ,64 bit) 32 bit (x86)

System hardware requirements Windows 8 32-bit

	Minimum configuration	Suggested configuration	Optimal configuration
Central processor	1 GHz (PAE, NX and SSE2 compatible)	2 GHz (PAE, NX and SSE2 compatible)	2 GHz multi-core processor
Memory	1 GB	2 GB	2 GB DDR3
Graphics Card	Microsoft Direct X9 graphics device with WDDM driver software	Direct X10 graphics device with WDDM driver software, with 128 MB of memory	Direct X11 graphics device with WDDM driver software, with 256 MB of memory
Residual hard disk space	16GB	More than 30 GB	64 GB SSD hard disk

System hardware requirements Windows 8 32-bit

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Tableta Windows 8

Architecture	ARM, X86AND X86-64
Central processor	NVIDIA Tegra/Qualcomm/Texas Instrucción TI ARM, AMD x86/x64 and Intel x86/x64
Memory	More 4G suggested
Graphics card	Adopt WDDM 1.2 or DirectX10 graphics device with higher driver software version
Hard disk/NAND Flash	10GB memory capacity or more after system installation
Push button	"n/off", "rotation lock", "Windows key", "volume up", "volume down"
Touch screen	At least 5-point capacitor-controlled touch screen, resolution greater than 1366 x 768, display must be 10 inches or larger
Display resolution	Resolution 1280 x 720
USB 2.0	At least one power socket
Internet connection	Wi-Fi and bluetooth 4.0
System boot firmware	UEFI
Other	Speaker, microphone, gyroscope, gravity sensor

4.5 Operation and communication of the mobile application

Intended use of the mobile application: The mobile application is intended to download and display PM 10 case data, will not monitor PM 10 and will not alter any PM 10 settings.
 Software installation
 Run the installation software, and you will see a window like the one below, click on the "install" button.
 The mobile application can be installed on a mobile phone with Android system, currently,

the mobile application is not suitable for mobile phones installed with Apple system.
 The mobile application has the following functions:
 - Connect to the PM10 via Bluetooth.
 - Download case data (date, time, duration of measurement and average heart rate)
 - View, store and review case data.
 - Data communication
 Start the software, turn on the device, then the software starts searching for the device, click on "PM10".
 Device information will be displayed, including: measurement duration, check time, heart rate.

Chapter 5. Troubleshooting and problem solving

If the device has an account of problems, please look at the following sheet for solutions first, if you are not included in the following topics and cannot solve any, please contact customer service.

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Chapter 6. Maintenance, transport and storage

6.1 Cleaning and disinfection

Switch off the device before cleaning. To disinfect the device, you can use medical alcohol and then air dry it. Or simply wipe it with a clean, dry cloth to clean it. Do not allow any liquid to enter the device.

6.2 Maintenance

- The maintenance personnel appointed by our company must not open the housing of the device in order to avoid damage to the internal components.
 - Any maintenance and upgrading of the equipment must be carried out by trained professionals authorised by the company.
 - Prevent any liquid from leaking into the device, as it will affect the safety and performance of the device.
 - The device should avoid violent shocks or impacts.
 - Do not place objects on the device. This may damage the touch screen.
 - If you do not use the device for a long time, charge the battery to the maximum every 3 months, otherwise it will cause permanent damage to the battery.
 - The device should not be held during use.

6.3 Transport and storage

The transport of the device adopts the general means of transport or follows the requirements of the contract. Avoid violent shocks, vibrations, rain and snow splashes during the transportation process.
 Store the packaged device in an environment with temperature -40 °C ~+55 °C, relative humidity not exceeding 95%, atmospheric pressure 500 hPa ~1060 hPa, no corrosion gas and in a well-ventilated room.

Chapter 7. Explanation of symbols

Icon	Description	Icon	Description
	Please refer to the instruction manual/brochure		TYPE BF
	Pulse rate (bpm)		Bluetooth
	Power button / function button		International protection
	USB		WEEE (2012/19/EU)
	Moisture limitation		Atmospheric pressure limitation
			Temperature limitation
	This article complies with Directive 93/42/EEC of 14 June 1993 concerning medical devices; including, as of 21 March 2010, amendments to Council Directive 2007/47/EEC.		Manufacturer
	Serial number		This side ABOVE
	Serial number		Keep dry
	Date of manufacture		Recyclable
	Fragile, handle with care		Authorised representative in the European Community
	Do not use this equipment in the MRI room.		

Capítulo 8. Packing list

No	Description	Quantity
1	Host	1
2	Cable USB	1
3	User manual	1

Chapter 9 Electromagnetic compatibility and interference

Guidance and Manufacturer's Declaration - Electromagnetic Emissionsfor all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emission

The PM10 is intended for use in the electromagnetic environment specified below. The customer of the user of the PM10 must ensure that it is used in such an environment.

Emissions testing	Compliance	Electromagnetic environment - guidance
Emission of RF CISPR 11	Group 1	The PM10 uses RF energy only for its internal operation. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Emission of RF CISPR 11	Class B	PM10 is suitable for use in all establishments, including domestic establishments and those directly connected to a low voltage power supply network supplying buildings used for domestic purposes.

Guidance and manufacturer's declaration: electromagnetic immunity

The PM10 is intended for use in the electromagnetic environment specified below. The customer or the user of the PM10 must ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Level of compliance	Electromagnetic environment - orientation
Electrostatic discharge (ESD)	±8 kV contact	±8 kV contact	Floors should be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%.
IEC 61000-4-4	±15 kV air	±15 kV air	
Power frequency (50/60 HZ)	30 A/m	30 A/m	The quality of the electricity power quality should be that of a typical commercial or hospital environment.
magnetic field IEC 61000 4-8 NOTE			

Guidance and manufacturer's declaration - electromagnetic electromagnetic immunity - for EQUIPMENT and SYSTEMS other than life-supporting

Guidance and manufacturer's declaration - electromagnetic immunity

The PM10 is intended for use in the electromagnetic environment specified below. The customer or the user of the PM10 must ensure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Level of compliance	Guide to the electromagnetic environment
Portable and mobile RF communications equipment should be used no closer to any part of PM10, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.			
RF radiated IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	Recommended separation distance 80 MHz to 800 MHz 800 MHz to 2.7 GHz
			$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
			$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$
			$d = \left[\frac{7}{E_1} \right] \sqrt{P}$

where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2 These guidelines may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
 Field strengths from fixed transmitters such as base stations of radiotelephones (cellular/wireless) and land mobile radios, amateur radio, AM and FM radio broadcasts, and television broadcasts cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the PM10 is used exceeds the RF compliance level applicable above, the PM10 should be observed to verify normal operation. If abnormal operation is observed, it may be necessary to take additional measures, such as reorienting or relocating the PM10.

Recommended separation distances between communications equipment communications equipment RF communications equipment and the EQUIPMENT or SYSTEM -- for non-life-supporting EQUIPMENT or SYSTEMS Recommended separation distances between portable and mobile RF communications equipment and the PM10

The PM10 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the PM10 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the PM10 as recommended below, depending on the maximum output power of the communications equipment.

Separation distance according to the frequency of the transmitter (m)

Rated maximum transmitter power output (W)	80 MHz to 800 MHz	800 MHz to 2.7 GHz
0.01		0.23
0.1	0.12	0.74
1	0.37	2.38
10	1.17	7.38
100	3.69	23.33
	11.67	

$$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$$

For transmitters with a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

CONTEC MEDICAL SYSTEM CO., LTD
 No. 112 Qinhuang West Street, Economic & Technical Development Zone, Qinhuangdao, Hebei Province, PEOPLE'S REPUBLIC OF CHINA

EC REP Legal representative inside EU / Representante legal en la UE: Shanghai International Holding Corp. GmbH (Europe) Eiffestrasse 80, 20537 Hamburg, Germany

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Imported by/Importado por: Grupo R. Queralto S.A.