
MANUAL DE INSTRUCCIONES

PULSIOXÍMETRO

CMS60F

Manual de instrucciones

ES

Instruction manual

EN

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
EN

NOTICE

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

WARNINGS

- The device may cause discomfort or pain if the device is used repeatedly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours
- For certain patients, there should be more care in the placing process. The device can not be clipped on the edema and tender tissue.
- The light (the infrared is invisible) emitted from the device is harmful to the eyes, so the user and the maintenance personnel should not stare at the light.
- Testee should not use nail polish or other makeup.
- Testee's fingernail should not be too long.
- Please refer to the correlative literature about the clinical restrictions and caution.
- This device is not intended for treatment.

1. SAFETY**1.1. Instructions for safe operation**

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance. It is recommended that the device should be inspected weekly at least. When there is obvious damage, stop using the device.
- Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
- The Oximeter cannot be used together with devices not specified in User's Manual. Only the accessory that is appointed or recommendatory by manufacture can be used with this device.
- This product is calibrated before leaving factory.

1.2. Warning

- Explosive hazard—DO NOT use the Oximeter in the environment with tinder such as anesthetic .
- Ensure that the environment in which the device is operated is not subject to any sources of strong electromagnetic interference, such as radio transmitters, mobile telephones, etc. Keep them far away high level electromagnetic radiation emitted from such devices may greatly affect the monitor performance.
- DO NOT use the Oximeter while the patient is being scanned by MRI or CT. N" The person who is allergic to rubber can not use this device.
- The disposal of scrap instrument and its accessories and packing (including battery, plastic bags, foams and paper boxes) should follow the local laws and regulations.
- Please check the packing before use to make sure the device and accessories are totally in accordance with the packing list, or else the device may have the possibility of working abnormally.
- Please choose the accessories and probe which are approved or manufactured by the manufacturer, or else it may damage the device.
- The device can only be matched with the compatible probe.
- Please don't measure this device with functional tester for the device's related information.
- The equipment uses lithium polymer battery, attention should be paid to that improper operation (put in a high temperature environment, short circuit etc) will lead to the risk of fire or explosion.
- Before using the equipment, ensure that the battery is full, so as to avoid influencing the normal use during the measurement for low energy, or bringing you unnecessary losses.
- No modification of this equipment is allowed.

1.3. Attention

- Keep the Oximeter away from dust, vibration, conosive substances, tinder, high temperature and moisture.
- If the Oximeter gets wet, please stop operating it.

- When it is carried from cold environment to warm or humid environment, please do not use it immediately.
- DO NOT operate keys on front panel with sharp materials.
- High temperature or high pressure steam disinfection of the Oximeter is not permitted. Refer to User Manual in the relative chapter (7.1) for instructions of cleaning and disinfection.
- DO NOT have the Oximeter immersed in liquid. When it needs cleaning, please wipe its surface with medical alcohol by soft material. Do not spray any liquid on the device directly.
- When cleaning the device with water, the temperature should be lower than 60°C.
- The fingers which are too thin or too cold may affect the measure accuracy, please clip the thicker finger such as thumb and middle finger deeply enough into the probe.
- The pulse Oximeter can be used to adult, child and infant. Whether the device is used to adult or children, it depends on the probe selected.
- The update period of data is less than 5 seconds, which is changeable according to different individual pulse rate.
- Please read the measure value when the waveform on screen is equably and steady-going. This measure value is optimal value, and the waveform at the moment is the standard one.
- If some abnormal conditions appear on the screen during test process, pull out the finger and reinsert to restore normal use.
- The device has normal life for three years since the first electrified use.
- The device has alarm function, users can check on this function according to chapter 6.1 as reference.
- The device has the function of limit alarm. When the measure data is beyond the highest or lowest limit, the device would start to alarm automatically if the alarm function is on.
- The device has alarm function. This function can either be paused, or closed for good. This function could be turned on through menu operation if you need. Please check the chapter 6.1 as reference.
- The sign "BATTERY" in the PCB board means where the lithium battery should be welded, and the symbol "+" indicates the battery anode.
- The device may not work for all patients. If you are unable to achieve stable readings, discontinue use.

2. OVERVIEW

The pulse oxygen saturation is the percentage of HbO₂ in the total Hb in the blood, so-called the O₂ concentration in the blood. It is an important bio-parameter for the respiration. A number of diseases relating to respiratory system may cause the decrease of SpO₂ in the blood, furthermore, some other causes such as the malfunction of human body's self-adjustment, damages during surgery, and the injuries caused by some medical checkup would also lead to the difficulty of oxygen supply in human body, and the corresponding symptoms would appear as a consequence, such as vertigo, impotence, vomit etc. Serious symptoms might bring danger to human's life. Therefore, prompt information of patients' SpO₂ is of great help for the doctor to discover the potential danger, and is of great importance in the clinical medical field.

The Pulse Oximeter is small in volume, convenient in operation and portable. It is only necessary for patient to put one finger into probe for diagnosis, and display screen will directly show the measure value of pulse oxygen saturation with the high veracity and repetition.

2.1. Features

- Operation is simple and convenient
- Small in volume, light in weight and convenient to carry.
- Low power consumption.

2.2. Major Application and Scope of Application

The Pulse Oximeter can be used in measuring the pulse oxygen saturation and pulse rate through finger. The product is suitable for being used in hospital, oxygen bar, community healthcare, physical care in sports (It can be used before or after doing sports and it is not recommended to use the device during the process of having sport) and etc.

The problem of overrating would emerge when the patient is suffering from toxicosis which caused by carbon monoxide, the device is not recommended to be used under this circumstance.

2.3. Environment Requirements

- Storage environment
 - a) Temperature: $-40^{\circ}\text{C}\sim 60^{\circ}\text{C}$
 - b) Relative humidity: $\leq 95\%$
 - c) Atmospheric pressure
- Operating environment
 - a) Temperature: $10^{\circ}\text{C}\sim 40^{\circ}\text{C}$
 - b) Relative humidity: $\leq 75\%$
 - c) Atmospheric pressure: $700\text{hPa}\sim 1060\text{hPa}$

3. PRINCIPLE

Principle of the Oximeter is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (Hb) and Oxyhemoglobin (HbO₂) in glow & near-infrared zones. Operation principle of the device is: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning & Recording Technology, so that two beams of different wavelength of lights can be focused onto human nail tip through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.


Figure 1

4. TECHNICAL SPECIFICATION

4.1. Main performance

- SpO₂ value display
- Pulse rate value display, bar graph display PI value display
- Pulse waveform display
- Low-power indication: low-power indication symbol appears before working abnormality which is due to low-power.
- Review function
- Screen brightness can be adjusted Pulse sound indication
- With alarm function
- With data record function of multiuser, continuous record 24 hours data for each user, and the record data can be uploaded to computer.
- It can be connected with adult, child, infant Oximeter probe With clock function
- With locking touch key function

4.2. Main parameters

- **Measurement of SpO₂**

Measuring range: 0 100% Accuracy:

When the SpO₂ measuring range is 70% 100%, the permission of absolute error is +2%;

Below 70% unspecified.

- **Measurement of pulse rate**

Measuring range: 30bpm 250bpm

Accuracy: +2 bpm or +2% (select larger)

- **Measurement of PI**

Range: 0 20%

- **Resolution**

SpO₂: 1%, Pulse rate: 1bpm

- **Measurement performance in weak filling condition:**

SpO₂ and pulse rate can be shown correctly when pulse-filling ratio is 0.4%. SpO₂ error is +4%; pulse rate error is +2 bpm or +2% (select larger).

- **Resistance to surrounding light:**

The deviation between the value measured in the condition of man-made light or indoor natural light and that of darkroom is less than +1%.

- **Power supply requirement: 3.6 V DC 4.2V DC.**

- **Optical sensor**

Red light (wavelength is 660nm, 6.65mW)

Infrared (wavelength is 880nm, 6.75mW)

- **Adjustable alarm range:**

SpO₂: 0~100%

Pulse Rate: 0bpm~254bpm

5. INSTALLATION

5.1. View of the front panel

Figure 2. Front view

Probe jack : It is used to connect Oximeter probe to measure the oxygen saturation and pulse rate.

USB port :It is used to connect a personal computer to export the trend data by data line.

5.2. Left view and Right view


Figure 3. Left view and right view

5.3 Probe connection


Figure 4. Installation

Inserting the Oximeter probe into the USB port of the pulse Oximeter as Figure 4. (The probe is limited to the one that is produced by our company; and can't be replaced with the similar one by other manufacturers).

5.4. Accessories

- A. A User Manual
- B. A power adapter
- C. A data line
- D. A disk (PC software)
- E. An Oximeter probe

6. OPERATING GUIDE

6.1. Application Method

6.1.1. Connect the probe well according to the method introduced in chapter 5.3, then put the finger into the probe. Refer to Figure 5.

- A. Long press the "power button", until the device is turned on.
- B. Do not shake the finger and keep the patient in a stable state during the process.
- C. The data can be read directly from the screen in the measure interface.

 **Fingernails and the luminescent tube should be in the same side.**

 **If the alarm function is on, the device will provide medium-priority alarm signal when probe or finger is out. Intermittent alarm will occur and the user interface presents "FINGER OUT". Medium priority indicating that prompt operator response is required.**

6.1.2. Pause alarm


Figure 5.

- A. Alarm including the alarm of measure data going beyond the limits, the alarm of low-voltage, the alarm of finger out.
- B. When alarm is on, short press the "power button" can pause the alarm, it can renew alarm after period of time, alarm pause time can be set by menu.
- C. Only alarm sound can be closed, the prompt information displayed can't be closed.

6.1.3. Menu operations

In the measuring interface, press the "menu button" can enter the main menu interface of Figure 6.


Figure 6. Main menu interface

A. Alarm setting

In the main menu interface, move the choice bar to "Alarm" item, then press "right button" to enter the alarm setting menu of Figure 7:

Alarm Menu		
High Limit		
SpO2	100	
PR	120	
Low Limit		
SpO2	83	
PR		
Alarm	off	
Pulse Sound	off	
OK	13:57	Back

Figure 7. Alarm setting menu

a. The high/low limit of alarm setting:

In alarm setting menu, you can set the high/low limit of alarm. Move the choice bar to the menu item that you want to set, and press "right button" to begin to set value. Press "menu button" for affirming.

⚠ If the alarm function is on, the device will provide medium-priority alarm signal when the data of SpO2 or pulse rate is beyond the limit. Intermittent alarm will occur and the value shows in yellow. Medium priority indicating that prompt operator response is required.

b. The alarm setting

In alarm setting menu, move the choice bar to the "Alarm" item, then set the alarm state by pressing "right button". Choose "on" to turn on the alarm, and choose "off" to turn off the alarm.

c. Pulse sound indication setting

In alarm setting menu, move the choice bar to the "Pulse Sound" item, then press "right button" to set state, choose "on" to turn on pulse sound, and choose "off" to turn off pulse sound.

B. Clock setting

In the main menu interface, move the choice bar to "Clock" item, then press the "right button" to enter the clock setting menu of Figure 8:

Clock Menu		
Year	2009	
Month	07	
Day	03	
Hour	14	
Minute	31	
Week	Fri.	
OK	13:57	Back

Figure 8. Clock setting menu

Move the choice bar to the menu item that you want to set, and press "right button" to begin to set. Press "menu button" for affirming. After resetting time, press "return button", and the prompt information whether to affirm setting will appear. The interface will return to main menu after choosing.

C. System setting

In the main menu interface, move the choice bar to the "System" item, then press the "right button" to enter the System setting menu as Figure 9:


Figure 9. System setting menu

- In "System Menu" interface, move the choice bar to the "Version" item, then press the "right button" to see the edition information of hardware and software.
- In "System Menu" interface, move the choice bar to the "Equipment ID" item, then press the "right button" to see the ID information of device. The ID of the device can be set by the PC software. Please refer to <SpO2Assistant user manual> for detail.
- Move the choice bar to the "Memory Manager" item, then press the "right button" to enter "Select ID" interface as Figure 10. Select the patient ID that you want to delete, then press menu button to delete the patient ID and storage data.


Figure 10. Select ID interface

- In "System Menu" interface, move the choice bar to the "Factory Default" item and press the "right button" to choose "Yes", then press menu button, and the device will resume leave-factory state.
- In "System Menu" interface, move the choice bar to the "Volume" item, then you can set the volume of voice.
- In "System Menu" interface, move the choice bar to the "Pause Time" item, then you can set the pause time of alarm. The default of pause time is 60 seconds, and press the "right button" to adjust the pause time of alarm.

g. In "System Menu" interface, move the choice bar to the "Key Sound" item, and press the "right button" to turn on/off the key sound.

D. Brightness adjustment

In the main menu interface, move the choice bar to the "Brightness" item, then press the "right button" to adjust the brightness of screen

E. Data storage setting

a. In the main menu interface, press the "up button" or "down button" to select "Record", then press "right button", here the "off" will turn to black display. Press "right button" again, and the black "off" will turn to black "on". Press menu button to enter the "Input ID" interface as Figure 11.

Figure 11. Input ID interface

b. In "Input ID" interface, press menu button to input new ID, here a glimmer cursor will appear in the ID input box and the red choice box will appear in the alphabet as Figure 12. Press direction key to choose the letter which is wanted to input and press menu button to input. You can input 4 letters most. After inputting the letters, move the choice box to the enter sign and press menu button to confirm input, here there is not a glimmer cursor in the ID input box. Press menu button again to begin recording and return to measure interface.

Figure 12. Input new ID interface

c. In "Input ID" interface,press direction key to choose the existent ID for data storage, here the ID input box will turn to blue,at the same time there will be a red choice box on the first ID in the ID listing as Figure 13. Press direction key to move the red choice box to the ID which is wanted to choose, then press menu button, here prompt information whether to erase the storedata for this ID will appear.Choose "Yes" to erase the ID data, at the same time begin to record and return to the measure interface;choose "No" to return the "Input ID" interface.

Figure 13. Choose existent ID

- d. Every patient ID could save 24-hour data,the device could save 16 patient ID.
- e. In the record state, whatever interfaces the device is in(measure interface, menu interface,clock interface) , if there is no operation in 30 seconds,the display screen will closed automatically.If pressing power button,the device will return the measure interface.
- f. After the display screen closes automatically, the pulse sound will close automatically too for saving electricity under record state.

F. Demo function setting

In the main menu interface, move the choice bar to the "Demo" item,then press the "right button" to turn on demo mode.

G. Review function setting

In the main menu interface,move the choice bar to the "Review" item,then press the "right button" to choose review mode.There are two review mode of value and trend. Under value review mode as Figure 14,press "left button" or "right button" to page up or page down,press "up button" or "down button" to page up or page down quickly.The displayed time on the right top is the total recording time.

ID:HHN	09-1-1 00:01:49		
Time	SpO2	PR	PI
01:56:30	99	077	23.1
01:56:31	99	078	12.1
01:56:32	99	079	23.4
01:56:33	99	080	23.0
01:56:34	99	081	23.1
01:56:35	99	082	12.0
01:56:36	99	083	23.1
01:56:37	99	083	12.1
01:56:38	99	083	23.2
01:56:39	99	082	12.1
01:56:40	98	082	23.2
OK	13:57	Back	

Figure 14. Value review mode

Under trend review mode as Figure 15, the white number on the left bottom is current recording time point, the middle blue number is SpO₂ value, green number is pulse rate value, the white number on the right bottom is PI value. Press "up button" or "down button" to page up or page down; press "left button" or "right button" to move the red triangle.


Figure 15. Trend review mode

H. Close the device

a In the main menu interface, move the choice bar to the "Power off" item, then press menu button to close the device. If the storage function has been opened, the prompt interface of "Recording..." will appear when closing the device. It means that the device is in the record state, can't be closed.

b In the state of boot-strap, long press power button could close the device too.

! If "finger out" appears and you don't turn on the record function and there's no operation in 10 minutes, the device will close automatically.

6.1.4. PC software operation

Please connect the device to the computer by USB data line, which is affiliated with the device, then double click "SpO₂ Assistant" icon to run the PC software. The functions such as uploading record data and change device ID could be carried out by the software. Please refer to <SpO₂ Assistant user manual> for detail.

! If the users choose to turn on the synchronizing display function on computer, it would probably take several seconds for the data to appear on the computer screen. (If there is no data on the computer screen, unplug USB data line, then repeat step t 6.1.4" again).

6.1.5. Charge

There are two kinds of charge method:

A. Connect the device to computer with data line, then the device should be in charge state.

B. Connect the device to power supply with power adaptor, then the device should be in charge state.

When pressing power button, if the device is closed, the charge icon will display, and when the battery status is full, the icon of full power will display. If there is no operation in 30 seconds, the display will close automatically, here short press power button to turn on display. When charging, long press power button to turn on the device, here the battery status icon on the right top will display dynamically. It means that the device is charging up. When the battery status is full, the icon of full power will display.

! If the alarm function is on, the device will provide high-priority alarm signal when the battery is in low power status. Intermittent alarm will occur. High priority indicating requires that operator responds immediately.

6.2. Attention for Operation

- A. Please check the device before use, and confirm that it can work normally.
- B. The finger should be in a proper position (see the attached illustration of Figure 5 for reference), or else it may result in inaccurate measure.
- C. The ray between luminescent tube and photoelectric receiving tube must get across subject's arteriole.
- D. The Oximeter probe should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
- E. Ensure nothing, such as a plaster, can impede the light passage., or else it may result in inaccurate measure of SpO₂, pulse rate and PI.
- F. Excessive ambient light may affect the measure result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight and etc.
- G. Exquisite action of the subject or extreme electrosurgical interference may also affect the accuracy.
- H. Testee can not use enamel or other makeup.
- I. Please clean and disinfect the device after operating according to the User Manual (7.1).

6.3. Clinical Restrictions

- A. As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measure will be more sensitive to interference.
- B. For those with a substantial amount of staining dilution drug (such as methylene blue, indigo green and acid indigo blue), or carbon monoxide hemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic hemoglobin, and some with icterus problem, the SpO₂ determination by this device may be inaccurate.
- C. The drugs like dopamine, procaine, prilocaine, lidocaine and butacaine may also be a major factor resulted in serious error of SpO₂ measure.
- D. As the SpO₂ value serves as a reference value for judgment of anemic anoxia and toxic anoxia, some patients with serious anemia may also report good SpO₂ measure.

7. MAINTAIN, TRANSPORTATION AND STORAGE

7.1. Cleaning and disinfecting

Using medical alcohol to disinfect the device, nature dry or clean it with a soft cloth.

7.2. Maintain

- A. Please clean and disinfect the device before use according to the User Manual (7.1).
- B. Please recharge the battery when the screen shows red aglitter "Battery Too Low" alarm information.
- C. Recharge the battery soon after the over-discharge. The device should be recharged ever six months when it is not regular used. It can extend the battery life following this guidance.
- D. The device needs to be calibrated once a year (or according to the calibrating program of hospital). It also can be performed at the state-appointed agent or just contact us for calibration.

7.3 Transportation and Storage

- A. The packed device can be transported by ordinary conveyance or according to transport contract. The device can not be transported mixed with toxic, harmful, corrosive material.
- B. The packed device should be stored in room with no corrosive gases and good ventilation.
Temperature: -40°C-60°C; Relative Humidity: ≤95%

8. TROUBLESHOOTING

Problem	Possible Cause	Solution
The SpO ₂ and pulse rate can not be displayed normally.	<ol style="list-style-type: none"> 1. The finger is not properly positioned. 2. The patient's SpO₂ is too low to be detected 	<ol style="list-style-type: none"> 1. Place the finger properly and try again. 2. Try again/Go to a hospital for a diagnosis if you are sure the device work properly.
The SpO ₂ and Pulse Rate are not displayed reliably.	<ol style="list-style-type: none"> 1. The finger is not placed inside deep enough. 2. The finger is shaking or the patient is moving. 	<ol style="list-style-type: none"> 1. Place the finger properly and try again. 2. Let the patient keep calm.
The device can not be turned on.	<ol style="list-style-type: none"> 1. The battery is drained away or almost drained away. 2. The malfunction of the device. 	<ol style="list-style-type: none"> 1. Please recharge the battery. 2. Please contact the local service center.
The display has turned off suddenly.	<ol style="list-style-type: none"> 1. The device is damaged. 2. The battery is drained away or almost drained. 	<ol style="list-style-type: none"> 1. Please contact the local service center. 2. Please recharge the battery.
The device can not be used for full time after charging.	<ol style="list-style-type: none"> 1. The battery is not fully charged. 2. The battery is broken. 	<ol style="list-style-type: none"> 1. Please recharge the battery. 2. Please contact the local service center.
The battery can not be full charged even after 10 hours charging time.	The battery is broken.	Please contact the local service center.

9. KEY OF SYMBOLS

Symbol	Description
	Refer to instruction manual/booklet
	The pulse oxygen saturation (%)
	Pulse rate (bpm)
	Perfusion index (%)
	The battery is full
	Lack of battery power (please change batteries in time for exact measure).
	Alarm sound indication off
	Alarm sound indication paused
	Alarm sound indication on

	Pulse sound indication off
	Pulse sound indication on
	Power on/off
	Type BF
IPX0	International Protection (the whole equipment)
IPX1	Internations Protection (only suitable for SpO2 probe)
	USB
HOLD ▶	Touch key switch
RST	Restoration orifice
SN	Serial number
	1. The finger clip falls off (no finger inserted) 2. Probe error 3. Signal inadequacy indicator
	WEEE (2002/96/EC)
	European Representative
	The item is compliant with Medical Device Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community

10. FUNCTION SPECIFICATION

Information	Display Mode
Pulse Oxygen Saturation (SpO2)	2-digit digital TFT display
Pulse Rate (PR)	3-digit digital TFT display
Pulse Intensity	bar-graph TFT display
PI (perfusion Index)	3- digit digital TFT display
SpO2 Parameter Specification	
Measuring rance	0~100% (the resolution is 1%)
Accuracy	70%~100%: ±2%, Below 70% unspecified
Average value	Calculate the average calue in every 4 measure calue. The deviation between average value and the true value dos not exceed 1%
Pulse Parameter Specification	
Measuring Range	30bpm~250bpm, (the resolution is 1bpm)
Accuracy	±2 bpm or ±2%

Average Value	Moving calculate the average pulse rate every 4 cardio-beat's cycle. The deviation between average value and true value doesn't exceed 1%
Perfusion Index Specification	
Range	Continuous bar graph display, the higher display indicates the stronger pulse.
Safety Type	Interior Battery, BF Type
Battery requirement	
Voltage 3.7 rechargeable lithium battery x1 (the red wire on the battery denotes anode, the black wire on the battery denotes cathode).	
Battery working life	
Charge and discharge 500 times.	
Power adapter	
Input Voltage	100 to 240 VAC, 50/60 Hz
Output Voltage	5V DC
Output current	1000mA
Dimensions and weight	
Dimensions	98 (L) x 58 (W) x 12 (H) mm
Weight	About 65g (with a lithium battery)

11. FACTORY DEFAULT

	Default	Unit
Brightness	3	
Key sound	on	
Pulse sound	off	
Alarm sound	off	
Alarm pause time	60	second
SpO ₂ alarm high limit	100	%
SpO ₂ alarm low limit	85	%
Pulse rate alarm high limit	120	bpm
Pulse rate alarm low limit	50	bpm

APPENDIX 1

Guidance and manufacturer's declaration electromagnetic emission for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration -- Electromagnetic emission		
The pulse oximeter (CMS60F) is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The pulse oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The pulse oximeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emission IEC 6100-3-3	Not applicable	

Guidance and manufacturer's declaration electromagnetic immunity for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration-electromagnetic immunity			
The pulse oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the pulse oximeter should assure that it is used in such an environment			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 KV contact ±9 KV air	±6 KV contact ±9 KV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration electromagnetic immunity for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration-electromagnetic immunity			
The pulse oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the pulse oximeter should assure that it is used in such an environment			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3V rms 150kHz to 80 MHz	3Vrms	Portable and mobile RF communication equipment should be used no closer to any part of the pulse oximeter including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Radiated RF IEC 61000-4-3	3V/M 80 MHz to 2.5Ghz	3V/m	<p>Recommended separation distance</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strength from fixed RD transmitters, as determined by and electromagnetic site survey, ^a should be less than the compliance level in each frequency range ^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
---------------------------	-----------------------	------	---

NOTE 1 At 80MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, object and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic site survey should be considered. If the measured field strength in the location in which the pulse oximeter is used exceeds the applicable RF compliance level above, the pulse oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the pulse oximeter.

^b Over the frequency range 150 KHz to 80 MHz, field strengths should be less the 3V/m

Recommended separation distances between portable and mobile RF communications for EQUIPMENT or SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distance between portable and mobile RF communications equipment and the pulse oximeter		
The pulse oximeter is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Pulse oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the pulse oximeter as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.1167	0.2334

0.1	0.3689	0.7378
1	1.1667	2.3334
10	3.6893	7.3786
100	11.6667	23.3334


For transmitters rated at a maximum output power not listed above, the recommended separation distanced in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 at 80mHZ and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objets and people.

APPENDIX 2

State	Alarm condition delay	Alarm signal generation delay
Low voltage alarm	1s	20 ms
SpO2 alarm	330 ms	20 ms
Pulse rate alarm	300 ms	20 ms
Probe error alarm	16 ms	20 ms

	<p>CONTEC MEDICAL SYSTEMS CO., LTD No. 112 Qinhuang West Street, Economic & Technical Development Zone, Qinhuangdao, Hebei Province, P.R.C.</p>
	<p>SHANGHAI INTERNATIONAL HOLDING CORP. GMBH (EUROPE) Eiffestrasse 80, 20537 Hamburg, Germany</p>
	<p>GRUPO R. QUERALTÓ S.A. Polígono Industrial El Pino, Calle Pino Albar, 24, 41016, Sevilla (Spain)</p>

MD CE 0123

Made in PRC