ri-charger[®]

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Instruction for use

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Contents

1. Introduction
 - 1.1. Important information before commissioning.
 - 1.2. Safety information and electromagnetic compatibility
 - 1.3. Packaging symbols
 - 1.4. Intended use
 - 1.4.1 Indication
 - 1.4.2 Contraindication
 - 1.4.3 Intended patient population
 - 1.4.4 Intended operators / users
 - 1.4.5 Required skills / operator training
 - 1.4.6 Environmental conditions
 - 1.5 Warnings / caution
 - 1.6 User responsibility
2. First commissioning of the device
 - 2.1 Scope of delivery
 - 2.2 Device function
3. Operation and function
 - 3.1 Attachment
 - 3.2 Symbol identification
 - 3.3 Commissioning
 - 3.4 Signalling the operating state:
4. Cleaning and disinfection
 - 4.1 General information
5. Technical specifications
6. Care instructions
7. General note
8. Spare parts and accessories
9. Accompanying documents on electromagnetic compatibility according to IEC 60601-1-2
 - 9.1 EMC (electromagnetic compatibility)
 - 9.2 Disposal
10. Warranty

1. Introduction

1.1 Important information before commissioning.

You have purchased a high-quality Riester charging station, which has been manufactured in compliance with Regulation (EU) 2017/745 and is subject to the strictest quality controls at all times. Read through these instructions for use carefully before using the device and keep them in a safe place. If you have any questions, we are available to answer queries at any time. You can find our address in these instructions for use. The address of our sales partner can be obtained on request. Please note that all instruments described in these instructions for use may only be used by appropriately trained personnel. The faultless and safe functioning of this device can only be guaranteed if Riester original parts and accessories are used.

1.2 Safety information and electromagnetic compatibility

Symbol	Symbol note
	Follow the instructions in the user manual.
	Medical device
	Type B applied part
	Class II protective devices
	Warning! The general warning symbol indicates a potentially dangerous situation that can lead to serious injuries.
	Caution! Important note in this manual. The caution symbol indicates a potentially dangerous situation that can lead to minor or moderate injuries. It can also be used to warn of unsafe practices
	Direct current (DC)
	Alternating current (AC)
	Date of manufacture YYMMDD (year, month, day)
	Manufacturer
	Manufacturer's serial number
	Lot/batch number
	Reference number
	Temperature for transportation and storage
	Relative humidity for transport and storage
	Air pressure for transportation and storage.
	Operating ambient air pressure
	CE marking
	Symbol for the marking of electrical and electronic equipment in accordance with Directive 2002/96/EC.
	Caution: used electrical and electronic equipment should not be treated as normal household waste and must be disposed of separately in accordance with national or EU guidelines.
	Non-ionising radiation
	Caution: (US) federal law restricts this device from being sold by or by order of a doctor (licensed physician).

1.3 Packaging symbols

Symbol	Symbol note
	Indicates that the contents of the transport package are fragile and should therefore be handled with care.
	Store in a dry place
	Indicates the correct position for transporting the package.
	Keep away from sunlight
	„Green Dot“ (country-specific)

The device meets the requirements for electromagnetic compatibility. Please note that under the influence of unfavourable field strengths, e.g. during the operation of mobile phones or radiological instruments, adverse effects on function cannot be excluded.

The electromagnetic compatibility of this device has been verified by tests according to the requirements of IEC 60601-1-2.

1.4 Intended use

The **ri-charger®** and **ri-charger® L** have been manufactured for the purpose of charging Riester handles with Riester batteries (**ri-accu®**).

The battery charger serves as an accessory for charging Riester diagnostic instruments that are equipped with **ri-accu®**.

1.4.1 Indication

The **ri-charger®** battery charger serves as an accessory for charging Riester diagnostic instruments that are equipped with **ri-accu®**.

The instruments serve the trained doctor or specialist as an aid in the detection, diagnosis, monitoring, treatment or alleviation of illnesses, injuries and disabilities.

1.4.2 Contraindication

There are no known contraindications.

The device is only intended for the use specified in the instructions for use.

1.4.3 Intended patient population

The device is intended for the user (doctors, nursing staff).

1.4.4 Intended operators/users

The charging stations are only used by doctors in clinics and medical practices.

1.4.5 Required skills/operator training

Since the charging stations are only used by doctors/nursing staff, the relevant qualifications are available.

1.4.6 Environmental conditions

The charging station is intended for use in rooms with a controlled environment.

The charging station must not be exposed to adverse/harsh environmental conditions.

1.5 Warnings / caution

Warning

The general warning symbol indicates a potentially dangerous situation that can lead to serious injuries.

Do not use in a magnetic resonance environment!

There is a risk of ignition of gases if the device is operated in the presence of flammable mixtures or mixtures of medicinal products and air or oxygen or nitrous oxide!

The device must not be operated in rooms in which flammable mixtures or mixtures of pharmaceuticals and air or oxygen or nitrous oxide are present, e.g. operating theatres.

Electric shock!

The housing of the **ri-charger® L** and **ri-charger®** may only be opened by authorised persons.

Damage to the device due to a fall or strong ESD influence!

If the device is not functioning, it must be returned to the manufacturer for repair.

The device must be used in a controlled environment.

The device must not be exposed to harsh environmental conditions.

Caution!

The caution symbol indicates a potentially dangerous situation that can lead to minor or moderate injuries. It can also be used to warn of unsafe practices.

The faultless and safe functioning of this charger can only be guaranteed if Riester original parts and accessories are used.

Old electronic devices must be disposed of in accordance with the institutional guidelines for the disposal of expired devices.

Cleaning frequency and sequence must comply with the cleaning regulations of non-sterile products in the respective facility. Cleaning/disinfection instructions in the instructions for use must be observed.

We recommend disconnecting the charging station's power adapter from the power supply before cleaning or disinfecting.

Clean and disinfect the charging station carefully so that no liquid penetrates the interior.

Never place the charging station in liquids!

The charging stations are delivered in a non-sterile condition. Do not use ethylene, oxide gas, heat, autoclaves or any other harsh methods to sterilise the device. The devices have not been approved for mechanical reprocessing or sterilisation. This leads to irreparable damage!

The patient is not the intended operator.

The product may only be used by qualified personnel.

Qualified personnel are doctors or nurses in hospitals, medical facilities, clinics and medical practices.

1.6 User responsibility

Caution!

User responsibility

You are responsible for the following:

The user must check the integrity and completeness of the charger/instruments before each use. All components must be compatible with each other.

Incompatible components can reduce performance. Never knowingly use a defective device.

Replace parts that are defective, worn, missing or incomplete.

Contact the nearest factory-approved service centre if repairs or replacements are required.

In addition, the user of the device bears sole responsibility for malfunctions resulting from improper use, incorrect maintenance, improper repair, damage or changes by persons other than Riester employees or authorised service personnel.

All serious incidents related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

If there are any problems with the product or use of the product, please contact your doctor immediately.

2. First commissioning of the device

2.1 Scope of delivery

10700 **ri-charger**[®] 2.5 V / 3.5 V / 230 V
with wall mount

- Instructions for use
- Wall mounting material

10701 **ri-charger**[®] 2,5 V / 3,5 V / 120 V
with wall mount

- Instructions for use
- Wall mounting material

10705 **ri-charger**[®]L 3,5 V / 230 V
with wall mount


- Instructions for use
- Wall mounting material

10706 **ri-charger**[®]L 3,5 V / 120 V
with wall mount

- Instructions for use
- Wall mounting material

2.2 Device function

The **ri-charger**[®] battery charger serves as an accessory for charging Riester diagnostic instruments that are equipped with **ri-accu**[®].

- 1 Shaft for battery handle with rechargeable battery
- 2 Removable ear specula storage compartment
- 3 Shaft for battery handle with rechargeable battery
- 4 Green LED (lights up when the battery is fully charged)
- 5 Orange LED (lights up during charging)
- 6 Orange LED (lights up during charging)
- 7 Green LED (lights up when the battery is fully charged)

3. Operation and function

3.1 Attachment

Wall mounting

Position the supplied wall mount on the wall and fix it in place with screws.

The charging station can now be pushed onto the wall mount using the receptacle in the base plate.

3.2 Symbol identification

	Applied part Type B
	Caution! Observe important information in these instructions.
	CE marking

3.3 Commissioning

Preparation for charging

Connect the power adapter to the socket. Make sure that the voltage stated on the type plate on the underside of the device matches the values of the mains connection. After connecting the charger, a self-test is performed. During the test, the diodes on the front of the device light up several times.

If you wish to charge AA handle instruments (handle $\varnothing 19$ mm), please leave the plastic reducers in the handle receptacles. If you wish to charge **ri-scope**[®] instruments (handle $\varnothing 28$ mm), simply pull the reducers upwards to remove them from the handle receptacles.

The **ri-charger**® from Riester can be used to charge 2.5 V and 3.5 V NiMH rechargeable batteries. It does not matter which charging compartment the instruments are inserted into. Of course, it is also possible to charge a 2.5 V and a 3.5 V instrument at the same time.

The **ri-charger**® L from Riester can be used to charge 3.5 V Li-ion rechargeable batteries. It does not matter which charging compartment the instruments are inserted into. Of course, it is also possible to charge two 3.5 V instruments at the same time.

Charging process with **ri-charger**® and **ri-charger**® L

The rechargeable battery is checked when it is inserted into the charging compartment.

If the battery is deeply discharged or defective, an attempt is made to regenerate the battery with a pulsed charging current.

Meanwhile, the yellow charging diode flashes.

If the battery voltage is above the deep discharge value, charging begins.

The yellow charging diode lights up.

Charging time

The end of the charging time is detected by evaluating the negative voltage difference that is caused by the battery cells heating up when fully charged. The yellow charging diode goes out and the green battery diode lights up.

The max. charging time for NiMH rechargeable batteries is limited to <650 minutes.

The max. charging time for Li-ion rechargeable batteries is limited to <480 minutes.

When the battery is fully charged (green diode lights up), the device switches to trickle charging.

Caution!

If an instrument with a fully charged battery (green diode lights up) is removed from the charging compartment for a few seconds without it being switched on, the yellow charging diode lights up when the unused battery is reinserted. However, this is not because of a charger malfunction, rather the described full battery detection.

When the battery is fully charged, no negative voltage difference can be detected.

The battery must first warm up slightly to signal „battery full“.

The time until the green full battery diode lights up again can vary greatly (15 minutes and longer). However, the additional temperature detection and timer ensure that the charging process is carefully monitored at all times.

We recommend placing the instruments back in the charging compartments after use to ensure that the instruments are always ready.

3.4 Signalling the operating state:

Operating state	LED 4 and 5	LED 6 and 7
No handle with rechargeable battery in the charging compartment	Off	Off
Charging active	LED 5 on, LED 4 off	LED 6 on, LED 7 off
Charging complete	LED 4 on, LED 5 off	LED 7 on, LED 6 off
Error: deeply discharged battery	LED 4 off, LED 5 on/ off @ 5Hz	LED 6 on/off @ 5Hz, LED 7 off
Error: temperature monitoring	LED 4 off, LED 5 on/ off @ 5Hz	LED 6 on/off @ 5Hz, LED 7 off
Error: precharge timeout	LED 4 off, LED 5 on/ off @ 5Hz	LED 6 on/off @ 5Hz, LED 7 off
Error: charge timeout	LED 4 off, LED 5 on/ off @ 5Hz	LED 6 on/off @ 5Hz, LED 7 off

4. Cleaning and disinfection

4.1 General information

The cleaning and disinfection of medical devices protects patients, users and third parties and preserves the value of medical devices. Due to the product design and the materials used, it is not possible to set a maximum limit for reprocessing cycles. The service life of a medical device depends on its function and use. Before returning defective products for repair, you should follow the instructions below.

Caution!

We recommend disconnecting the charging station's power adapter.

Clean and disinfect the charging station with care. The exterior of the charging station can be cleaned with a damp cloth until it is visually clean. Only use disinfectant products in accordance with the manufacturer's instructions. Only disinfectants that have been deemed effective in accordance with the national guidelines should be used. After disinfection, wipe the instruments with a damp cloth to remove any remaining disinfectant before cleaning or disinfection.

Never place the charging station in liquids!

The charging station is delivered in a non-sterile condition. Do not use ethylene, oxide gas, heat, autoclaves or any other harsh methods to sterilise the device.

The devices must not be mechanically serviced or sterilised. This can lead to irreversible damage!

If a reusable device shows signs of material deterioration, it should no longer be reused and should be disposed of/claimed according to the procedure described in the Disposal/Warranty sections.

5. Technical specifications

Technical specifications	
Medical device:	Medical device for power supply (charging of battery-operated instruments)
Electric shock:	Class II insulation devices
Model	ri-charger®L for Li-ion rechargeable batteries ri-charger® for NiMH rechargeable batteries
Power adapter	Input: 100 V-240 V AC / 50-60 Hz/0.6 A Output: 9 V DC/2 A/18 W
Classification	Applied part type B
Operating conditions	0°C to + 40°C, 10% to 85% relative humidity
Storage and transport conditions	-5°C to + 50°C, 10% to 85% relative humidity
Air pressure	700 to 1050 hPa
Dimensions	186 x 87 x 74 mm
Weight	500 g

6. Care instructions

7. General note

The cleaning and disinfection of the medical device serves to protect the patient, the user and third parties and to preserve the value of the medical devices.

Due to the product design and the materials used, there is no defined upper limit for possible reprocessing cycles. The service life of medical devices is determined by their function and careful handling. Defective products must have undergone the entire reprocessing process before being returned for repair.

8. Spare parts and accessories

10711 **ri-charger®** AA type handle insert, black

10710 **ri-charger®** wall mount with screws and dowels

10712 **ri-charger®** cover for ear specula compartment

9. Accompanying documents on electromagnetic compatibility according to IEC 60601-1-2

The device meets the requirements for electromagnetic compatibility. Please note that under the influence of unfavourable field strengths, e.g. during the operation of mobile phones or radiological instruments, adverse effects on function cannot be excluded.

The electromagnetic compatibility of this device has been verified by tests according to the requirements of IEC 60601-1-2.

9.1 EMC (electromagnetic compatibility)

- 9.1.1** Follow the instructions below during the installation and operation of the device:
- 9.1.2** To avoid electromagnetic interference with the operation of the device, do not use the device at the same time as other electronic devices.
- 9.1.3** To avoid electromagnetic interference with the operation of the device, do not use or stack the device near, on or under other electronic devices.
- 9.1.4** Do not use the device in the same room as other electronic devices, e.g. life-supporting devices that have a significant impact on the life of the patient and treatment outcomes, or other measuring instruments or treatment devices that use little electrical power.
- 9.1.5** Do not use cables or accessories that are not specified for the device, as this can increase the emission of electromagnetic waves from the device and reduce the device's electromagnetic immunity.

Caution!

Medical electrical devices are subject to special precautions with regard to electromagnetic compatibility (EMC).

Portable and mobile radio frequency communication devices may affect medical electrical equipment. The ME device is intended for use in an electromagnetic environment and for professional facilities such as industrial areas and hospitals.

The user of the device should ensure that it is used within such an environment.

Warning!

The ME device may not be stacked, arranged or used directly next to or with other devices. If the ME device must be operated near or stacked with other devices, the ME device and the other ME devices must be monitored to ensure proper operation in this arrangement. This ME device is only intended for use by medical professionals. This ME device is intended for use in professional healthcare facilities. This device may cause radio interference or impair the operation of nearby devices. It may become necessary to take appropriate corrective measures, such as redirecting or rearranging the ME device or shield.

The rated ME device does not exhibit any basic performance features in the sense of IEC 60601-1 that would present an unacceptable risk to patients, operators or third parties should the power supply fail or malfunction.

Warning!

Portable RF communications equipment (radios), including accessories such as antenna cables and external antennas, should not be used in closer proximity than 30 cm (12 inches) to parts and cables of the **ri-charger® L** and **ri-charger®** specified by the manufacturer. Failure to comply may result in reduced device performance.

Guidance and manufacturer's declaration - electromagnetic emissions

The chargers are intended for use in the electromagnetic environment specified below. The customer or user of the charger should ensure that it is used in such an environment.

Guidance and manufacturer's declaration - electromagnetic emissions

The charging station is intended for use in the electromagnetic environment specified below. The customer of the user of the charging station should ensure that it is used in such an environment.

Emissions testing	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The charging station uses RF energy for its internal function only. Therefore, the RF emissions are very low and are unlikely to cause interference in nearby electronic devices.
RF emissions CISPR 11	Class B	The charging station is intended for use in all facilities, including residential areas and those directly connected to a public supply network that also supplies buildings used for residential purposes.
Harmonic emissions IEC 61000-3-2	Pass	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Pass	

Guidance and manufacturer's declaration - electromagnetic immunity

The charging station is intended for use in the electromagnetic environment specified below. The customer or user of the charging station should ensure that it is used in such an environment.

Immunity testing	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Con: ± 8 kV Air: ± 15 kV	Con: ± 8 kV Air: ± 15 kV	Floors should be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%.
Fast transient electrical disturbances / bursts IEC 61000-4-4	5/50 ns, 100 kHz; ± 2 kV	5/50 ns, 100 kHz; ± 2 kV	The quality of the supply voltage should be that of a typical business or hospital environment
Surge IEC 61000-4-5	1,2/50 (8/20) μ s LtL: $\pm 1,0$ kV LtG: $\pm 2,0$ kV	1,2/50 (8/20) μ s LtL: $\pm 1,0$ kV LtG: $\pm 2,0$ kV	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short interruptions and voltage fluctuations on power supply input lines IEC 61000-4-11	0% UT for 0.5 cycle (1 phase) 0% UT for 1 cycle 70% UT for 25/30 cycles (50/60 Hz)	0% UT for 0.5 cycle (1 phase) 0% UT for 1 cycle 70% UT for 25/30 cycles (50/60 Hz)	The quality of the supply voltage should be that of a typical business or hospital environment.
	0% UT for 250/300 cycles (50/60 Hz)	0% UT for 250/300 cycles (50/60 Hz)	
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz	30 A/m 50 Hz	Mains frequency magnetic fields should be at a level characteristic of a typical location in a typical commercial hospital environment.
Hinweis U_T is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity

The charging station is intended for use in the electromagnetic environment specified below. The customer or user of the charging station should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	Pass	Pass	<p>Portable and mobile RF communications equipment should be used no closer to any part of the charging station, including the cables, than the recommended distance calculated using the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d = 1.2\sqrt{P}$ 150 KHz to 80 MHz $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2,7 GHz</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	<p>Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be below the compliance level in each frequency range. In the vicinity of devices marked with the following symbol, interference may occur:</p>
Proximity fields from RF wireless communications equipment	Pass	Pass	<p>Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and the recommended distance is indicated in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be below the compliance level in each frequency range.</p> <div style="text-align: center;">  </div>

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. The electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

- a.)
Field strengths from fixed transmitters, such as base stations for radio (mobile/cordless phones) and land mobile radios, amateur radio, AM and FM radio broadcasts and television broadcasts, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location where the charging station is used exceeds the RF compliance level specified above, the charging station should be monitored to verify normal operation. If abnormal performance is observed, additional measures may be required, such as re-orienting or relocating the charging station.
- b.)
With a frequency range over 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.

Recommended distances between portable and mobile RF communications equipment and the charging stations.

The charging station is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the charging station can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the charging station as recommended below, according to the maximum output power of the communications equipment.

Rated output power of the transmitter (W)	Separation distance according to the frequency of the transmitter (m)		
	150 kHz to 80 MHz $d = 12\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,7 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. The electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

9.2 Disposal

The used medical device must be disposed of in accordance with current medical practices or local regulations on the disposal of infectious biological medical waste.

Batteries and electrical/electronic devices may not be treated as household waste and must be disposed of in accordance with local regulations.

If you have any questions about the disposal of products, please contact the manufacturer or their representative.

10. Warranty

This product has been manufactured to the strictest quality standards and has undergone a thorough final quality check before leaving our factory. We are therefore pleased to issue a warranty of 2 years from the date of purchase on all defects that are demonstrably due to material or manufacturing defects. A warranty claim does not apply in case of improper handling. All defective parts of the product will be replaced or repaired free of charge within the warranty period. This does not apply to wear parts. We offer an additional 5-year warranty on the R1 shock-proof calibration, which is required for CE certification. A warranty claim can only be guaranteed if this warranty card has been filled in and stamped by the dealer and is enclosed with the product. Please note that all warranty claims must be made within the warranty period. Of course, we will be happy to carry out inspections or repairs for a fee after the warranty period has expired. You are also welcome to request a preliminary cost estimate from us free of charge. In the event of a warranty claim or repair, please send the Riester product with the completed warranty card to the following address:

Rudolf Riester GmbH
Dept. Repairs RR
Bruckstr. 31
72417 Jungingen
Germany

Serial number or batch number
Date, stamp and signature of the specialist dealer

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