ri-pen[®]

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Table of Contents

- 1. Important information prior to use**
 - 1.1. Safety symbols**
 - 1.2. Packaging symbols**
 - 1.3. Contraindications**
 - 1.4. Intended patient population**
 - 1.5. Intended operator/user**
 - 1.6. Required skills/training**
 - 1.7. Environmental conditions**
 - 1.8. Warnings / caution**
 - 1.9. Scope of delivery**
- 2. Examination lamp/device function**
 - 2.1. Purpose**
 - 2.2. Operational readiness**
 - 2.3. Specifications**
 - 2.4. Caution**
 - 2.5. Disposal**
 - 2.6. Care instructions**
 - 2.7. Cleaning and disinfection**
 - 2.8. Caution**
 - 2.9. Maintenance**
- 3. EMC requirements**
- 4. Disposal**
- 5. Warranty**

1. Important information prior to use

You have purchased a high-quality Riester diagnostic penlight, which was manufactured in accordance with Regulation (EU) 2017/75 on medical devices and is subject to the strictest quality controls at all times. The excellent quality guarantees reliable diagnoses. Please read the instructions carefully before use and keep them for future reference. If you have any questions, we are available at any time, and our contact information is provided at the end of this IFU. Contact information for distribution and sales partners can be obtained upon request. Please note all instruments described in this IFU should only be used by appropriately trained personnel. The safe functioning of Riester diagnostic instruments is only guaranteed if Riester original parts and accessories are used.

1.1.Safety symbols

Symbol	Note on symbol
	Follow the instructions in the user manual.
	Type B applied part
	Medical device
	Class II protective devices
	Warning! The general warning symbol indicates a potentially dangerous situation that can lead to serious injuries.
	Caution! Important note in this manual. The caution symbol indicates a potentially dangerous situation that can lead to minor or moderate injuries. It can also be used to warn of unsafe practices.
	Direct current
	Alternating current
	Date of manufacture YYYY-MM-DD / (Year-Month-Day)
	Manufacturer
	Lot / batch number
	Reference number
	Temperature for transportation and storage
	Relative humidity for transportation and storage
	CE-marking
	Symbol for the marking of electrical and electronic equipment in accordance with Directive 2002/96/EC. Caution: Used electrical and electronic equipment should not be treated as normal household waste, but should be disposed of separately in accordance with national and EU regulations.
	Non-ionizing radiation

1.2 Packaging symbols

Symbol	Note on symbol
	Indicates that the contents of the transport package are fragile and should therefore be handled with care.
	Store in a dry place
	Indicates the correct position for transporting the package.
	Keep away from sunlight
	"Green Dot" (country-specific)

Warning:

Please note the proper and safe functioning of our instruments is only guaranteed if both the instruments and their accessories are exclusively from Riester. The use of other accessories may result in increased electromagnetic emissions or reduced electromagnetic immunity of the device and may lead to incorrect operation.

1.3. Caution /contraindications

- There may be a risk of ignition of gases if the instrument is used in the presence of flammable mixtures or mixtures of pharmaceuticals.
- The diagnostic penlight must never be placed in liquids.
- Extended exposure to light during an eye exam with the diagnostic penlight may damage the retina.
- Only use Riester or Riester-approved accessories/consumables.
- Cleaning frequency and sequence must comply with the cleaning regulations for non-sterile products in the respective facility.
- Cleaning/disinfection instructions in this IFU must be observed.
- The product may only be used by trained personnel.

1.4. Intended patient population

The device is intended for adults and children.

1.5. Intended operator/user

Diagnostic penlights are intended exclusively for use by doctors and clinicians (trained medical professionals) in clinics and medical practices.

1.6. Required skills/training

Only trained medical professionals should operate the diagnostic penlights, as they have the required skills and qualifications.

1.7. Environmental conditions

The device is intended for use in rooms with a controlled environment. The device must not be exposed to adverse/harsh environmental conditions.

1.8. Warnings/caution

Warning

The general warning symbol indicates a potentially dangerous situation that can lead to serious injuries.

Do not use in a magnetic resonance environment!

There is a risk of ignition of gases if the device is operated in the presence of flammable mixtures or mixtures of medicinal products and air or oxygen or nitrous oxide!

The device must not be operated in rooms in which flammable mixtures or mixtures of pharmaceuticals and air or oxygen or nitrous oxide are present, e.g. operating rooms.

Damage to the device due to a fall or strong ESD influence!

If the device is not functioning, it must be returned to the manufacturer for repair.

The device must be used in a controlled environment.

The device must not be exposed to harsh environmental conditions.

Caution!

The caution symbol indicates a potentially dangerous situation that can lead to minor or moderate injuries. It can also be used to warn of unsafe practices.

The faultless and safe functioning of the examination lamps can only be guaranteed if Riester original parts and accessories are used.

Old electronic devices must be disposed of in accordance with the institutional guidelines for the disposal of expired devices.

Cleaning frequency and sequence must comply with the cleaning regulations of non-sterile products in the respective facility. Cleaning/disinfection instructions in the instructions for use must be observed.

We recommend removing the (rechargeable) batteries from the battery handle before cleaning or disinfecting.

Clean and disinfect the diagnostic penlights carefully so that no liquid penetrates the interior.

Never place the devices in liquids!

The instruments are delivered in a non-sterile condition. Do not use ethylene oxide gas, heat, autoclaves or any other harsh methods to sterilise the device.

The devices have not been approved for mechanical reprocessing or sterilisation. This leads to irreparable damage!

The patient is not the intended operator.

The product may only be used by qualified personnel.

Qualified personnel are doctors or nurses in hospitals, medical facilities, clinics and medical practices.

Caution!

User Responsibility

It is your responsibility to:

The user must check the integrity and completeness of the instruments before each use. All components must be compatible with each other.

Incompatible components can result in degraded performance.

Never knowingly use a defective device.

Replace parts that are defective, worn, missing or incomplete.

Contact the nearest factory-approved service centre if repairs or replacements are required.

In addition, the user of the device bears sole responsibility for malfunctions resulting from improper use, incorrect maintenance, improper repair, damage or changes by persons other than Riester employees or authorised service personnel.

All serious incidents related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is domiciled.

If there are any problems with the product or use of the product, please contact your doctor immediately.

1.9. Scope of delivery

ri-pen® diagnostic penlight , 1 piece

Art.no.: 5075-526 black

Art.no.: 5074-526 silver

Art.no.: 5077-526 red

Art.no.: 5071-526 blue

Art.no.: 5078-526 purple

Art.no.: 5072-526 green

ri-pen® diagnostic penlight, 6 pieces

Art.no.: 5070 assorted colours

Art.no.: 5075 black

Art.no.: 5074 silver

Art.no.: 5077 red

Art.no.: 5071 blue

Art.no.: 5078 purple

Art.no.: 5072 green

2. ri-pen diagnostic penlight

Device function:

1) On-off button

2) Pocket Clip

3) Battery Compartment

4) Lamp head, LED

2.1. Purpose

The Riester ri-pen was manufactured to illuminate and examine various parts of the body, including pupil examinations.

2.2. Operational readiness

The ri-pen® diagnostic penlight is equipped with an LED lamp. Twist the silver LED base at the bottom of the device and insert two standard AAA batteries with the positive terminal towards the LED.

2.3. Specifications:

Technical data for the bulb

Light field	LED 3V 5300 Kelvin Max. 40,000 lux Approx. 40 mm Ø/distance 15 cm
Classification	Applied part
Batteries	2 x AA alkaline batteries, 1.5 V
Ambient temperature:	0° to +40°
Relative humidity:	30% to 70% non-condensing
Transport and storage temperature:	-10 ° to +55 °
Relative humidity:	10% to 95% non-condensing

2.4. CAUTION:

- If you do not use the device for a long time or take it with you while travelling, please remove the batteries from the handle.
- New batteries should be inserted when the light intensity of the instrument becomes dim and could impair the examination.
- For optimal light output, we recommend use of new high-quality batteries when changing the battery.
- Make sure that no liquid or moisture penetrates the handle.

2.5. Disposal:

Please note that batteries must be disposed of specially. Information about this can be obtained from your municipality or from your environmental consultant.

Powering on and off

An on-off button is located on the handle, above the pocket clip. Press the button down to turn the device on; press the button again to turn the device off.

2.6. Care instructions

General note

The cleaning and disinfecting of medical devices serves to protect the patient, user, and third parties, and also to maintain the overall value (integrity & functionality) of the medical devices.

Due to the product design and materials used, a defined limit for the maximum possible number of reprocessing cycles cannot be determined. The service life of medical devices is defined by their designed and intended function and appropriate, careful handling.

Before returning for repair, defective products must have undergone the prescribed cleaning and disinfection process.

2.7. Cleaning and disinfection

In order to avoid possible cross-contamination, the examination lamps must be cleaned and disinfected regularly.

The diagnostic penlight can be cleaned on the outside using a damp cloth (if necessary, moistened with alcohol) until they are visually clean. Wipe with disinfectant (e.g. disinfectant Bacillol AF from Bode Chemie GmbH (time 30s)) ac-

ording to the instructions of use of the respective manufacturer of disinfectant. Only disinfectants with proven effectiveness according to national guidelines should be used. After disinfecting, wipe the diagnostic penlight with a damp cloth to remove potential residue.

Please make sure that the cloth is moistened but not wet, so no moisture penetrates the openings of the examination lamp.

Make sure that glass and lenses are only cleaned with a dry and clean cloth.

Caution! 

The diagnostic instruments are not sterile devices; they cannot be sterilized.

Caution! 

Never place the examination lamps in liquids! Make sure that no liquids penetrate the housing interior! The article is not approved for machine reprocessing and sterilisation.

This can lead to irreparable damage!

If a reusable device shows signs of material deterioration, it should no longer be used and should be disposed of according to the procedures described in the Disposal/Warranty section.

2.9. Maintenance

The instruments and their accessories require no special maintenance.

If an instrument needs to be tested for any reason, please contact the Riester office directly, or contact an authorized Riester dealer in your area, the details of which we will provide you upon request.

CAUTION! 

There may be a risk of ignition if the device is operated in the presence of flammable mixtures of pharmaceuticals and air, oxygen, nitrous oxide or anaesthetic gases! Safety information according to international standard IEC 60601-1 „Medical electrical equipment - Part 1: General requirements for basic safety and essential performance“: It is not permitted to open the battery handle in the vicinity of the patient or to touch batteries and the patient simultaneously.

3. Electromagnetic compatibility according to IEC 60601-1-2, 2014, Ed. 4.0

Caution: 

Medical electrical equipment is subject to special precautions regarding electromagnetic compatibility (EMC).

Portable and mobile radio frequency communication devices can affect medical electrical equipment. The ME device is intended for operation in a home health care electromagnetic environment and for professional facilities such as industrial areas and hospitals.

The user of the device should ensure that it is operated within such an environment.

Warning: 

The ME device may not be stacked, arranged or used directly next to or with other devices. When use close to or stacked with other devices is required, the ME device and the other ME devices must be monitored to ensure intended operation within this configuration. This ME device is intended for use by medical professionals only. This device may cause radio frequency interference or interfere with the operation of nearby devices. It may become necessary to take appropriate corrective measures, e.g. re-positioning, re-configuring or shielding the ME device.

The rated ME device does not show basic performance features in terms of EN60601-1, the malfunctioning of which, or rather present an unacceptable risk to patients, operators or third parties in the event that it should present an outage or the power supply should disconnect.

Warning: 

Portable RF communications equipment (radios), including accessories, such as antenna cables and external antennas, should not be used in closer proximity than 30 cm (12 inches) to parts and cables of the examination lamp specified by the manufacturer. Failure to comply may result in a reduction of the device's performance features.

Guidelines and manufacturer's declaration - electromagnetic emissions		
The ri-pen instrument is intended for use in the electromagnetic environment specified below. The customer or user of the ri-pen should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions RF emissions pursuant to CISPR 11	Group 1	The ri-pen uses RF energy exclusively for internal functions. Therefore, its RF emissions are very low and unlikely to interfere with nearby electronic devices.
RF emissions RF emissions according to CISPR 11	Class B	The ri-pen is intended for use in all establishments, including residential areas and those directly connected to a public supply network that also supplies buildings used for residential purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Emissions of voltage fluctuations, flicker IEC 61000-3-3	Not applicable	

Guidelines and manufacturer's declaration – electromagnetic immunity

The ri-pen instrument is intended for use in the electromagnetic environment specified below. The customer or user of the ri-pen should ensure that it is used in such an environment.

Immunity testing	IEC 60601 test level	Compliance	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Con: ± 8 kV Air: $\pm 2, 4, 8, 15$ kV	Con: ± 8 kV Air: $\pm 2, 4, 8, 15$ kV	Floors should be wood, concrete or ceramic tile. If the floor is covered with synthetic material, the relative humidity must be at least 30%.
Fast transient electrical disturbances / bursts IEC 61000-4-4	5/50 ns, 100 kHz, ± 2 kV	Not applicable	The quality of the supply voltage should be that of a typical commercial or hospital environment.
Impulse voltage IEC 61000-4-5	± 0.5 kV voltage Phase-to-phase conductor ± 2 kV voltage Outer conductor to earth	Not applicable	The quality of the supply voltage should be that of a typical commercial or hospital environment.
IEC 61000-4-11 Voltage dips, short interruptions and voltage variations according to IEC 61000-1-11	$\leftarrow 0\%$ UT 0.5 period at 0, 45, 90, 135, 180, 225, 270 and 315 degrees 0% UT 1 period and 70% UT 25/30 periods Single-phase: at 0 degrees (50/60 Hz)	Not applicable	The quality of the supply voltage should be that of a typical commercial or hospital environment.
Magnetic field with energy-efficient rated frequencies IEC 61000-4-8	30A/m 50/60 Hz	30A/m 50/60 Hz	Mains frequency magnetic fields should be at a level characteristic of a typical location in a typical commercial hospital environment.

Note: U_i is the AC source. Mains voltage before the application of the test level.

Guidelines and manufacturer's declaration – electromagnetic immunity

The ri-pen instrument is intended for use in the electromagnetic environment specified below. The customer or user of the ri-pen should ensure that it is used in such an environment.

Immunity testing	IEC 60601 test level	Compliance	Electromagnetic environment - guidance
<p>Conducted RF disturbances pursuant to IEC 61000-4-6</p>	<p>3 Vrms 0.5 MHz to 80 MHz 6 V in ISM frequency bands Between 0.15 MHz and 80 MHz 80% AM at 1 kHz</p>	<p>Not applicable</p>	<p>Portable and mobile RF communications equipment should not be used closer to any part of the non-contact ri-pen, including the cables, than the recommended distance, which is calculated using the equation applicable to the transmitter frequency. Recommended separation distance:</p> <p>$d = 1.2 \times P$ 80 MHz to 800 MHz $d = 2.3 \times P$ 800 MHz to 2.7 GHz</p> <p>Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended distance in metres (m). Field strengths of fixed RF transmitters, determined by an electromagnetic site survey, should be lower than the compliance standard in each frequency range. Interference may occur in the vicinity of devices marked with the following symbol:</p> 

Emitted RF	3 V/m		
IEC 61000-4-3	80 MHz to 2.7 GHz	10 V/m 27 V/m	
Proximity fields of wire-less RF commu-nications equipment	380 – 390 MHz	28 V/m	
	27 V/m, PM 50%, 18 Hz	9 V/m	
	430 – 470 MHz	28 V/m	
	28 V/m; (FM ±5 kHz, 1 kHz sine)	28 V/m	
	PM, 18 Hz	9 V/m	
	704 – 787 MHz		
	9 V/m, PM 50%, 217 Hz		
	800 – 960 MHz		
	28 V/m, PM 50%, 18 Hz		
	1700 – 1990 MHz		
28 V/m, PM 50%, 217 Hz			
2400 - 2570 MHz			
28 V/m, PM 50%, 217 Hz			
5100 - 5800 MHz			
9 V/m, PM 50%, 217 Hz			

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. The electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

a: Field strengths of fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM broadcast and television transmission, cannot be theoretically predicted with accuracy. To evaluate the electromagnetic environment based on fixed RF transmitters, an electromagnetic assessment should be considered. If the measured field strength at the location where the ri-pen is used exceeds the above-mentioned RF degree of compliance, the ri-pen should be observed to ensure normal operation.

If abnormal performance is observed, additional measures may be required, such as reorienting or moving the ri-pen.

b: With a frequency range of 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.

Recommended distances between portable and mobile RF communications equipment and the ri-pen.

The ri-pen is intended for use in an electromagnetic environment in which RF emissions are controlled. The customer or user of the ri-pen can help to avoid electromagnetic interference by observing the minimum distance between portable and mobile RF communications equipment (transmitters) and the ri-pen in accordance with the maximum output power of the communication equipment.

Rated maximum output power of the transmitter (W)	Separation distance according to the frequency of the transmitter (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,7 GHz
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended distance in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
 Note 1: At 80 MHz and 800 MHz, the separation distance applies to the higher frequency range.

Note 2: These guidelines may not apply in all situations. The electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

4. Disposal

The used medical device must be disposed of in accordance with current medical practices or local regulations on the disposal of infectious biological medical waste.

Batteries and electrical/electronic devices may not be treated as household waste and must be disposed of in accordance with local regulations.

If you have any questions about the disposal of products, please contact the manufacturer or their representative.

5. WARRANTY

This product was manufactured to the highest quality standards and subjected to a thorough final inspection before leaving our factory.

We are pleased to issue a warranty of 2 years from the date of purchase on all defects traceable to material or manufacturing defects. A warranty claim is excluded from cases of improper handling or use.

All defective parts will be replaced or repaired free of charge within the warranty period.

A warranty claim can only be made if the product is accompanied by this warranty card, which is filled out in full and stamped by the dealer. Please note that warranty claims must be made within the warranty period.

After the warranty period, we are happy to provide device repairs, maintenance, quality and functionality checks; however, there would be a fee for such services. We offer free, no-obligation quotes for requested services.

In case of warranty coverage or repair, we ask you to return the Riester product with the completed warranty card to the following address:

Rudolf Riester GmbH
Repairs dept. RR
Bruckstr. 31
D-72471 Jungingen
Germany

Serial number or batch number:

Date, stamp and signature of the specialist dealer

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