- DE** Gebrauchsanweisung  
Laryngoskope  
Warm und Kaltlicht
- GB** Instructions  
Laryngoscopes  
Standard and Fiber Optic
- FR** Mode d'emploi  
Laryngoscopes  
Lumière standard et à fibre optiques
- ES** Instrucciones para el uso  
Laringoscopios  
de luz estándar y de fiber óptica
- IT** Istruzioni per l'uso  
Laringoscopi  
standard e a fibre ottiche
- RU** инструкция по эксплуатации  
ларингоскопы стандартные и  
фиброоптические

## DEUTSCH

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## Table of Contents

<b>1.</b>	<b>Introduction</b>
<b>1.1</b>	<b>Important information prior to use</b>
<b>1.2</b>	<b>Safety symbols</b>
<b>1.3</b>	<b>Packaging symbols</b>
<b>1.4</b>	<b>Purpose</b>
<b>1.4.1</b>	<b>Indications</b>
<b>1.4.2</b>	<b>Contraindications</b>
<b>1.4.3</b>	<b>Intended patient population</b>
<b>1.4.4</b>	<b>Intended operators/users</b>
<b>1.4.5</b>	<b>Required skills/operator training</b>
<b>1.4.6</b>	<b>Environmental conditions</b>
<b>1.5</b>	<b>Warnings/caution</b>
<b>2.</b>	<b>First use</b>
<b>2.1</b>	<b>Scope of delivery</b>
<b>2.2</b>	<b>Device function</b>
<b>3.</b>	<b>Operation and function</b>
<b>3.1</b>	<b>Symbol identification</b>
<b>3.2</b>	<b>Startup</b>
<b>3.3</b>	<b>Inserting and removing batteries and rechargeable batteries</b>
<b>3.4</b>	<b>Attachment of laryngoscope blades</b>
<b>4.</b>	<b>Care instructions</b>
<b>4.1</b>	<b>General information</b>
<b>4.2</b>	<b>Cleaning and disinfection</b>
<b>4.2.1</b>	<b>Laryngoscope blades/handles</b>
<b>4.2.2</b>	<b>Point of use:</b>
<b>4.2.3</b>	<b>Containment and transport:</b>
<b>4.2.4</b>	<b>Preparation for decontamination:</b>
<b>4.2.5</b>	<b>Blades:</b>
<b>4.2.6</b>	<b>Handle</b>
<b>4.2.7</b>	<b>Pre-cleaning</b>
<b>4.2.8</b>	<b>Cleaning: manually</b>
<b>4.2.9</b>	<b>Disinfection: manually</b>
<b>4.2.10</b>	<b>Cleaning and disinfection: Automated WD</b>
<b>4.2.11</b>	<b>Maintenance, inspection and testing:</b>
<b>4.2.12</b>	<b>Packaging:</b>
<b>4.2.13</b>	<b>Sterilisation:</b>
<b>4.2.14</b>	<b>Storage:</b>
<b>4.2.15</b>	<b>Additional information</b>
<b>4.3</b>	<b>Test procedure for blade and handle</b>
<b>5.</b>	<b>Technical specifications</b>
<b>5.1</b>	<b>Technical data for the bulb</b>
<b>6.</b>	<b>Spare parts and accessories</b>
<b>7.</b>	<b>Maintenance/accuracy check/calibration/applied standards</b>
<b>7.1</b>	<b>Applied standards</b>
<b>8.</b>	<b>Disposal</b>
<b>9.</b>	<b>Warranty</b>

**1. Introduction****1.1 Important information prior to use**

You have purchased a high-quality Riester product, which has been manufactured in compliance with Regulation [EU] 2017/745 and is subject to the strictest quality controls at all times. The excellent quality guarantees reliability in use.

These instructions for use are valid for laryngoscopes with fibre optics, standard illumination and ri-dispo (disposable) blades. Should you have any questions, please contact us or your Riester representative at any time.

Our contact details are listed on the last page of this user manual. We will gladly provide you with the address of Riester sales and distribution representatives upon request. Please note that all instruments described in this user manual are only suitable for use by appropriately trained personnel.

Before each use, check the status of the internal electrical power source by powering on the device.

Please also note the proper and safe functioning of our instruments is only guaranteed if both the instruments and their accessories are exclusively from Riester.

## 1.2 Safety symbols

Symbol	Note on symbol
	Meaning of the symbol on the outer packaging: Caution, please observe the operating instructions
	The operator is obliged to read the instructions of the operating manual
	Application part type B
<b>MD</b>	Medical device
	Protection class II
	<b>Warning!</b> The general warning symbol indicates a potentially dangerous situation that can lead to serious injuries.
	<b>Caution!</b> The caution symbol indicates a potentially dangerous situation that can lead to minor or moderate injuries. The symbol may also indicate unsafe practices.
	Date of manufacture YYYY-MM-DD / (Year-Month-Day)
	Manufacturer
<b>LOT</b>	Manufacturer's batch number
	Temperature requirements for transport and storage
	Relative humidity for transport and storage
CE	CE Mark
<b>S/S</b>	Stainless Steel
	Stainless steel
	Not sterile
	Latex-free
	<b>Caution:</b> Used electrical and electronic equipment should not be treated as normal domestic waste but should be disposed of separately in accordance with national or EU directives
	Recycling symbol lithium ion rechargeable battery
	Month of manufacture of the battery

### 1.3 Packaging symbols

Symbol	Note on symbol
	Fragile. The package should be handled with care.
	Keep the package from getting wet
	This way up. The symbol indicates the correct positioning for transporting the package.
	Keep away from sunlight
	„Green Dot“ (country-specific)

### 1.4 Purpose

The laryngoscopes are made for endotracheal intubation.

The Riester battery handles are used to supply the lamps with energy and to hold a blade.

#### 1.4.1 Indications

The laryngoscopes (med. laryngeal mirrors) from Riester are inserted into the oral cavity as far as the pharynx, they are used for endotracheal intubation.

#### 1.4.2 Contraindications

Any other use, or use going beyond this, is considered inappropriate. The manufacturer is not liable for any damages resulting from such. The user bears the sole risk.

#### 1.4.3 Intended patient population

The device is intended for all patients.

#### 1.4.4 Intended operator/user

The device is for use trained and licensed clinicians in hospitals, medical institutions, clinics, doctor's offices.

#### 1.4.5 Required skills/operator training

The user must have the qualifications of a doctor.

Laryngoscopes are a day-to-day tool of the anaesthetist in the operating room.

#### 1.4.6 Environmental conditions

The device is intended for use in rooms with a controlled environment.

The device must not be exposed to adverse/harsh environmental conditions.

### 1.5 Warnings/caution

The device is intended for use in rooms with a controlled environment.

The device must not be exposed to adverse/harsh environmental conditions.

If there is a suspicion that liquid or condensation has penetrated the handle, it must not be charged under any circumstances.

Never place the handles in liquids.

Make sure that no liquids penetrate the housing interior. Batteries and rechargeable batteries must never be cleaned with water or aqueous solutions.

At exposed places that remain illuminated, bulbs can generate enough heat to burn human tissue.

Use only cleaned/disinfected/sterilised laryngoscopes to limit the risk of cross-contamination.

Remove the batteries and bulbs before cleaning, disinfecting, or sterilising the laryngoscope system.

Please do not allow excess fluid to seep into the electronic contacts. The batteries must be removed before cleaning and sterilisation.

The sterilization guidelines listed by Rudolf Riester GmbH are set out as procedures with correspondingly compatible specific materials. The sterilization must be carried out according to the approved hospital protocol. Rudolf Riester GmbH cannot guarantee sterility. This is to be validated by the hospital or the manufacturer of the sterilization equipment.

Ultrasonic cleaning is strictly prohibited.

Do not immerse in bleach, Betadine or potassium hydroxide solutions. Doing so will severely damage instruments. Avoid metal to metal contact during and after soaking.

The batteries must be inserted in the correct position.

Attention!

It is your responsibility to:  
Before each use, check the status of the internal electrical power source by powering on the device.

Make sure that an open battery handle cover is used when inserting the ri-accu®/ri-accu® L into the rechargeable battery handle.

With certain rechargeable cells, it is possible for the voltage to drop rapidly during use and this can lead to a rapid decline in illumination.

New batteries should be inserted when the light intensity of the instrument becomes weaker and could thereby adversely affect its application

For optimal light output, we recommend that you always use new high-quality batteries when changing the battery (as described in 3.2.).

Note that for the warm light or cold light blades, the handles intended for each respectively should be used. (Cold light handles are marked with a green ring.)

For disposable blades (ri-dispo), please use our handles for the fibre optic laryngoscopes.

If you do not use the device for a long time or take it with you while travelling, please remove the batteries and rechargeable batteries from the handle.

Cleaning frequency and procedure must comply with the cleaning regulations of non-sterile products in the respective facility.

The product may only be used by trained personnel.

The function of the laryngoscope blades and handles should always be tested after cleaning/disinfection/sterilization and before use. To test, place the laryngoscope blade on the handle and click it to the „ON“ position. If the light on the blade does not light up, or flickers, please check the bulb/batteries and the electrical contacts. Make sure that there is always an adequate supply of spare bulbs and spare batteries at hand. If the problem persists, please contact the manufacturer.

All serious incidents related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is resident.


## 2. First use

### 2.1 Scope of delivery

Item no.: 12230 to 12225	1 laryngoscope blade 1 user manual
Item no.: 8040 to 8110	1 set (case, 1 handle, 3 blades) 1 user manual
Item no.: 8120 to 8158	1 set (case, 1 handle, 1 battery, 3 blades) 1 Instructions for use of laryngoscopes

Item no.: 7040	1 ri-charger with accessories 1 ri-charger instructions for use 1 set (case, 1 handle, 3 blades) 1 user manual
Item no.: 7041 to 7042	1 set (case, 1 handle, 1 battery, 3 blades) 1 Instructions for use of laryngoscopes 1 ri-charger with accessories 1 ri-charger instructions for use
Item no.: 7050	1 set (case, 1 handle, 3 blades) 1 user manual
Item no.: 7051 to 7052	1 set (case, 1 handle, 1 battery, 3 blades) 1 Instructions for use of laryngoscopes 1 ri-charger with accessories 1 ri-charger instructions for use
Item no.: 7060	1 set (case, 1 handle, 3 blades) 1 user manual
Item no.: 7061 to 7062	1 set (case, 1 handle, 1 battery, 3 blades) 1 Instructions for use of laryngoscopes 1 ri-charger with accessories 1 ri-charger instructions for use
Item no.: 7070	1 set (case, 1 handle, 3 blades) 1 user manual
Item no.: 7071 to 7072	1 set (case, 1 handle, 1 battery, 3 blades) 1 Instructions for use of laryngoscopes 1 ri-charger with accessories 1 ri-charger instructions for use

## 2.2 Device function

1. Laryngoscope blade
2. Handle

Note that for the warm light or cold light blades, the handles intended for each respectively should be used. [Cold light handles are marked with a green ring.]

## 3. Operation and function

### 3.1 Symbol identification

Green dot

Green ring

Only combine laryngoscopes with a green dot with handles with a green ring.

### 3.2 Startup

#### 3.2.1 Battery handles product range

a. Battery handle type C 28 mm diameter, 2.5 V for standard laryngoscopes, fibre optic laryngoscopes and ri-dispo. To operate these battery handles, you need 2 standard alkaline batteries type C 1.5 V [IEC standard designation LR14] or a rechargeable battery 2.5 V or 3.5 V [ri-accu®/ri-accu® L]. The handle with the Riester battery [ri-accu®/ri-accu® L] can only be charged in the Riester ri-charger® charger. Important:

The handle with Riester ri-accu® L can only be charged in the Riester ri-charger® L.

This is different if you have purchased a ri-accu® L (10694) with plug-in charger from Riester. This can only be charged with our new plug-in charger (10707,10708)

b. Battery handle type AA Ø 19 mm, 2.5 V for standard laryngoscopes, fibre optic laryngoscopes and ri-dispo.

To operate this battery handle, you need 2 standard batteries type AA Mignon 1.5 V [IEC standard designation LR6] or a rechargeable battery [ri-accu®/ri-accu® L] 2.5 V.

Important:

The handle with Riester ri-accu® L can only be charged in the Riester ri-charger® L

c. Rechargeable battery handle type C 2.5 V or 3.5 V for charging in a 230 V or 120 V socket

Fibre optic laryngoscopes and ri-dispo.

The handle is available as a 2.5 V or 3.5 V version and can be ordered for 230 V or 120 V operation. Please note that the handle can only be operated with the ri-accu®/ri-accu® L from Riester.

When using the new ri-accu® L 10692, make sure that no insulation is attached to the spring of the battery handle cover. When using the old ri-accu® L 10692 insulation must be attached to the spring (risk of short circuit).

new ri-accu® L

old ri-accu® L

### 3.2.2 Notes on electromagnetic compatibility

There are currently no indications that electromagnetic interactions with other devices can occur when the devices are used as intended. Nevertheless, under the increased influence of unfavourable field strength, e.g. when operating radio telephones and radiological instruments, interference cannot be completely ruled out.

## 3.3 Inserting and removing batteries and rechargeable batteries

### 3.3.1 Handle types 3.2 a and b

Unscrew the handle cover on the lower part of the handle.

Insert either 2 batteries [as described in 3.2.] or a rechargeable battery from Riester [ri-accu®/ri-accu® L] 2.5 V or 3.5 V into the handle casing so that the positive poles point towards the upper part of the handle. On the rechargeable battery [ri-accu®/ri-accu® L] you will also find an arrow next to the plus sign, to indicate the direction of insertion into the handle. Screw the handle cover firmly onto the handle. Remove the batteries or the ri-accu®/ri-accu® L rechargeable battery by first unscrewing the battery handle cover. The handle must be shaken downwards a little so that the batteries or the rechargeable battery can be removed. Before using the rechargeable battery handles, the ri-accu®/ri-accu® L battery must be inserted into the handle and charged in the Riester ri-charger®/ri-charger® L charger. The ri-charger® charger comes with an additional user manual that must be observed.

Important:

The handle with Riester ri-accu® L can only be charged in the Riester ri-charger® L.

The batteries must be inserted in the correct position.

### 3.3.2 Handle type 3.2 C.

Before using the socket handle for the first time, it should be charged in the socket for a maximum of 24 hours.

Please note whether you are using a new ri-accu® L or an old one, see pictures above for a description. The socket handle [only with NiMH batteries] must not be charged for longer than 24 hours.

Unscrew the handle cover on the lower part of the handle. Depending on which handle you have purchased for which voltage [see 3.2], insert the respective rechargeable batteries into the handle sleeve. With the 2.5 V batteries, make sure that you insert the battery into the handle with the plus side in the direction of the upper part of the handle; next to the plus sign you will also find an arrow that shows you the direction in which to insert it into the handle. With old NiMH 3.5 V rechargeable batteries, it doesn't matter in which direction you use them. Please note whether you are using a new ri-accu® L or an old one, see pictures above for a description.

Screw the handle cover firmly back onto the handle.

Turn the lower part of the handle counter-clockwise. The plug contacts become visible. Round contacts are for 230 V mains operation, flat contacts are for 120 V mains operation. Now plug the handle base into the socket for charging.


The handle must never be in the socket when replacing the battery!

Technical data: Optionally 230 V or 120 V

### 3.4 Attachment of laryngoscope blades

#### 3.4.1 On/off function

Attach the desired laryngoscope blade onto the latch on the upper part of the handle and pull the front part of the blade upwards until the back of the blade is on the same line as the upper part of the handle. A clicking sound can be heard. When the blade clicks into place, the lamp is switched on. The lamp switches off as soon as the blade is disengaged from its latch.

Note that for the warm light or cold light blades, the handles intended for each respectively should be used. [Cold light handles are marked with a green ring.]

At exposed places that remain illuminated, bulbs can generate enough heat to burn human tissue.

#### 3.4.2 Replacing the light bulb

Warm light blade [vacuum lighting]:

Unscrew the bulb from the front part of the blade and screw in a new one.

Cold light blade [Xenon and LED lighting]:

- Unscrew the upper part of the handle off of the handle casing
- Unscrew the bulb mount by the knurl and pull it out.
- Unscrew the bulb.
- Screw in a new bulb.
- Screw the bulb mount into the upper part of the handle.
- Screw the upper part of the handle into the handle casing.

#### 3.4.3 Removing/inserting the light guide ri-modul

By laterally unlatching/latching the light guide at the hook-in connection of the blade.

### 4. Care instructions

#### 4.1 General information

The cleaning and disinfecting of the medical devices serves to protect the patient, the user and third parties and to maintain the value of the medical devices.

The product design and materials used make it impossible to define an upper limit on max. feasible treatment cycles. The service life of medical devices is determined by their function and careful handling.

Before return for repair, defective products must have undergone the prescribed reprocessing procedure.

Remove the batteries and bulbs before cleaning, disinfecting, or sterilising the laryngoscope system.

If a reusable device shows signs of material deterioration, it should no longer be used and should be disposed of/claimed according to the procedures described in the disposal/warranty sections.

#### 4.2 Cleaning and disinfection

##### 4.2.1 Laryngoscope blades/handles

###### Cleaning procedure

Standard blade: Bulb must be removed before cleaning, disinfection and sterilization.

F.O. handle: Batteries and bulb must be removed before cleaning, disinfection and sterilization.

Standard handle: Batteries must be removed before cleaning, disinfection and sterilisation.

Do not allow excess fluid to seep into electrical contacts.

Do not immerse in bleach, Betadine or potassium hydroxide solutions.

Doing so will severely damage instruments. Avoid metal to metal contact during and after soaking.

Ultrasonic cleaning is strictly prohibited.

#### 4.2.2 Point of use:

No particular requirements

#### 4.2.3 Containment and transport:

It is recommended to transport the contaminated medical devices in a closed container.

It is recommended to reprocess the medical devices immediately after use (2 hours between soiling and pre-cleaning was validated). A prolonged intermediate storage of medical devices used with impurities such as blood residues can lead to corrosion damage.

#### 4.2.4 Preparation for decontamination:

##### 4.2.5 Blades:

Personal protective equipment (gloves, water repellent protective skirt, face protection mask or protective glasses)

1) Disassemble the medical device (not ri-integral) (see Figure 1 example).


Figure 1

##### 4.2.6 Handle

Personal protective equipment (gloves, water repellent protective skirt, face protection mask or protective glasses)

1) Disassemble the medical device (see Figure 1 example).

2) Unscrew the cover of the batteries.

3) Remove the batteries.

4) Unscrew and remove the bulb from the head (not standard version).


Fig. 2

##### 4.2.7 Pre-cleaning

The laryngoscope system should be rinsed under cool running tap water (minimum drinking water quality), until all visible soil is removed (1 minute was validated).

Rinse the hard-to-reach areas with a cleaning gun for at least 1 minute (1 minute was validated).

**Note:** If the connector and spring mechanism on the bulb are contaminated, rinse this area with a cleaning gun for at least 1 minute (1 minute was validated).

##### 4.2.8 Cleaning: manually

- **Equipment:** enzymatic cleaner (neodisher MediZym, Dr. Weigert # 404033 was validated), cleaning brush (Interlock # 09098 was validated), tap water/flowing water ( $20 \pm 2^\circ\text{C}$  /  $68^\circ\text{F} \pm 35,6$ ) (minimum drinking water quality), tank/basin for cleaner.

**Note:** the bulb may be cleaned and disinfected with a cotton ball dampened with alcohol

1. The cleaner/cleaning solution is to be prepared in accordance with the manufacturer's instructions for the detergent (neodisher Medizym; 0.5% and  $40^\circ\text{C}$  /  $104^\circ\text{F}$  was validated).

2. Soak the medical device completely in the cleaning solution.

3. Brush the hard-to-reach areas of the medical device with a softbrush for at least 1 minute (1 minute was validated). Pay attention to the critical hard-to-reach areas where it is not possible to assess the cleaning efficacy.

4. Handle: Move the spring mechanism 10 times during the exposure time.
5. Make sure that all surfaces are completely moistened with cleaning solution.
6. Total exposure time in the cleaning solution in accordance with the manufacturer`s instructions (20 minutes was validated).
7. Remove the medical device from the cleaning solution.
8. Rinse the medical device completely under running tap water (at least drinking water quality) for minimum 1 minute (1 minute was validated) to remove the cleaning solution completely. Check for cleanliness; if dirt is visible, repeat the steps above.

#### 4.2.9 Disinfection: manually

**- Equipment:** Non-protein-fixing VAH-listed instrument disinfectant based on aldehydes (Johnson & Johnson, CIDEX OPA solution # 20391) was validated [active ingredient: ortho-phthalaldehyde 0.55%/100 g) (20°C / 68°F) demineralised water [demineralised water, according to KRINKO/BfArM recommendation free of facultative pathogenic microorganisms) (20 ± 2°C / 68°F ± 35,6), disinfectant tank, tank for demineralised water, lint-free gauze, medical compressed air (according to the European Pharmacopoeia).

1. Pour the disinfectant solution into a disinfectant tank (CIDEX OPA is a ready-to-use solution whose concentration must be confirmed using test strips, see manufacturer`s instructions).
2. Immerse the medical device completely in the disinfectant solution.
3. Make sure that all surfaces are completely moistened with disinfectant solution.
4. Brush the hard-to-reach areas during the exposure time for at least 1 minute (1 minute was validated).
5. Exposure time of the disinfectant solution in accordance with the manufacturer`s instructions (CIDEX OPA, for high-level disinfection 12 minutes was validated).
6. Remove the products from the disinfectant solution.
7. Place the instruments in a tank of demineralised water for at least 1 minute (1 minute was validated).
8. Repeat step 5 two times with fresh demineralised water to remove the disinfectant solution completely.
9. Wipe with a lint-free gauze and/or dry with medical compressed air.

#### 4.2.10 Cleaning and disinfection: Automated WD

**- Equipment:** Washer disinfecter in accordance with EN ISO 15883-1 + -2 with thermal program (temperature 90-95°C / 194- 203°F), enzymatic cleaner (neodisher MediZym, Dr. Weigert # 404033 was validated), lint-free gauze, medical compressed air (as per the European Pharmacopoeia).

1. Place the medical device in a suitable tray or place on the load carrier so that all internal and external surfaces are cleaned and disinfected.
2. Close the WD and start the program; the program parameters are shown in the table below.

Program step	Water	Dosage	Time	Temperatures
Pre-rinse	Cold		5 min	
Dosage cleaner		According to manufacturer's instructions (0.2% was validated)		According to manufacturer's instructions
Cleaning	Deionised water		According to manufacturer's instructions (5 min was validated)	[40°C was validated]
Rinse	Deionised water		2 min	
Disinfection	Deionised water		A0 value >30001 (e.g. 5 min, 90°C)	
Drying			15 min	Up to 120°C
<sup>1</sup> Public authorities may adopt other implementing provisions (parameters for disinfection performance) within their field of responsibility				

3. At the end of the program, remove the medical device.
4. Check for dryness; if necessary, wipe with lint-free gauze and/or dry with medical compressed air.
5. After removal from the WD check the products for cleanliness. If dirt is still visible, clean the products manually. Afterwards the automated cleaning procedure must be carried out again.

#### 4.2.11 Maintenance, inspection and testing:

1. Assemble the medical device
2. All medical devices must be visually checked for cleanliness, dryness and damage (e.g. cracks, fractures, corrosion, mobility, bending etc.), if necessary, using an illuminated magnifier [3-6 Dptr.]
3. Damaged medical devices must be separated out and not be used again.

#### 4.2.12 Packaging:

- **Equipment:** Sterilisation fleece (e.g. steriCLIN # 3FVLI330114 was validated), sealing machine (e.g. Hawo, Hm 880 DC-V)

1) A suitable method (sterile barrier system) must be used for packaging the products. Packaging according to DIN EN ISO 11607 (single) or DIN 58953-9 (sets)

- **Single:** Use a sterile barrier system (e.g. see-through sterilisation pouches or sterilisation fleece) according to DIN EN ISO 11607, which is intended by the manufacturer for steam sterilisation. The products can be packed single or double. The packaging must be large enough so that the seal is not under tension (for validation it was double wrapped in sterilisation fleece).

- **Sets:** Use a sterile barrier system (sterilisation fleece) according to DIN EN ISO 11607, which is intended by the manufacturer for steam sterilisation. The medical devices can be packed single or double in sterilisation fleece or in a sterile container according to DIN 58953-9. The weight of the sterile container must not exceed 10 kg.

- **Note:** After the sealing process, the seal must be visually inspected for possible faults. If the packaging is damaged the device must be repackaged and resealed.

#### 4.2.13 Sterilisation:

- **Equipment:** Steam steriliser according to DIN EN 285 or DIN EN 13060 with type B process

Note: Do not exceed temperature of 135°C / 275°F and pressure of 28 psi. Pre-vacuum process, 134°C / 273,2°F and sterilisation time at least 3 min or 132°C / 269,6°F at least 3 min (longer holding times are possible) [132°C / 269,6°F, 90 seconds was validated].

1. Place the packaged medical device in the steriliser chamber.
2. Start the sterilisation program.
3. At the end of the sterilisation program remove the device and let it cool down.
4. Check the packaging for damage or moisture penetration. If a package fails this check it must be considered non-sterile. The products must be repackaged and sterilised again.

#### 4.2.14 Storage:

Storage (at least protected against dust and moisture) and storage time according to user specifications.

#### 4.2.15 Additional information

The processing of medical devices should be done with validated processes.

#### 4.3 Test procedure for blade and handle

The function of the laryngoscope blades and handles should always be tested after cleaning/disinfection/sterilization and before use. To test, place the laryngoscope blade on the handle and click it to the „ON“ position. If the light on the blade does not light up, or flickers, please check the bulb/batteries and the electrical contacts. Make sure that there is always an adequate supply of spare bulbs and spare batteries at hand. If the problem persists, please contact the manufacturer.

The sterilization guidelines listed by Rudolf Riester GmbH are set out as procedures with correspondingly compatible specific materials. The sterilization must be carried out according to the approved hospital protocol. Rudolf Riester GmbH cannot guarantee sterility. This is to be validated by the hospital or the manufacturer of the sterilisation equipment.

#### 5. Technical specifications

Operating conditions: 0°C / 32°F to 40°C / 104°F with a relative humidity of 15% to 90% (non-condensing)

Storage conditions: -20°C / -4°F to 50°C / 122°F with a relative humidity of 15% to 85% (non-condensing)

Air pressure: 620 – 1060 hPa

#### 5.1 Technical data for the bulb

Bulb, small, warm light	2.7 V 300 mA average lifespan 20 h
Bulb, large, warm light	2.7 V 300 mA average lifespan 20 h
Bulb, cold light XL	2.5 V 680 mA average lifespan 20 h
Bulb, cold light XL	3.5 V 700 mA average lifespan 20 h
Bulb, cold light LED	2.5 V 120 mA average lifespan 20,000 h
Bulb, cold light LED	3.5 V 280 mA average lifespan 20,000 h

## 6. Spare parts and accessories

Item no.:

- 12320 Bulb, LED 2.5 V for F.O. laryngoscopes + ri-dispo, e-xam
- 12321 Bulb, LED 3.5 V for F.O. laryngoscopes + ri-dispo, e-xam
- 11428 Pack of 6 bulbs XL 2.5 V, F.O. laryngoscopes
- 11429 Pack of 6 bulbs XL 3.5 V, F.O. laryngoscopes
- 12320 Bulb, LED 2.5 V for F.O. laryngoscopes + ri-dispo, e-xam
- 12321 Bulb, LED 3.5 V for F.O. laryngoscopes + ri-dispo, e-xam
- 11428 Pack of 6 bulbs XL 2.5 V, F.O. laryngoscopes
- 11429 Pack of 6 bulbs XL 3.5 V, F.O. laryngoscopes
- 11380 Pack of 6 bulbs 2.7 V, small, without spring, for blade ri-standard Miller no. 00 to 1
- 11381 Pack of 6 bulbs 2.7 V, large, without spring, for blade ri-standard Miller no. 2-4, McIntosh No. 0-
- 12314 Battery handle type C, rechargeable, LED 3.5 V, for blades with F.O.
- 12316 Battery handle type AA rechargeable, LED 3.5 V, for blades with F.O.
- 12318 Battery handle type AA rechargeable, XL 3.5 V, for blades with F.O.
- 12300 Battery handle type C rechargeable, XL 2.5 V, for blades with F.O.
- 12301 Battery handle type C rechargeable, XL 3.5 V, for blades with F.O.
- 12302 Battery handle type AA rechargeable, XL 2.5 V, for blades with F.O.
- 12311 Battery handle type C, non-rechargeable, LED 2.5 V, for blades with F.O.
- 12312 Battery handle type AA, non-rechargeable, LED 2.5 V, for blades with F.O.
- 12317 Battery handle type AA short, non-rechargeable, LED 2.5 V, for blades with F.O.
- 12305 Battery handle type C non-rechargeable, XL 2.5 V, for blades with F.O.
- 12306 Battery handle type AA non-rechargeable, XL 2.5 V, for blades with F.O.
- 12309 Battery handle type AA short, non-rechargeable, XL 2.5 V, for blades with F.O.
- 12303 Battery handle type C, rechargeable for blades without F.O.
- 12304 Battery handle type AA rechargeable, for blades without F.O.
- 12307 Battery handle type C, non-rechargeable, for blades without F.O.
- 12308 Battery handle type AA non-rechargeable, for blades without F.O.
- 12310 Battery handle type AA non-rechargeable, short, for blades without F.O.
- 12265 Light guide for blade 12260 Macintosh No. 0
- 12266 Light guide for blade 12261 Macintosh No. 1
- 12267 Light guide for blade 12262 Macintosh No. 2
- 12268 Light guide for blade 12263 Macintosh No. 3
- 12269 Light guide for blade 12264 Macintosh No. 4
- 12273 Light guide for blade 12272 Macintosh No. 5
- 12280 Light guide for blade 12274 Miller no. 00
- 12281 Light guide for blade 12275 Miller no. 0
- 12282 Light guide for blade 12276 Miller no. 1
- 12283 Light guide for blade 12277 Miller no. 2
- 12284 Light guide for blade 12278 Miller no. 3
- 12285 Light guide for blade 12279 Miller no. 4
- 10694 Li-ion rechargeable battery, 3.5 V, ri-accu®L for plug-in charger for battery handle type C
- 10680 ri-accu® 2.5 V NiMH for battery handle type AA
- 10681 ri-accu® 2.5 V NiMH for battery handles type C
- 10682 ri-accu® 3.5 V NiMH for battery handles type C
- 10690 Li-ion battery, 3.5 V, ri-accu®L for battery handle type AA and ri-charger® L
- 10691 Li-ion battery, 3.5 V, ri-accu®L for battery handle type C and ri-charger® L
- 10685 4 alkaline AA batteries (Mignon), type E 91, alkaline, 1.5 V
- 10686 2 alkaline C batteries, 1.5 V, MV 1400, LR14

### Charging station and accessories

- 10707 Plug-in charger for Li-ion battery, 3.5 V, ri-accu®L for battery handle type C, with EU plug
- 10708 Plug-in charger with Li-ion battery, 3.5 V, ri-accu®L for battery handle type C, with EU plug
- 10700 ri-charger 2.5 V / 3.5 V / 230 V with wall mount
- 10701 ri-charger® 2.5 V / 3.5 V / 120 V with wall mount
- 10705 ri-charger® L 3.5 V/230 V with wall mount
- 10706 ri-charger® L 3.5 V/120 V with wall mount
- 10710 ri-charger® wall mount with screws and dowels
- 10711 ri-charger® AA type handle insert, black
- 10712 ri-charger® cover for ear speculum compartment
- 11389 Pouch 24 x 20 cm black, laryngoscopes
- 11393 Hard case for up to 12 blades and 1 handle
- 11385 Case for laryngoscope sets ri-modul Miller
- 11386 Case for laryngoscope sets ri-modul Miller baby
- 11392 Case for laryngoscope sets ri-integral, ri-modul and ri-standard Macintosh baby
- 11391 Case for laryngoscope sets ri-integral, ri-modul and ri-standard Macintosh

11387 Case for laryngoscope sets ri-integral and ri-standard Miller  
11388 Case for laryngoscope sets ri-integral and ri-standard Miller baby  
99204 GA, laryngoscopes

## 7. Maintenance/accuracy check/calibration/applied standards

The instruments and their accessories require no special maintenance. If an instrument needs to be tested for any reason, please send it to us or an authorized Riester dealer in your area, the details of which we will provide you with upon request.

### 7.1 Applied standards

Applied standard: ISO 7376

Blades and handles from manufacturers that meet these standards are compatible.

The following combinations were tested:

F.O. Riester blade (ri-integral, ri-modul) with handles from the following manufacturers;

- Heine F.O. handle
- Rusch F.O. handle
- Vital Signs F.O. handle
- Riester F.O. handle
- Penlon F.O. handle
- Timesco F.O. handle
- Kawe F.O. handle
- MD F.O. handle
- Welch Allyn F.O. handle

F.O. Riester handles with blades from the following manufacturers;

- Riester F.O. blade
- MD Maxlite blade
- Heine F.O. blade
- MD F.O. blade

## 8. Disposal

### Caution!

The used medical device must be disposed of in accordance with current medical practices or local regulations on the disposal of infectious biological medical waste.

Batteries and electrical/electronic devices may not be treated as domestic waste and must be disposed of in accordance with local regulations.

If you have any questions about the disposal of products, please contact the manufacturer or their representative.

## 9. Warranty

This product has been manufactured under the strictest quality standards and has undergone a thorough final quality check before leaving our factory. We are therefore pleased to be able to provide a warranty of 2 years from the date of purchase on all defects, which can verifiably be shown to be due to material or manufacturing faults. A warranty claim does not apply in the case of improper handling. All defective parts of the product will be replaced or repaired free of charge within the warranty period. This does not apply to wearing parts. For r1 shock-proof, we grant an additional warranty of 5 years for the calibration, which is required by CE-certification. A warranty claim can only be granted if this Warranty Card has been completed and stamped by the dealer and is enclosed with the product. Please remember that all warranty claims have to be made during the warranty period. We will, of course, be pleased to carry out checks or repairs after expiry of the warranty period at a charge. You are also welcome to request a provisional cost estimate from us free of charge. In case of a warranty claim or repair, please return the Riester product with the completed Warranty Card to the following address:

**Rudolf Riester GmbH**  
**Dept. Repairs RR**  
**Bruckstr. 31**  
**72417 Jungingen**  
**Germany**

**Serial number or batch number, date,**  
**stamp and signature of the specialist dealer**

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