# MANUAL DE INSTRUCCIONES

## Silla de ruedas BOLONIA

Manual de instrucciones	ES
Instruction manual	EN
Manuel D'Utilisation	FR
Manuale D'Istruzioni	IT
Anweisungen	DE
Instruções	PT
Gebruiksaanwijzing	NL
Bruksanvisning	SWE
Instrukcja obsługi	PL
Brugsanvisning	DK

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**NOTICE**

Notice to the user and/or patient: any serious incident which has occurred in connection with the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.


**INDICATIONS**

The wheelchair is used as a mobility tool for the handicapped, patients and the elderly who have difficulties in elderly people who have difficulties in action.

Beginners should be accompanied by assistants to avoid accidents due to improper handling.

This wheelchair is designed for frequent or occasional use by a single occupant weighing up to 110 kg and can be used indoors and outdoors. The wheelchair can be manually moved forward and backward. All Mobiclinic wheelchairs are manufactured and CE compliant.


**FOLDING AND UNFOLDING THE WHEELCHAIR:** To fold the wheelchair, just lift the center edge of the seat upholstery from the bottom. To unfold, tilt the wheelchair slightly to one side to lift the wheels off the floor. Then press down on one or both seat rails. Before folding the chair, completely remove the footrests or make sure they are fixed and in an upright position.



**LOCKING THE WHEELS:** Push the wheel locks forward until they engage in the locked position. Do not attempt to sit down or stand up without having activated the brakes. Always apply both brakes simultaneously.

**ADJUSTING THE LENGTH OF THE FOOTREST:** The position of the footrest (standard-flat or elevated) can be adjusted to suit the length of the user's legs. The footrest should be adjusted to support the weight of the user's feet and lower legs in a position to support the weight of the thighs. If the footrest adjustment is too long, the thighs will be pressed under the front of the seat upholstery. If the footrest adjustment is too short, the user's knees will rise and the buttocks will have to support excessive weight. Both situations lead to an increased risk of pressure ulcers. The user's knees and hips should be at approximately the same level. The lower edge of the footrest should be at least 7 cm above the floor to allow sufficient clearance when going over ramps or slopes. To do this, loosen the adjustment screw with a wrench and fold the foot plates inward or outward to achieve the desired position. Then tighten the screw carefully.


**WARNINGS:**

1. Make sure the brakes are applied before getting in or out of the wheelchair.
2. Do not stand on the footrests when getting out of the chair, otherwise the chair may tip over.

3. Keep the wheelchair away from fire.

**USING THE REAR SUPPORT PEDALS AND BRAKES:** Use one of the rear support pedals to lift the front wheels, for example, when climbing a curb. To do this, push the rear support pedals down with one foot. Do not lift the small front wheels by pushing down on the handles, as this could damage the wheelchair.

**CLIMBING A CURB:** Approach the curb from the front. The caregiver should use the rear support pedals to raise the front wheels and lower them onto the curb once the curb is raised. Finally, the caregiver should push the wheelchair forward, lifting it slightly to raise the curb if necessary.

**LOWERING A CURB:** Place the front wheels in front of the curb. The caregiver should use the support pedals to raise the front wheels and tilt the occupant back slightly. While keeping the wheels up, slightly lower the wheelchair toward the curb.

**Warning:** When going down a curb, the front wheels should be raised to prevent the occupant from falling off. These are recommendations only; they may not be suitable for all cases.

### IMPORTANT ASPECTS TO TAKE INTO ACCOUNT

The user or caregiver should check the safety of the wheelchair at regular intervals.

1. Check the handgrips and non-skid tires and make sure that they are solid and secure.
2. For a correct fit, check the locking devices and make sure that the large wheels are locked when they are activated.
3. Look at all bolts, nuts and joints to check that they remain solid.
4. If the wheelchair is equipped with pneumatic tires, check that the tire pressure is adequate.
5. Check that the adjusted length of the footrest is appropriate.
6. If other conditions affecting safety are present, stop using the wheelchair.
7. Remember to engage the wheelchair lock before transferring to or from the wheelchair.
8. Do not lean forward out of the wheelchair unless both feet are on the floor. This is especially important for users with heavy casts with elevated foot rests.

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FOSHAN KAIYANG MEDICAL TECHNOLOGY GROUP CO., LTD.
SUNGO EUROPE B.V.
GRUPO R. QUERALTÓ S.A.