MANUAL DE INSTRUCCIONES

DOPPLER

Detector fetal 2Mhz

Detector vascular 8Mhz

Last revision date: 21 de marzo de 2023
Última revisión: March 21st, 2023

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
Fetal detector 2Mhz

Vascular detector 8Mhz

Please read instructions before use

2.782.116.01(QUE)TSS/1.0 1.4.01.34.439 2019.02

NOTICE
Notice to the user and/or patient: any serious incident which has occurred in connection with the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

ATTENTION
This user manual is written and compiled in accordance with Council Directive MDD93/42/EEC for medical devices and harmonised standards. In the case of software modifications and updates, the information contained in this document is subject to change without notice. The manufacturer makes no warranties of any kind with respect to this material, including but not limited to the implied warranties of merchantability or fitness for a particular purpose. The manufacturer assumes no liability for any errors that may appear in this document, for incidental or consequential damages in connection with furnishing, performance or use of this material. No part of this document may be photocopied, reproduced or translated into another language without the manufacturer's prior written consent. The information contained in this document is subject to change without notice.

MANUFACTURER'S LIABILITY
The manufacturer shall only be held liable for effects on safety, reliability and performance of the equipment if: Assembly operations and repairs are carried out by persons authorised by the manufacturer and the device is used in accordance with its instructions for use.

WARNING
This device is not intended for treatment. The intended use is for fetal heart rate detection. If the result is not credible, please use other methods such as stethoscopy to verify immediately.

1.SAFETY GUIDE
This unit is internally equipped with power; the degree of discharge protection is type B applied part. Type B applied part means that these patient connections comply with the permissible leakage current, dielectric strength of IEC 60601-1.

- 1.1 Safety precautions**
- The warning messages must be observed. To avoid the possibility of injury, observe the following precautions during operation of the device.
 - The device is not explosion proof and should not be used in the presence of flammable anaesthetics.
 - Do not dispose of batteries in fire as they may explode.
 - Do not attempt to recharge normal dry cell batteries, they may leak and cause a fire or explosion.
 - Do not touch the signal input or output connector and the patient simultaneously.
 - Auxiliary equipment connected to the analogue and digital interfaces shall be certified according to the respective IEC standards (e.g. IEC 950 for data processing equipment and IEC60601-1 for medical equipment). In addition, all configurations shall comply with the valid version of the IEC60601-1-1 system standard. If you connect the additional equipment to the signal input or output connector, you configure a medical system and are therefore responsible that the system complies with the requirements of the valid version of the IEC 60601-1-1 standard. If in doubt, please consult our service department or your local distributor.
 - The pocket fetal detector is a tool to assist in professional healthcare and should not be used in place of a normal fetal monitor.
 - Please use the probe provided by the manufacturer. Use of other accessories may damage the device and affect the performance and safety of the device.
 - Do not pull the probe cable more than 1 metre, otherwise the probe may separate from the detector connector.
 - If in doubt about the measurements, consult a physician.
 - Keep the equipment out of the reach of children.
 - Do not make any modifications to the device.
 - The device should only be serviced by authorised personnel. The device must not be serviced or repaired while in use.
 - The battery can be replaced by the patient.
 - Use conductive gel. In case of allergies, contact your physician.
 - Keep the equipment clean. Avoid vibrations.
 - Do not use high temperature sterilisation processes, gamma radiation or electron beam irradiation.
 - Electromagnetic Interference: Make sure that the environment where the device is to be used is not subject to sources of electromagnetic interference sources, such as radio transmitters, mobile phones, etc. Keep them away.
 - The user should check the equipment for visible evidence of damage that may affect patient safety or monitoring capability prior to use. The recommended inspection interval is once per month or less. If damage is evident, replacement before use is recommended.
 - The following safety checks should be performed once every two years or as specified in the institution's test and inspection protocol by a qualified person.
 - Inspect the equipment for mechanical and functional damage.
 - Check relevant safety labels for legibility.
 - Verify that the device is functioning properly as indicated in the instructions for use.
 - Approve the patient leakage current according to IEC 60601-1/1998: Limit μA (B). The leakage current must not exceed the limit. The data must be recorded in an equipment log. If the device does not work properly or any of the above tests are missing, the device must be repaired.
 - The battery must be disposed of in accordance with local regulations after use.
 - The battery must be removed from its compartment if the device is not to be used for an extended period of time.
 - The device may only be used if the battery cover is closed.
 - The battery should be stored in a cool, dry place.
 - If using a rechargeable battery, to ensure capacity and life, please fully charge the batteries before the first use, usually the batteries should be continuously charged for more than 14 hours or charges according to the guide on the battery.
 - Be sure to install the batteries correctly.
 - The valid period of this product is five years.
 - After the service life, please dispose of the products according to local regulations.
 - The device should not be used with defibrillator or high frequency surgical unit.
 - Please choose the accessories authorized by our company or the device may be damaged.
 - Please keep the probe away from sharp objects.
 - Do not use strong solvents such as acetone.
 - Never use abrasive products such as aluminium scouring pad or metal polish.
 - Do not allow any liquid to enter the product and do not immerse any part of the device in any liquid.


2. INTRODUCTION

2.1. General information
The Pocket Fetal Detector is a hand-held obstetric unit that can be used in the hospital, in clinics or at home for daily self-monitoring. It contains ultrasonic signal transmitting and receiving components, analogue signal processing unit, FHR calculation unit, LCD display control unit, etc. It has three modes of operation: FHR, FHR calculation unit and FHR control unit. It has three modes of operation: real-time fetal heart rate display mode, averaged fetal heart rate display mode and manual mode. It also has an audio output that can be connected to headphones or recorders with audio input.


2.2 Features

- Aesthetic design, portable and easy to use.
- The probe has a flexible structure, which facilitates its use and provides comfort to pregnant women during the operation of the product.
- Battery status indicator.
- Built-in speaker.
- Headphone output
- Interchangeable probe
- Probe inspection
- Automatic shutdown

3. OVERVIEW AND CONFIGURATION

3.1. Display
The LCD display is as follows:

WARNING: Do not attempt to connect probes not authorised by the manufacturer.

3.2. Buttons

There are three buttons (POWER, MODE and START/STOP) and a volume control button on the pocket fetal detector. The main functions are as follows:

3.2.1 Power button

Function: On/Off
On: Press the button once.
Off: Press the button and hold for 3 seconds to turn off.

3.2.2 Mode button

Function: Mode selection, press once to enter the next mode of operation while in use. As the fetal detector has a memory function, when powering on the device will enter the selected mode before powering off automatically.

3.2.3 Start/stop button

Function:
1. In mode 1 and 2, press the button to turn on/off the backlight.
2. In mode 3, the button is to start/stop the operation.

3.2.4 Volume Control Indicator


Volume adjustable direction indicator. From left to right, it means that the sound level is from high to low.

3.3 Headphone Input

Headphone socket: A socket for audio output and can be connected to headphones or recorder with audio input for recording.

3.4 Probe

The meaning of CD--IP-- on the label is as follows:

- C: The operating mode of the probe is a continuous wave.
- D: The shape of the structure for the probe is cell type.
- 2.0/8.0: The frequency of the probe is 2 MHz or 3 MHz.
- IPX4: Degree of protection against liquids. It is protected from splashing water.

4. GENERAL OPERATION

4.1 Fetal heart rate monitoring

- Switch on by pressing the power button. It will perform the automatic test when the machine is turned on. After the automatic test, the LCD display is like the figure on the front panel.
- Finding the position of the fetus. At first, feel the position of the fetus by hand. Find the best direction to inspect the fetal heart. Apply a generous amount of gel on the probe faceplate, place the probe faceplate in the best position to detect the fetal heart. Adjust the probe for optimal audio signal by angling the probe around. Adjust the volume according to requirements.
- Calculation of the fetal heart rate: The result will be displayed on the LCD screen

CAUTION

- Place the probe on the best detection position to obtain better detection effect.
- Do not place the probe in the position where it has strong placental blood sound or strong umbilical sound.
- If the pregnant woman adopts the horizontal position and the fetal position is normal, place the probe in the position of the lower midline of the umbilicus to obtain the clearest fetal heartbeat sound.
- Do not measure the fetal heart rate unless you have heard the audible fetal sound.

4.2. Mode selection

4.2.1 Real-time fetal heart rate display mode (mode 1)

At the time of heart rate signal detection, the LCD will display the flashing heart symbol and real-time heart rate display simultaneously.

4.2.2 Fetal heart rate averaging display mode (mode 2)

Used to obtain more stable heart rate readings. In this mode, the fetal heart rate is an average of 8 beats. The LCD display shows the flashing heart symbol when the heart rate is displayed.

4.2.3 Manual Mode (Mode 3)

When entering Mode 3, the system will automatically count the audible beats, the fetal heart rate is displayed in the format "--" and the LCD display flashes the heart symbol. Press the backlight control button to stop calculating. The unit will automatically calculate the derived fetal heart rate as an average over the calculation time and display the result. If you measure the heart rate again, press the backlight control button to start. Press it again, it will stop calculating. This value remains until you repeat the measurement or change the mode.

4.3 Probe operation

4.3.1. Probe inspection

When the probe is moved away from the device, the LCD will flash "--" and the probe frequency indication data will disappear. At this time the probe must be reconnected. After successful reconnection, the LCD will stop flashing and display the probe frequency data.

4.3.2 Replacing the probe


The manufacturer has connected a probe to the device when packing. If users need to replace it with another probe, turn off the device, then remove the probe from the package. Disconnect the probe plug from the probe socket. Then connect the replacement probe plug to the socket.

Note: place the temporarily unused probe carefully and avoid dropping it. When the device is not used for a long time, the user is recommended to connect the plug of a probe to the detector socket and put the probe in the package. Then store the device with the probe in the box.

4.3.3 Removing and inserting the probe

1. Remove the probe

Hold the main unit with one hand and grasp the probe handle with the other hand to pull out the probe.

2. Insert the probe: This is the opposite of removing the probe. Grasp the main unit with one hand and hold the top of the probe with the other hand, then push the probe into the probe holder.

4.4 Low Power Inspection

When operating normally, the LCD displays the battery status and the number in the status box indicates the remaining battery. When the battery is low, the battery displays 0 in the status box to remind the user of the need to replace or charge the battery (rechargeable batteries only).

4.5. Replacing the battery


4.5.1 Removing the battery

On the rear panel, open the battery compartment, remove the batteries.

4.5.2 Battery replacement

Insert two AA size batteries into the battery compartment (for battery orientation, refer to the instructions inside the compartment). Close the battery compartment.

CAUTION: If the device is not to be used for an extended period of time, remove the batteries from the device.

5. PRODUCT SPECIFICATIONS

Product name: Fetal pocket detector

Classification:

- Anti-electrical shock type: internally powered equipment.
- Anti-electrical shock rating: CF type applied part.
- Harmful liquid proof rating: IPX4
- Operating system: Continuous operation equipment.
- EMC: Group 1 Class B
- Suitable range of use: Suitable for use after the 12th week of pregnancy.
- Size: 135 mm (length) x 92 mm (width) x 29 mm (height)
- Weight: Approximately 245 g (including batteries).
- Display: LCD display 38 mm x 28 mm

Environment

1. Operation:
 - Temperature: 5°C~40°C
 - Humidity: ≤80%
 - Atmospheric pressure: 70 kPa~106 kPa
2. Transport and storage
 - Temperature: -10°C~55°C
 - Humidity: ≤93%
 - Atmospheric pressure: 50 kPa~106 kPa

Fetal Heart Rate performance

- Fetal heart rate measurement range: 50~240 BPM (BPM: Beats per minute)
- Resolution: 1 BPM
- Accuracy: ± 2 BPM
- Power consumption: <1 W
- Auto power off: After 1 minute without signal.
- Battery type: 2 x 1.5 V DC batteries (size AA LR6).

Probe

- Normal frequency: 2.0 or 8.0 MHz
- Operating frequency: 2.0 MHz±10%
- P-: <1.5 MPa
- Iob: <20 mW/cm²
- Ispta: <100 mW/cm²
- Ultrasound output intensity: Isata <20 mW
- Mode of operation: Continuous wave Doppler
- Effective area of transducer radiation: ≤157 mm²

6. MAINTENANCE

6.1 Maintenance

This device is breakable and should be handled with care. The gel should be wiped off the probe after use. These precautions will prolong the life of the unit. The user should check the equipment for visible evidence of damage that may affect patient safety or the capability of the pocket fetal detector prior to use. The recommended inspection interval is once per month or less. If damage is evident, replacement prior to use is recommended. The safety of the device should be checked periodically to ensure adequate protection of the patient from current leakage. This should include current leakage measurement. The recommended testing interval is once every two years or as specified by the institution's testing and protocol. Heart rate accuracy is controlled by the device and cannot be adjusted by the user. If the fetal heart rate result is not credible, use other methods such as stethoscopy to verify immediately or contact the distributor or manufacturer for assistance.

6.2 Cleaning

Before cleaning turn off the power and remove the batteries. Keep the outer surface of the device clean and free from dust and dirt and wipe the outer surface (including the display) of the unit with a soft dry cloth. If necessary, clean the cover with a soft cloth soaked in a mixture of soap and water and wipe it dry with a clean cloth immediately. Wipe the probe with a soft cloth to remove any remaining ultrasound coupling gel. Clean with soap and water only.

CAUTION:

- Do not use strong solvents such as acetone.
- Never use abrasives such as aluminum scouring pads or metal polish.
- Do not allow any liquids to enter the product and do not immerse any parts of the device.
- Avoid pouring liquids into the device.
- Do not allow any cleaning mixture to remain on the surface of the device.

NOTE: Dry the surface of the probe with 70% ethanol and allow to dry or wipe with a clean, dry cloth.

6.3. Disinfection

Clean the equipment case, probe, etc. as above and then wipe the probe with a wipe soaked in alcohol (70% ethanol). Wipe the probe with a clean, dry cloth to remove any remaining moisture.

WARNING: Never attempt to sterilize the equipment with low temperature steam or any other method.

7. TROUBLESHOOTING

Problem	Possible reasons	Solutions
Weak sound	1. Volume is too low. 2. Battery is running low. 3. The gel hasn't been applied.	1. Adjust the volume. 2. Replace the batteries. 3. Apply the gel.
Noise	1. The probe is too close to the main unit. 2. There is interference in the output signal. 3. The battery is running low.	1. Put more distance between the probe and the main unit. 2. 2. Keep it away from the signal. 3. Replace the batteries.
Low sensitivity	1. The position of the probe is not correct. 2. The gel hasn't been applied	1. Adjust the position of the probe. 2. Smear the gel.

APPENDIX 1

Relevance of the Internal Fetal Monitor

Modern medicine thinks that fetal heart rate is the essential thing to identify fetal health, recording changes in heart rate can observe fetal hypoxia, fetal distress, umbilical cord around the neck and other symptoms. The fetal home monitor monitors changes in fetal heart rate, listening primarily to the sound of the fetal heart.

The fetal heart rate changes most obviously in the following three periods:

1. 30 minutes after she has stood up.
2. 60 minutes after she has eaten.
3. 30 minutes before going to bed.

In the above three periods, due to the change in the fetal state of the pregnant woman, the activity of the digesting food needs to provide the fetus with more oxygen, relatively the oxygen to the fetus becomes less. It is easy to create symptoms such as fetal anoxia. Checking the heart rate at this time can visualize the healthy state of the fetus.

APPENDIX 2

Guidance and manufacturing declaration - electromagnetic emissions - for all EQUIPMENT AND SYSTEMS

Guidance and Declaration of Manufacture - Electromagnetic Emissions		
The pocket foetal detector is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in these environments.		
Emission Testing	Compliance	Electromagnetic environment-guidance
CISPR 11 RF Emissions	Group 1	The pocket fetal detector uses RF energy only for its internal function. Therefore, its RF emissions are very low and do not cause interference with nearby electronic equipment.
CISPR 11 RF Emissions	Class B	The pocket foetal detector is suitable for use in any environment, including domestic and other environments connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guidance and Declaration of Manufacture - electromagnetic immunity - for all EQUIPMENT AND SYSTEMS

Guidance and Declaration of Manufacture - Electromagnetic Immunity			
The pocket foetal detector is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in these environments.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic-guided environment
Electrostatic discharge (ESD) IEC61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be made of wood, cement or ceramic tile. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Electric frequency magnetic field (50 Hz) IEC 61000-4-8	3 A/m	3 A/m	Magnetic power frequency fields should be at levels characteristic of a typical hospital or commercial environment.

Guidance and Declaration of Manufacture - electromagnetic immunity - for NON-LIFE ASSISTANCE EQUIPMENT AND SYSTEMS

Guidance and manufacturing specifications - electromagnetic immunity			
The pocket foetal detector is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in these environments.			
Immunity test	IEC 60601 test level	Compliance	Electromagnetic environment-guidance
IEC 61000-4-3 radiated RF	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the fetal pocket detector, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rate of the transmitter in watts (W) according to the transmitter manufacturer and, d is the recommended separation distance in metres (m). The extent of the field from fixed RF transmitters, as determined in an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol:</p>
NOTE 1: The highest frequency range applies at 80 MHz and 800 MHz.			
NOTE 2: These guidelines may not apply to all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a. The extent of fields from fixed transmitters, such as base stations for radio (cellular/wireless), telephones and land mobile radios, amateur radio, AM and FM radio broadcasting and TV broadcasting cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey must be considered. If the measured field extent at the site where the fetal detector is used exceeds the applicable RF compliance level, the detector should be tested to verify normal operation. If abnormal performance is observed, additional measures such as reorientation or relocation of the detector may be required.			
b. Above the frequency range 150 kHz to 80 MHz the extended fields must be less than 3 V/m.			

Recommended separation distances for portable and mobile RF communications equipment and EQUIPMENT or SYSTEMS - for NON-LIFE ASSISTANCE EQUIPMENT or SYSTEMS

Recommended separation distances for portable and mobile RF communications equipment and the Pocket Fetal Detector		
The pocket foetal detector is intended for use in an electromagnetic environment where radiated RF disturbances are controlled. The customer or user of the pocket fetal detector can help prevent electromagnetic interference by maintaining a minimum distance between the portable and mobile RF communications equipment (transmitters) and the pocket fetal detector as recommended below, according to the maximum output power of the communications equipment		
Energía de salida Rated maximum output power of the transmitter (W)	Separation distance according to transmitter frequency	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz
0.01	0.01167	0.2334
0.1	0.3689	0.7378
1	1.1667	2.3334
10	3.6893	7.3766
100	11.6667	23.3334
For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		
NOTE 1: The separation distance for the higher frequency range applies to 80 MHz and 800 MHz.		
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

APPENDIX 3

Symbols key

	CF type		Recycle
	Read instructions before use		This side up
	This device is compliant with the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.		Fragile, handle with care
	WEEE (2002/96/EC)		Keep dry
	European representative		Transport and storage temperature limit
	Serial number		Transport and storage humidity limit
	Manufacturer		Transport and storage atmospheric pressure limit
	Manufacturing date		Earphone socket

CONTEC MEDICAL SYSTEMS CO., LTD
 No. 112 Qinhuang West Street, Economic & Technical Development Zone, Qinhuangdao, Hebei Province, P.R.C.

SHANGHAI INTERNATIONAL HOLDING CORP. GMBH (EUROPE)
 Eiffestrasse 80, 20537 Hamburg, Germany

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GRUPO R. QUERALTÓ S.A.