MANUAL DE INSTRUCCIONES PULSIOXÍMETRO PX-01

Manual de instrucciones **ES**

Instruction manual **EN**

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Manuale D'Istruzioni **IT**

Anweisungen **DE**

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Bruksanvisning **SWE**

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AVISO

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NOTICE

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

INDICATIONS

- In all types of patients to detect episodes of hypoxaemia that are difficult to detect clinically.
- Situations of respiratory failure of any aetiology.
- Emergency tracheal intubation and connection to mechanical ventilation.
- As a heart rate monitor in any emergency procedure.

This manual is written in accordance with the MDD93/42/EEC directive for medical devices and harmonised standards. The manual contains instructions for the current pulse oximeter. In case of software modifications and updates, the information contained herein is subject to change without notice. Please read these instructions carefully before using this Mobiclinic PX-01 pulse oximeter. These instructions describe the operating procedures that must be followed as meticulously as possible. Failure to follow these instructions may result in abnormal measurements, equipment damage and personal injury. The manufacturer is NOT responsible for safety, performance, reliability or any abnormal monitoring or personal injury or damage to equipment due to user negligence in complying with the operating instructions. The manufacturer's warranty does not cover such damage. This product is a medical device and may be used repeatedly.

1. SAFETY**1.1 Instructions for safe use**

Check the main unit and all accessories periodically to ensure that there is no visible damage that could affect patient safety or monitoring. It is recommended to check the equipment at least once a week for continuous use. Please do not use the equipment if any damage is detected. Maintenance should ONLY be performed by qualified personnel. Users must not perform maintenance without supervision. The oximeter must not be used in conjunction with other equipment not specified in the user manual. Please use equipment recommended by the manufacturer.

1.2 Warnings

- Explosion Hazard - DO NOT use the Oximeter in an environment with flammable gas such as some anaesthetic agents.

- DO NOT use the oximeter while the patient is undergoing a CT or X-ray examination.
- Be careful with the carrying cord and do not break it during use to avoid damage to the equipment. Please do not use the carrying cord if you are allergic to its fabric.
- Please do not use the product if you are allergic to the rubber of the pulse oximeter cushions.
- Please dispose of the equipment, accessories and packaging (including the plastic bag, foam and cardboard) in accordance with applicable law.

1.3 Precautions

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperatures and humidity.
- If the oximeter gets wet, please stop using it.
- When moving from a cold to a warm and humid environment, please do not use the equipment immediately.
- Do not press the keys on the front panel with sharp materials.
- High temperature and high pressure steam are not allowed to disinfect the oximeter. Refer to the User Manual for cleaning and disinfection.
- Do not immerse the equipment in any liquid. When cleaning is required, wipe with a cloth lightly moistened with a mild disinfectant solution. Do not spray any type of liquid directly on the equipment.
- When cleaning the equipment with water, the temperature of the water should be below 60°C.
- Fingers that are too thin or too cold can affect the patient's normal SpO₂ and pulse rate reading, use a thicker finger such as thumb or middle finger and insert it deep enough into the sensor.
- Do not use the device on newborns or infants.
- The product is suitable for adults and children aged 5 years and older (patient weight should be between 15 kg and 110 kg).
- The device may not work for all patients. If you are unable to get a stable reading, discontinue use.
- The upload period is less than 5 seconds, which may vary from patient to patient.
- The waveform is normalised. Please read the measured value when the waveform on the display is stabilised and constant. At that point you will be able to obtain the optimal value.
- If any abnormal condition appears on the display during the test, remove your finger and reinsert it to restore normal use.
- The device has a useful life of three years from first use.

- The pulse oximeter does not have a low voltage alarm function, it only displays it on the screen. Do not use in situations where an alarm is required.
- The batteries should be removed if the device is to be stored for more than one month.
- A flexible circuit connects the two parts of the device. Do not bend or pull the connections.

2. FEATURES

The oxygen saturation pulse is the percentage of HbO₂ in total blood Hb, also called blood O₂ concentration. It is an important bioparameter for respiration. This pulse oximeter has been developed for the purpose of measuring SpO₂ more easily and accurately. It can simultaneously measure heart rate. The finger pulse oximeter is small and with low power consumption, very convenient to carry and use. Simply place your finger on the sensor of the pulse oximeter and the SpO₂ value will appear on the display immediately.

2.1 Classification

Class II B, (MDD93/42/EEC Annex IX Regulation 10)

2.2 Features

- The operation of the product is simple and convenient.
- The device is small, light (total weight is 50g with batteries) and easy to carry.
- Two AAA batteries are required for operation. The batteries allow for approximately 30 hours of operation.
- The device will automatically switch off if no signal is detected for 5 seconds.
- The orientation of the display can be changed automatically and is easy to view.

2.3 Environment of use

The pulse oximeter can detect SpO₂ and rhythm through the patient's finger and indicates the pulse intensity through the plethimographic waveform. This device is suitable for use in the home, hospitals (primary care), family practice, sports medicine (can be used before or after sports, but not recommended during sports), and others.

WARNING: The device is not suitable for continuous patient monitoring.

WARNING: Higher values may occur when the patient is suffering from carbon monoxide toxicosis. Use of the device in these circumstances is not recommended.

2.4 Environmental Conditions

Transport and Storage

- Temperature: -40° C ~ 60° C
- Humidity: ≤ 95 %.
- Pressure: 500 hPa ~ 1060 hPa


Operation

- Temperature: 10° C ~ 40° C
- Humidity: ≤75 %.
- Pressure: 700 hPa ~ 1060 hPa

3. PRINCIPLE AND PRECAUTIONS

3.1. Measurement principle

A data processing experience formula is established by taking the use of Lambert Beer's Law according to the absorption characteristics of the spectrum of reducing haemoglobin (Hb) and oxyhaemoglobin (HbO₂) in the brightness and near-infrared zones. The working principle of the device is the oxyhaemoglobin photoelectric inspection technology is adopted in disagreement with the capacitance pulse scanning and recording technology, so that two light beams of different wavelength can be focused on the tip of human fingernail through a perspective gripper finger type sensor. Then the measured signal can be obtained by a photosensitive element, the information acquired through which will be displayed on the screen through processing in electronic circuits and microprocessor.

3.2. Precautions before use

- The finger must be properly positioned, otherwise an inaccurate measurement may occur.
- The path between the SpO₂ sensor and the receiving photoelectric tube must be interrupted by a patient artery.
- The SpO₂ sensor must not be placed on an extremity with a strangulated arterial channel or pressure cuff or with an intravenous line for drug administration.
- Ensure that the optical path is free of any obstructions, such as rubber debris, otherwise venous pulsation and inaccurate SpO₂ measurement may result.
- Excessive ambient light may affect the measurement result. This includes fluorescent lights, dual ruby lights, infrared heaters, direct sunlight, etc.
- Excessive patient activity or electrosurgical interference.
- The patient should not wear nail polish or other make-up.

3.3. Clinical constraints

- As the measurement is taken based on arterial pulse, substantial pulsation of the patient's blood flow is necessary. For a patient with a weak pulse due to shock, low body/environment temperature, major bleeding or the use of a vasoconstrictor drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.
- For those with a significant amount of staining drug dilution (such as methylene blue, indigo green and indigo blue acid), or carbon monoxide in haemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic acid in haemoglobin, and some with stroke, SpO₂ determination by this monitor may be inaccurate.
- Drugs such as dopamine, procaine, prilocaine, lidocaine and butacaine can induce a serious error in SpO₂ measurement.
- The SpO₂ value only serves as a reference to assess anaemic anoxia and toxic anoxia, some patients with severe anaemia may give good SpO₂ measurements.

4. TECHNICAL SPECIFICATIONS

Information	Display mode
SpO ₂	LCD
Heart rate	LCD
Pulse intensity	LCD bar graph
Pulse waveform	LCD
Specification of SpO ₂ parameters	
Measurement range	0%~100% (the resolution is 1%)
Accuracy	70%~100% ±2%, unspecified below 70%.
Optical sensor	Red light (wavelength 660 nm) Infrared (wavelength 880 nm)


Specification of pulse parameters	
Measuring range	30 ppm ~250 ppm (resolution is 1 ppm)
Accuracy	± 2 ppm or $\pm 2\%$, whichever is higher
Pulse intensity	
Range	Continuous display of the bar graph, the larger the display, the stronger the pulse.
Battery requirements	
Alkaline batteries x2 x 1.5 V (AAA)	
Battery life	
Batteries can operate continuously for 20 hours	
Size and weight	
Size	61 mm x 36 mm x 32 mm
Weight	Approx. 57 g (with battery)

5. ACCESSORIES

Hanging lanyard / Two AAA batteries / One instruction manual

6. INSTALLATION

6.1 Front panel view

6.2 Batteries

1. Refer to the figure below. Insert two AAA batteries in the correct direction.
2. Replace the cover.

WARNING: Be careful when inserting the batteries as incorrect insertion may damage the device.

6.3 Attaching Hanging Cord

1. Pass the end of the rope through the hole.
2. Pass the other end of the cord through the first end and tighten.

7. OPERATION GUIDE

- a) Insert the batteries according to the polarity indication. Replace the cover.
- b) Open the clamp.
- c) Place the patient's finger inside the device between the two rubber pads (make sure that the finger is well positioned) and release the clamp.
- c) Press the start button on the front panel.
- d) Use your index, middle or ring finger to take the measurement. A low finger temperature or a finger that is too thick or too thin may cause an incorrect measurement.
- e) Do not shake the finger and keep the patient at rest during the measurement.
- f) Observe the information directly from the display.
- g) During the measurement, the device may change the viewing direction depending on the direction of the hand. There are four different orientations.
- h) During the measurement, a long press of the start button changes the brightness of the display.

WARNING: With the finger placed for measurement inside the instrument, make sure that the fingernail is facing upwards.

8. REPAIR AND MAINTENANCE

1. Replace the batteries when the low battery warning lights up on the display.
2. Clean the surface of the equipment after each use. Moisten the equipment first with alcohol and allow it to air dry or wipe it dry with a clean, dry cloth.
3. Remove the batteries from the equipment if you are not going to use it for a long period of time.
4. The storage conditions of the equipment are: ambient temperature between -20°C and 55°C and relative humidity not exceeding 95%.
5. Maintain the equipment in a proper way to ensure correct operation.
6. It is recommended that the device be calibrated periodically.

WARNING: The equipment cannot be sterilised by high pressure steam.

WARNING: Do not immerse the equipment in any liquid.

WARNING: It is recommended that the equipment be kept in a dry environment. Moisture may reduce the life of the device or even damage it.

9. PROBLEM SOLVING

Error	Possible reason	Solution
SpO ₂ and rate are not displayed on the screen.	<ol style="list-style-type: none"> 1. The finger is not positioned correctly. 2. Patient's SpO₂ is too low to be detected. 	<ol style="list-style-type: none"> 1. Position the finger correctly and try again. 2. Try again and go to a hospital if you are sure that the equipment is working.
SpO ₂ and rate are displayed unstable.	<ol style="list-style-type: none"> 1. The finger is not placed far enough in. 2. The patient moves or shakes the finger. 	<ol style="list-style-type: none"> 1. Position the finger correctly and try again. 2. Keep the patient still
The device does not turn on	<ol style="list-style-type: none"> 1. The batteries have run out or are about to run out. 2. Batteries are not inserted correctly. 3. The equipment is not working properly 	<ol style="list-style-type: none"> 1. Replace the batteries. 2. Contact technical support.
Display turns off suddenly	<ol style="list-style-type: none"> 1. The device switches off automatically if no signal is detected for 5 seconds. 2. The batteries are almost empty. 	<ol style="list-style-type: none"> 1. Normal. 2. Replace batteries.

10. SYMBOL KEY

	Type BF
	Please read instructions
%SpO ₂	Peripheral Oxygen Saturation
PRbpm	Pulse rate (beats per minute)
	Low battery (replace batteries to avoid inaccurate readings)
	No finger inserted or insufficient signal
	Battery positive terminal
	Battery negative terminal
	Exit standby mode / Change display brightness
	WEEE (2002/96/EC)
IP22	Ingress Protection
	The device complies with the medical device directive 93/42/EEC of the EEC.
	Manufacturer
	Date of manufacture

	Importer
	Transport and storage temperature limit
	Transport and storage humidity limit
	Atmospheric pressure transport and storage limit
	Upward
	Fragile
	Keep dry
	Recyclable

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9. PROBLEM SOLVING

Error	Possible reasons	Solution
SpO ₂ and rate are not displayed on the screen.	<ol style="list-style-type: none"> 1. The finger is not positioned correctly. 2. Patient's SpO₂ is too low to be detected. 	<ol style="list-style-type: none"> 1. Position the finger correctly and try again. 2. Try again and go to a hospital if you are sure that the equipment is working.
SpO ₂ and rate are displayed unstable.	<ol style="list-style-type: none"> 1. The finger is not placed far enough in. 2. The patient moves or shakes the finger. 	<ol style="list-style-type: none"> 1. Position the finger correctly and try again. 2. Keep the patient still
The device does not turn on	<ol style="list-style-type: none"> 1. The batteries have run out or are about to run out. 2. Batteries are not inserted correctly. 3. The equipment is not working properly 	<ol style="list-style-type: none"> 1. Replace the batteries. 2. Contact technical support.
Display turns off suddenly	<ol style="list-style-type: none"> 1. The device switches off automatically if no signal is detected for 5 seconds. 2. The batteries are almost empty. 	<ol style="list-style-type: none"> 1. Normal. 2. Replace batteries.

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CONTEC MEDICAL SYSTEMS CO., LTD
GRUPO R. QUERALTO S.A.