# MANUAL DE INSTRUCCIONES PULSIOXÍMETRO PX-02

Manual de instrucciones **ES**

Instruction manual **EN**

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Manuale D'Istruzioni **IT**

Anweisungen **DE**

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Bruksanvisning **SWE**

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**AVISO**

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**NOTICE**

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**INDICATIONS**

- In all types of patients to detect episodes of hypoxaemia that are difficult to detect clinically.
- Situations of respiratory failure of any aetiology.
- Emergency tracheal intubation and connection to mechanical ventilation.
- As a heart rate monitor in any emergency procedure.

This manual is written in accordance with the MDD93/42/EEC directive for medical devices and harmonised standards.

The manual contains instructions for the current pulse oximeter. In case of software modifications and updates, the information contained herein is subject to change without notice. Please read these instructions carefully before using this Mobiclinic PX-02 pulse oximeter. These instructions describe the operating procedures that must be followed as meticulously as possible. Failure to follow these instructions may result in abnormal measurements, equipment damage and personal injury. The manufacturer is NOT responsible for safety, performance, reliability or any abnormal monitoring or personal injury or damage to equipment due to user negligence in complying with the operating instructions. The manufacturer's warranty does not cover such damage. This product is a medical device and may be used repeatedly.

**WARNINGS**

- Explosion Hazard - DO NOT use the Oximeter in an environment containing flammable gas such as some anesthetic agents.
- DO NOT USE the oximeter while the patient is undergoing a CT or X-ray test.
- Do not take the information displayed on the device as the sole basis for clinical diagnosis. This device is used as an auxiliary diagnostic and should only be used in conjunction with a physician's advice, clinical manifestations and symptoms.
- Users should not repair the device themselves.
- Discomfort or pain sensation may occur if the device is used uninterruptedly, especially in the microcirculation areas of patients. It is recommended that the sensor not be applied to the same finger for more than 2 hours.

- For certain patients, careful selection of the application site should be made. The equipment cannot be placed on oedema or sensitive tissue.
- The light (infrared is invisible) emitted by the device is harmful to the eyes. The user or maintenance personnel must avoid staring into the light.
- The device contains components with silicone, PVC, TPU, TPE and ABS whose biocompatibility has been tested according to and complies with the requirements of ISO 10993-1. The device is not recommended for use by persons who are allergic to any of the components.
- Be careful with the carrying cord and do not break it during use to avoid damage to the equipment. Please do not use the carrying rope if you are allergic to its fabric.
- Please dispose of the equipment, accessories and packaging (including plastic bag, foam and cardboard) according to the applicable law. Packaging materials should be kept away from children.
- The device must not be used with accessories not approved or recommended by the manufacturer or damage to the operator, patient or pulse oximeter may result.
- Inspect the device prior to use for visible damage that could impair user safety or equipment performance. If damage is observed, replace the part before use.
- There are no functional tests to analyse the accuracy of the pulse oximeter.
- Some bump tests or patient simulators may be used to verify that the device is functioning properly.
- When using the device, keep it away from equipment that can generate strong electrical or magnetic fields.
- Measurement accuracy may be affected by interference from electrosurgical equipment.
- If several devices are used on the patient simultaneously, there may be a danger of leakage current.
- Higher values may occur when the patient suffers from carbon monoxide toxicosis. Use of the device in these circumstances is not recommended.
- The device is not intended for treatment.
- Avoid servicing the device when it is in use.

## 1. OVERVIEW

Oxygen saturation pulse is the percentage of HbO<sub>2</sub> in total Hb in blood, also called blood O<sub>2</sub> concentration. It is an important bioparameter for respiration. This pulse oximeter has been developed for the purpose of measuring SpO<sub>2</sub> more easily and accurately. It can simultaneously measure heart rate. The finger pulse oximeter is small and with low power consumption, very

convenient to carry and use. Simply place your finger on the sensor of the pulse oximeter and the SpO<sub>2</sub> value will appear on the display immediately.

### 1.1 Features

Easy to use; small, light and portable; low power consumption.

### 1.2 Environment of use

The finger pulse oximeter can detect SpO<sub>2</sub> and rhythm through the patient's finger and indicates the pulse intensity through the plethimographic waveform. This equipment is suitable for use in the home, hospitals (primary care), family practice, sports medicine (can be used before or after sports, but not recommended during sports), and others.

### 1.3 Environmental conditions

Transport and storage

- Temperature: -40° C ~ 60° C
- Humidity: ≤ 95 %.
- Pressure: 500 hPa ~ 1060 hPa

Operation

- Temperature: 10° C ~ 40° C
- Humidity: ≤ 75 %.
- Pressure: 700 hPa ~ 1060 hPa

### 1.4. Precautions

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and humidity.
- If the oximeter gets wet, please stop using it.
- When moving from a cold to a warm and humid environment, please do not use the equipment immediately.
- Do not press the keys on the front panel with sharp materials.
- High temperature and high pressure steam are not allowed to disinfect the oximeter. Refer to the User Manual for cleaning and disinfection.
- Do not immerse the equipment in any liquid. When cleaning is required, wipe with a cloth lightly moistened with a mild disinfectant solution. Do not spray any type of liquid directly on the equipment.
- When cleaning the equipment with water, the temperature of the water should be below 60°C.
- Fingers that are too thin or cold can affect the patient's normal SpO<sub>2</sub> and pulse rate reading, use a thicker finger such as thumb or middle finger and insert it deep enough into the sensor.
- Do not use the device on newborns or infants.
- The product is intended for use on adults and children 5 years and older (patient weight should be between 15 kg and 110 kg).
- The device may not work for all patients. If you are unable to get a stable reading, discontinue use.
- The upload period is less than 5 seconds, which may vary from patient to patient.

- The waveform is normalised. Please read the measured value when the waveform on the display is stabilised and constant. At that point you will be able to obtain the optimal value.
- If any abnormal condition appears on the display during the test, remove your finger and reinsert it to restore normal use.
- The device has a useful life of three years from first use.
- The pulse oximeter does not have a low voltage alarm function, it only shows it on the display. Please replace the battery when it is exhausted.
- The instrument does not have an alarm function. Do not use it in situations where it is required.
- The batteries must be removed if the device is to be stored for more than one month.
- A flexible circuit connects the two parts of the device. Do not bend or pull the connections.


### 1.5. Clinical Restrictions

- As the measurement is taken based on the arterial pulse, substantial pulsation of the patient's blood flow is necessary. For a patient with a weak pulse due to shock, low body/environment temperature, major bleeding, or the use of a vasoconstrictor drug, the SpO<sub>2</sub> waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.
- For those with a significant amount of staining drug dilution (such as methylene blue, indigo green and indigo blue acid), or carbon monoxide in haemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic acid in haemoglobin, and some with stroke, SpO<sub>2</sub> determination by this monitor may be inaccurate.
- Drugs such as dopamine, procaine, prilocaine, lidocaine and butacaine can induce a serious error in SpO<sub>2</sub> measurement.
- The SpO<sub>2</sub> value only serves as a reference to assess anaemic anoxia and toxic anoxia, some patients with severe anaemia may give good SpO<sub>2</sub> measurements.

## 2. MEASUREMENT PRINCIPLE


A data processing experience formula is established by taking the use of Lambert Beer's Law according to the absorption characteristics of the spectrum of reducing haemoglobin (Hb) and oxyhaemoglobin (HbO<sub>2</sub>) in the brightness and near-infrared zones. The working principle of the device is The oxyhaemoglobin photoelectric inspection technology is adopted in disagreement with the capacitance pulse scanning and recording technology, so that two light beams of different wavelengths can be focused on the tip of human fingernail through a perspective gripper finger type sensor. Then the

measured signal can be obtained by a photosensitive element, the information acquired through which will be displayed on the screen through processing in electronic circuits and microprocessor.

### 3. INSTALLATION

#### 3.1 Front panel view

#### 3.2. Battery

1. Insert two AAA batteries in the correct direction.
2. Replace the cover.

**WARNING:** Be careful when inserting the batteries as incorrect insertion may damage the device.

### 3.3 Attaching the hanging cord

1. Pass the end of the cord through the hole.
2. Pass the other end of the cord through the first end and tighten it.

### 3.4 Accessories

Hanging cord / Two AAA batteries / One instruction manual

## 4 . OPERATION GUIDE

- a) Insert the batteries according to the polarity indication. Replace the cover.
- b) Open the clamp.
- c) Place the patient's finger inside the unit between the two rubber pads (make sure that the finger is well positioned) and release the clamp.
- d) Press the start button on the front panel.
- e) Do not shake the finger and keep the patient at rest during the measurement.
- f) Observe the information directly from the display.

**WARNING: With the finger placed for measurement inside the device, make sure that the fingernail is facing upwards.**

## 5. MAINTENANCE

### 5.1 Cleaning and disinfection

The device must be switched off and the batteries removed before cleaning and must not be immersed in any liquid.

Use 75% isopropanol to clean the surface and allow it to air dry or use a soft cloth. Do not spray any liquid directly on the device and prevent liquid from entering the device.

### 5.2 Maintenance

- A. Periodically inspect the device for damage. Weekly inspection is recommended.
- B. Clean and disinfect before use.
- C. Replace batteries when the low battery indicator appears.
- D. Remove the batteries if the device is not to be used for an extended period of time.
- E. It is not necessary to calibrate the device.

## 6. TROUBLESHOOTING

Error	Possible reason	Solution
SpO2 and rate are not displayed on the screen.	<ol style="list-style-type: none"> <li>1. The finger is not positioned correctly.</li> <li>2. The finger is trembling or the patient is moving.</li> <li>3. The device is being used in the wrong environment.</li> <li>4. The patient's SpO2 is too low to be detected.</li> </ol>	<ol style="list-style-type: none"> <li>1. Position the finger correctly and try again.</li> <li>2. Instruct the patient to hold still.</li> <li>3. Use the pulse oximeter in the indicated setting.</li> <li>4. Try again and go to a hospital if you are sure that the equipment is working.</li> </ol>
SpO2 and rate are displayed unstable.	<ol style="list-style-type: none"> <li>1. Finger is not placed all the way down.</li> <li>2. The patient moves or shakes the finger.</li> </ol>	<ol style="list-style-type: none"> <li>1. Position the finger correctly and try again.</li> <li>2. Keep the patient still</li> </ol>
The device does not turn on	<ol style="list-style-type: none"> <li>1. Batteries are dead or about to be dead.</li> <li>2. Batteries are not inserted correctly.</li> <li>3. The equipment is not working properly</li> </ol>	<ol style="list-style-type: none"> <li>1. Replace the batteries.</li> <li>2. Contact technical support.</li> </ol>
Display turns off suddenly	<ol style="list-style-type: none"> <li>1. The device switches off automatically if no signal is detected for 5 seconds.</li> <li>2. The batteries are almost empty.</li> </ol>	<ol style="list-style-type: none"> <li>1. Normal.</li> <li>2. Replace batteries.</li> </ol>

## 7. SYMBOLS KEY

	Type BF
	Read instructions
%SpO2	Peripheral oxygen saturation

PRbpm	Heart rate (beats per minute)
	Low Battery (replace batteries to avoid inaccurate readings)
	No finger inserted or inadequate signal
	Positive battery terminal/Battery positive terminal
	Negative battery terminal/Battery negative terminal
	Exit standby mode / Change display brightness
	WEEE (2002/96/EC)
IP22	Degree of protection IP
	The device complies with the EEC Medical Device Directive 93/42/EEC.
	Manufacturer
	Date of manufacture
	Importer
	Transport and storage temperature limit
	Transport and storage humidity limit

	Atmospheric pressure limit for transport and storage
	Upwards
	Fragile
	Keep dry
	Recyclable

## 8. TECHNICAL SPECIFICATIONS

SpO <sub>2</sub>	
Display range	0%~99%
Measuring range	0%~100%
Accuracy	70%~100%;±2%;/ 0%~69%; not specified
Resolution	1%
Heart rate	
Display range	30 bpm~250 bpm
Measuring range	30 bpm~250 bpm
Accuracy	±2 bpm or ±2%, whichever is greater

Resolution	1 bpm
Accuracy in low perfusion conditions	Low perfusion 0.4%: SpO <sub>2</sub> : ±4%; PR: ±2 bpm or ±2%, whichever is greater.
Light interference	In normal light conditions and environments the variation of SpO <sub>2</sub> is ≤1%.
Optical sensor	
Red light	Wavelength: approx. 660 nm, output power: <6.65 mW
Infrared light	Wavelength: approx. 905 nm, output power: <6.75 mW
Protection class	Device with internal power supply, applied part type BF
IP degree of protection	IP22
Voltage	DC 2.6 V ~ 3.6 V
Current	≤25 mA
Power supply	1.5V Alkaline batteries x2 (AAA)
Battery life	Two batteries can operate continuously for 24 hours
Size and weight	
Size	61 x 36 x 32 mm
Weight	Approx. 60g (with batteries)

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
## 2. MEASUREMENT PRINCIPLE

A data processing experience formula is established by taking the use of Lambert Beer's Law according to the absorption characteristics of the spectrum of reducing haemoglobin (Hb) and oxyhaemoglobin (HbO<sub>2</sub>) in the brightness and near-infrared zones. The working principle of the device is The oxyhaemoglobin photoelectric inspection technology is adopted in disagreement with the capacitance pulse scanning and recording technology, so that two light beams of different wavelengths can be focused on the tip of human fingernail through a perspective gripper finger type sensor. Then the measured signal can be obtained by a photosensitive element, the information acquired through which will be displayed on the screen through processing in electronic circuits and microprocessor.

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CONTEC MEDICAL SYSTEMS CO., LTD
GRUPO R. QUERALTÓ S.A.