MANUAL DE INSTRUCCIONES PULSIOXÍMETRO PX-03

Manual de instrucciones **ES**

Instruction manual **EN**

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Manuale D'Istruzioni **IT**

Anweisungen **DE**

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Bruksanvisning **SWE**

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AVISO

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NOTICE

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

USER INSTRUCTIONS

Please read these instructions carefully before using the Mobiclinic PX-03 pulse oximeter.

These instructions describe operating procedures that must be followed as closely as possible. Failure to follow these instructions may result in abnormal measurements, equipment damage and personal injury. The manufacturer is NOT responsible for safety, performance, reliability or any abnormal monitoring or personal injury or damage to the equipment due to user negligence in complying with the operating instructions. The manufacturer's warranty does not cover such damage. This product is a medical device and may be used repeatedly.

WARNINGS

- Discomfort or pain may occur if the device is used continuously, especially in the microcirculation areas of patients. It is recommended that the sensor not be applied to the same finger for more than 2 hours.
- For special patients, careful selection of the area of application should be made. The equipment cannot be placed on oedema or soft tissue.
- The light (infrared is invisible) emitted by the device is harmful to the eyes. The user or maintenance personnel must avoid staring into the light.
- The patient may not wear nail polish or other make-up.
- The patient's nails must not be too long.
- Please refer to existing documentation regarding clinical limitations and precautions.
- The equipment is not designed for treatment.

1. SAFETY**1.1 Instructions for safe use**

- Check the main unit and all accessories periodically to ensure that there is no visible damage that could affect patient safety or monitoring. It is recommended to check the equipment at least once a week for continuous use. Please do not use the equipment if any damage is detected.
- Maintenance must ONLY be carried out by qualified personnel. Users must

not perform maintenance tasks without supervision.

- The oximeter may not be used in conjunction with other equipment not specified in the user manual. Please use equipment recommended by the manufacturer.

1.2 Warnings

- Explosion Hazard - DO NOT use the oximeter in an environment with flammable gas such as some anaesthetic agents.

- DO NOT use the oximeter during a CT or scanner test.

- Be careful with the carrying cord and do not break it during use to avoid damage to the equipment. Please do not use the transport cord if you are allergic to its fabric.

- Please do not use the product if you are allergic to the rubber of the pulse oximeter cushions.

- Please dispose of the equipment, accessories and packaging (including plastic bag, foam and cardboard) in accordance with applicable law.

- For paediatric use, always use the device under adult supervision.

1.3 Precautions

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperatures and humidity.

- If the oximeter gets wet, please stop using it.

- When moving from a cold environment to a warm and humid environment, please do not use the equipment immediately.

- Do not press the keys on the front panel with sharp materials.

- High temperature and high pressure steam are not allowed to disinfect the oximeter. Refer to the User Manual for cleaning and disinfection.

- Do not immerse the equipment in any liquid. When cleaning is required, wipe with a cloth lightly moistened with a mild disinfectant solution. Do not spray any type of liquid directly on the equipment.

- When cleaning the equipment with water, the temperature of the water should be below 60°C.

2. INDICATIONS

The oxygen saturation pulse is the percentage of HbO₂ in the total Hb in the blood, also called blood O₂ concentration. It is an important bioparameter for respiration. Many respiratory diseases can cause hypoxia, which can be life-threatening. Thus, SpO₂ monitoring is indispensable in clinical care. The traditional method of measuring SpO₂ is to analyse a patient's blood sample to obtain the partial pressure of oxygen and calculate SpO₂ with a blood gas analyser. This method is complicated and outdated. In order to measure SpO₂

more easily and accurately, the finger pulse oximeter has been developed. The device can simultaneously measure heart rate. The finger pulse oximeter is small and with low power consumption, very convenient to carry and use. Simply place your finger on the sensor of the pulse oximeter and the SpO₂ value will appear on the display immediately. In clinical practice the tolerance is less than $\pm 2\%$ in the range 70 % to 100 %.

2.1 Features

- The device can accurately measure the SpO₂ value and has an audible pulse rate indication and limit alarm, as when the measured data exceeds the upper or lower limits, the alarm will automatically warn.
- The device will automatically switch off if no signal is detected for 5 seconds.
- Specially designed for paediatric use.
- The PX-03 pulse oximeter includes a lithium battery and has a low power consumption operation. When the Low Battery warning is displayed, the device may not operate normally.
- Lightweight and portable, the total weight with batteries is approximately 50 g.
- Classification: Class IIa, MDD 93/42/EEC, Annex IX, Regulation 10.

2.2 Environment of use

- The finger pulse oximeter can detect SpO₂ and rhythm through the patient's finger and indicates the pulse intensity through the bar graph. This device is suitable for use in the home, hospitals (primary care), family practice, sports medicine (can be used before or after sports practice, but not recommended during sports practice), and others.
- The device is not suitable for continuous patient monitoring.

2.3 Environmental conditions

Transport and storage

- Temperature: $-40^{\circ}\text{C} \sim 60^{\circ}\text{C}$
- Humidity: $\leq 95\%$.
- Pressure: 500 hPa \sim 1060 hPa

Operation

- Temperature: $10^{\circ}\text{C} \sim 40^{\circ}\text{C}$
- Humidity: 30 % \sim 75 %.
- Pressure: 700 hPa \sim 1060 hPa

3. PRINCIPLE AND PRECAUTIONS

3.1 Principle of measurement

Pulse oximetry measurement is the measurement using a multifunctional oxyhaemoglobin meter to transmit some narrow bands of the light spectrum

through blood samples and measure the attenuation of the spectrum with different wavelengths following the characteristic that RHB, O2Hb, Met Hb and COHb have to absorb light of different wavelengths. In this way the O2Hb saturation of different fractions is determined. The O2Hb saturation is called "fractional" O2Hb saturation.

Fractional O2Hb saturation = $[O_2Hb / (RHb + O_2Hb + Met\ Hb + COHb)] \times 100$

On the opposite side, the pulse oximeter measures the functional O2Hb saturation:

Functional O2Hb saturation = $[O_2Hb / (RHb + O_2Hb)] \times 100$

This Mobiclinic PX-03 SpO2 oximeter transmits light of only two wavelengths, red light (wavelength 660 nm) and infrared light (wavelength 940 nm), to differentiate HbO2 from HbR. One side of the sensor contains two LEDs and the other side contains a photoelectric detector. The SpO2 oximeter measures the HbO2 saturation in blood with the plethimographic light when the pulse beats. The result is quite accurate when the HbO2 saturation is around 70 % ~ 100 %.


Figure 1. Principle of operation

3.2 Pre-use precautions

- The equipment must be fully checked for normal use before use.
- The finger must be properly positioned (see illustrations in this manual), otherwise inaccurate measurement may occur.
- The path between the SpO2 sensor and the receiving photoelectric tube must be interrupted by a patient artery.
- The SpO2 sensor must not be placed on an extremity with a strangulated arterial channel or pressure cuff or with an intravenous line for drug administration.
- Ensure that the optical path is free of any obstructions, such as rubber debris,

otherwise venous pulsation and inaccurate SpO₂ measurement may result.

- Excessive ambient light may affect the measurement result. This includes fluorescent lights, dual ruby lights, infrared heaters, direct sunlight, etc.
- Excessive patient activity or extreme electrosurgical interference may affect accuracy.
- The patient may not wear nail polish or other make-up.
- Please clean and disinfect the equipment after use according to the user manual.

3.3 Clinical constraints

- Since the measurement is taken based on arterial pulse, substantial pulsation of the patient's blood flow is necessary. For a patient with a weak pulse due to shock, low body/environment temperature, major bleeding, or the use of a vasoconstrictor drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.
- For those with a significant amount of staining drug dilution (such as methylene blue, indigo green and indigo blue acid), or carbon monoxide in haemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic acid in haemoglobin, and some with stroke, SpO₂ determination by this monitor may be inaccurate.
- Drugs such as dopamine, procaine, prilocaine, lidocaine and butacaine can induce a serious error in SpO₂ measurement.
- The SpO₂ value only serves as a reference to assess anaemic anoxia and toxic anoxia, some patients with severe anaemia may give good SpO₂ measurements.

4. TECHNICAL SPECIFICATIONS

A. Display: LCD

Pulse measuring range: 30 ppm ~ 250 ppm

Pulse waveform presentation: bar graph display

B. Power supply:

Rechargeable lithium battery included

Adaptable range: 3,6 V ~ 4,2 V

C. Operating current: ≤ 30 mA

D. Resolution:

SpO₂ 1 %

Heart rate 1 ppm

E. Accuracy:

SpO₂ 70 % to 100 % ±2 digits, below 70 % not specified.

Heart rate ±2 ppm or ±2 % (whichever is greater)

F. Measurements at low perfusion: SpO₂ and pulse values can be displayed

properly when pulse saturation is around 0,4 %.

SpO2 accuracy ± 4 %.

Heart Rate Accuracy ± 2 ppm or ± 2 % (select wider)

G. Ambient light resistance:

The difference between a value measured in natural ambient light conditions and a dark room is less than ± 1 %.

H. Auto power off function:

The equipment shall automatically shut off after 5 seconds in which no finger is detected to be positioned for measurement.

5. ACCESSORIES

- A. A carrying rope
- B. Rechargeable lithium battery with charging accessories
- C. A user's manual in English and Spanish

6. COMMISSIONING

6.1 Front panel view

- As the measurement is taken on the basis of the arterial pulse, a substantial pulsation is necessary.

Figure 2. Front view

6.2 Inserting the battery

1. See figure 3 and insert the included lithium battery, observing the polarity.
 2. Replace the cover.
- Please be careful when inserting the battery, as incorrect insertion may damage the equipment.

Figure 3. Fitting the battery

7. OPERATION GUIDE

Open the clamp and place the patient's finger inside the device between the two rubber pads (make sure that the finger is correctly positioned) and with the fingernail facing upwards as shown in Figure 4.

Figure 4. Put your finger in position

- a) Press the start button on the front panel.
 - b) Use your index, middle or ring finger to take the measurement.
A low finger temperature, a finger that is too thick or too thin may cause an incorrect measurement.
 - c) Do not shake or shake the finger and keep the patient at rest during the measurement.
 - Observe the information directly from the display.
 - During measurement, the orientation of the display can be changed by quickly pressing the start button. There are six different orientations.
 - During the measurement, a long press of the start button changes the brightness of the display.
 - Setting the pulse sound
- After switching on the device, the pulse sound starts. If you press the button the sound is switched off and the pulse sound indicator icon disappears. If you

press the button again, the sound turns on again and the pulse sound indicator icon reappears.

- Alarm setting

The alarm has two functions, the measurement alarm and the low battery alarm: the measurement alarm occurs when the measurement data exceeds limits (SpO₂ below 90 % or the pulse rate is not between 50 bpm and 120 bpm) and the low battery alarm occurs when the battery power is low.

If the alarm function is enabled, when the measurement data exceeds the normal value range, the audible alarm will occur and the measurement value will flash. If you press the button, the alarm is suspended and the alarm sound indicator icon disappears, but the value continues to flash. The alarm function will resume in 30 seconds.

- Charging

Place the rechargeable lithium battery in the charger respecting the polarity. Connect one end of the adapter to the mains socket and connect the other end to the charger with the data line. During charging, a red light will be displayed and when charging is complete, the light will be green.

8. CLEANING AND DISINFECTION

- Use alcohol to disinfect the equipment and allow it to air dry or wipe it with a soft cloth.

- Clean the surface of the device before use.

9. REPAIR AND MAINTENANCE

1. Replace the batteries when the low battery warning lights up on the display.
2. Clean the surface of the device after each use. Dampen the equipment first with alcohol and let it air dry.
3. Remove the batteries from the equipment if you are not going to use it for a long period of time.
4. The storage conditions of the equipment are: ambient temperature between - 20° C and 55° C and relative humidity not exceeding 95 %.

Maintain the equipment in a proper manner to ensure proper operation.

The equipment cannot be sterilised by high pressure steam.

Do not immerse the equipment in any liquid.

It is recommended to keep the equipment in a dry environment. Humidity may reduce the lifetime of the device or even damage it.

10. ERRORS AND PROBLEMS

Error	Possible reason	Solution
SpO2 and rate are not displayed on the screen	<ol style="list-style-type: none"> 1. The finger is not positioned correctly. 2. The finger is trembling or the patient is moving. 3. The device is being used in the wrong environment. 4. The patient's SpO2 is too low to be detected. 	<ol style="list-style-type: none"> 1. Position the finger correctly and try again. 2. Instruct the patient to hold still. 3. Use the pulse oximeter in the indicated environment. 4. Try again and go to a hospital if you are sure that the equipment is working.
SpO2 and frequency displayed unsteadily	<ol style="list-style-type: none"> 1. The finger is not placed all the way down. 2. The patient moves or shakes the finger. 	<ol style="list-style-type: none"> 1. Position the finger correctly and try again. 2. Hold the patient still
The device does not switch on	<ol style="list-style-type: none"> 1. The batteries have run out of power or are about to run out. 2. The batteries are not inserted correctly. 3. The equipment is not working properly 	<ol style="list-style-type: none"> 1. Replace the batteries. 2. Contact technical support.
The screen turns off suddenly	<ol style="list-style-type: none"> 1. The device switches off automatically if no signal is detected for 5 seconds. 2. The batteries are almost empty. 	<ol style="list-style-type: none"> 1. Normal. 2. Replace batteries.

11. SPECIFICATION OF FUNCTIONS

On-screen information	Mode of presentation
Oxygen pulse saturation SpO2	2 digital digits LCD display
Heart rate bpm	3 digit digital LCD display

Pulse intensity (bar graph)	Graphic bar LCD display
SpO2 Parameter specification	
Measuring range	35%~99% (1% resolution)
Accuracy	±2% between 70%~99%. Below 70% not specified.
Optical sensor	Red light (wavelength 660 nm) Infrared (wavelength 940 nm) Photosensor

Rhythm Parameter specification	
Measuring range	30 ppm ~ 250 ppm, (resolution 1 ppm)
Accuracy	±2 ppm or ±2% (select the wider)

Pulse intensity	
Range	Continuous bar graph, the highest display corresponds to the strongest pulse.

Battery Requirement	
Rechargeable lithium battery 3.6 V x 1	

Autonomy

Autonomy of 6 hours
Charges and discharges more than 300 times

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CONTEC MEDICAL SYSTEMS CO.,LTD
en la UE: Shanghai International Holding Corp.
GmbH (Europe)
Imported by/Importado por: Grupo R. Queralto S.A.