# MOBICLINIC

Spanish Leading Company Since 1985

**ES** - Tensiómetro aneroide
**EN** - Aneroid sphygmomanometer
**IT** - Sfigmomanometro aneroide
**FR** - Tensiomètre aneroide
**PT** - Esfigmomanômetro aneroide
**DE** - Aneroides sphygmomanometer
**NL** - Aneroïde sfygmomanometer
**SWE** - Aneroid sphygmomanomet
**PL** - Sfigmomanometr aneroïdowy
**DK**- Aneroid sphygmomanometer

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### ANEROID SPHYGMOMANOMETER TA-01   **EN**

**Notice:**

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

**Indications:**

The aneroid blood pressure monitor is for professional medical or home use. It is a non-invasive blood pressure measurement system designed to measure the systolic and diastolic blood pressure of an adult individual.

**Technical data:**

\*Built to EN ISO 81060-1 standards.
\*Scale grading: 2mmHg.
\* Measurement precision: 0-300mmHg.
\*Measurement precision: ±3mmHg.
\*The lowest value on the scale is "0", the aneroid needle in the "0" position must be within a tolerance of ±3mmHg.

Aneroid sphygmomanometers offer the advantage of being sturdy and easy to transport, thus ensuring top reliability both in professional and domestic use. MOBICLINIC equipment complies with the European regulations in force in conformity with Regulation (UE) 2017/745 for medical devices.

**Pressure measurement:**

- Place the inflatable cuff on the right or left forearm, making sure that its lower border is 1-2cm above the hollow of the elbow.
- Close the cuff so that its marking is positioned on the artery and the diaphragm rests on the inside of the biceps.
- In order to obtain an optimal reading, the diaphragm must be placed on the artery; the slightly flexed arm must rest on a flat surface and the musculature must be totally relaxed.
- Check that the deflating valve on the bulb is tightly closed by rotating it anti-clockwise, insert the stethoscope chestpiece under the cuff, then attach

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the ear-pieces of the phonendoscope to the ears and the start measuring by inflating the cuff to its maximum pressure which should not exceed the expected maximum by more than 20-30mmHg.

- Open the exhaust valve a little by rotating it clockwise and allow the air to escape as uniformly as possible at a rate of 2-3mmHg per second.
- When the pressure inside the cuff drops lower than the arterial peak, the typical tone of the artery will be heard at every beat through the stethoscope.
- When you heard the first one, read on the pressure gauge scale the value that corresponds to the maximum (or systolic) pressure. Listen carefully and when you hear the last tone; read the value on the pressure gauge scale; this is the minimum (or diastolic) pressure.
- After reading the pressure values, let the air out of the cuff, if the reading is doubtful, repeat the operation after allowing 10-15 minutes to elapse.

**Cautions:**

- The equipment must be protected against shocks.
- Never inflate to more than 300mmHg.
- Do not damage the rubber parts with cutting or sharp objects.
- Do not expose the equipment to excessively high or low temperatures.

**Maintenance and care:**

- The appliance should be checked at least every two years, if the find the equipment's error that exceeds the technical data, must be recalibrated.
- To clean your sphygmomanometer use a clean soft cloth.
- Do not immerse the equipment in liquids to clean it.
- Do not hot sterilize.

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