MANUAL DE INSTRUCCIONES

Andador / Rollator VICTORIA

Manual de instrucciones

ES

Instruction manual

EN

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FR

Manuale D'Istruzioni

IT

Anweisungen

DE

Instruções

PT

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NL

Bruksanvisning

SWE

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PL

Brugsanvisning

DK

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NOTICE

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

Please read these instructions carefully to get to know your new rollator before use. The rollator is used by people who need support when walking due to lack of balance, disability or during a post-operative period. These models can support a weight of up to 100 kg.

1. INSTRUCTIONS FOR USE

1. To activate the brakes press the handlebars upwards.
2. Place both feet together in the middle of the rollator frame.
3. Disengage the brakes by lifting the brake handlebars.
4. Grasp both handlebars of the rollator.
5. Move the rollator forward to a distance that is comfortable for you.
6. Move forward into the rollator frame with your weaker leg.
7. Move forward with your more agile leg while leaning on the handlebars. Both feet should be between the rollator wheels. Continue walking repeating steps 5 and 7.

2. TIPS AND PRECAUTIONS

- Always keep a part of your body, such as your toes, inside the rollator frame.
- When standing up from a seated position, push yourself up from the seat rather than pulling on the handlebars of the rollator.
- Take short steps to turn and keep your hips parallel to the front bar of the rollator. Avoid turning your body to the opposite side of the rollator.
- The rollator must be braked if you are seated. Do not use it as a wheelchair. It can only be used to rest temporarily while being braked.
- Do not apply force on one side only or the rollator will tip over.
- Do not go up or down stairs or steps with the rollator.
- Do not adjust the rollator during use. Assembly and disassembly of the rollator should be done only when the rollator is not in use.
- Do not walk on escalators with a rollator.
- Do not change the direction of the rollator suddenly when using the rollator on a slope.
- It must be used on level streets. Do not use on bumpy roads. Avoid hitting the rollator against obstacles.
- Do not walk with a rollator uphill with an incline of more than 7 degrees. Do not walk a rollator downhill with a slope of more than 15 degrees. When walking uphill or downhill, the lateral inclination or slope gradient should not exceed 3.5 degrees.
- Clean the rollator regularly with a clean, dry cloth.

3. INDICATIONS

It is designed to assist walking for the elderly, disabled and people who have difficulty walking. Make sure all parts are tight and all four wheels are in contact with the ground before use.

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GUANGDONG DAYANG MEDICAL TECHNOLOGY CO., LTD
EC REP **SUNGO EUROPE B.V.**
GRUPO R. QUERALTÓ S.A.