# MANUAL DE INSTRUCCIONES

## Silla de ruedas MUSEO

Manual de instrucciones	ES
Instruction manual	EN
Manuel D'Utilisation	FR
Manuale D'Istruzioni	IT
Anweisungen	DE
Instruções	PT
Gebruiksaanwijzing	NL
Bruksanvisning	SWE
Instrukcja obsługi	PL
Brugsanvisning	DK

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**NOTICE**

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

This wheelchair is designed for frequent or occasional use by a single occupant up to 100kg and can be used indoors and outdoors. The wheelchair can be manually moved forwards and backwards. All Mobiclinic wheelchairs are manufactured to CE standards.

**1. FOLDING AND UNFOLDING THE WHEELCHAIR**

To fold the wheelchair, just lift the centre edge of the seat upholstery from underneath. To unfold, tilt it slightly to one side to raise the wheels away from the floor. Then press down on one or both seat rails. Before folding the wheelchair, completely remove the footrests or make sure they are fixed and in an upright position.

**2. LOCK THE WHEELS**

Push the locks forward until they engage in the locked position. Do not attempt to sit down or stand up without having activated the brakes. Always apply both brakes simultaneously.

**3. ADJUSTING THE LENGTH OF THE FOOTREST**


The position of the footrest (standard-flat or elevated) can be adjusted to suit the length of the user's legs. The footrest should be adjusted to support the weight of the user's feet and lower legs in a position to support the weight of the thighs.

If the footrest adjustment is too long, the thighs will be pressed under the front of the seat upholstery. If the footrest adjustment is too short, the user's knees will rise and the buttocks will have to bear excessive weight. Both situations lead to an increased risk of pressure ulcers. The user's knees and hips should be at approximately the same level.

The lower edge of the footrest must be at least 7 cm above the ground to allow sufficient clearance when going over ramps or slopes. To do this, loosen the adjusting screw with a spanner and fold the foot plates inwards or outwards to achieve the desired position. Then tighten the screw carefully.

Make sure the brakes are applied before getting in or out of the wheelchair.

Do not stand on the footrests when getting out of the wheelchair, otherwise the wheelchair may tip over.

Keep the wheelchair away from fire.

#### 4. USE THE REAR SUPPORT PEDALS AND BRAKES.

Use one of the rear support pedals to raise the front wheels, for example, when climbing a kerb. To do this, push down on the rear support pedals with one foot. Do not lift the small front wheels by pushing down on the handles, as this could damage the wheelchair.

#### 5. CLIMBING A KERB

Approach the kerb facing the kerb. The caregiver should use the rear support pedals to raise the front wheels and lower them onto the kerb once the kerb has been raised. Finally, the caregiver should push the wheelchair forward, lifting it slightly to raise the kerb if necessary.

#### 6. LOWERING A KERB

Position the front wheels in front of the kerb line. The caregiver should use the support pedals to raise the front wheels and tilt the occupant slightly backwards. While keeping the wheels up, lower the wheelchair slightly towards the kerb.

#### WARNING

When going down a kerb, the front wheels should be raised to prevent the occupant from falling off. These are recommendations only; they may not be appropriate in all cases.

#### 7. IMPORTANT ASPECTS TO BE TAKEN INTO ACCOUNT

The user or caregiver should check the safety of the wheelchair at regular intervals.

1. Check the handles and non-slip rubbers and make sure that they are solid and secure.
2. For correct adjustment, check the wheel locks and ensure that the large wheels are locked when engaged.
3. Check all bolts, nuts and joints to ensure that they remain sound.
4. If the wheelchair is equipped with tyres, check that the tyre pressure is correct.
5. Check that the adjusted length of the footrest is correct.
6. If other conditions affecting safety are present, discontinue use of the wheelchair.

Remember to activate the wheelchair lock before transferring the patient from the wheelchair. Do not lean forward out of the wheelchair unless both feet are on the floor. This is especially important for users with heavy casts with elevated foot rests.

#### 8. INDICATIONS

- It is designed for the elderly, the patient, the frail and the disabled, as a mobile chair.
- Check the rear wheel axle and front fork before using the wheelchair.
- Make sure they are tightened securely.

#### 9. CONTRAINDICATIONS

- Do not drive on a bumpy road. Avoid hitting the wheelchair against obstacles.

- When driving the wheelchair on inclined terrain, the user must wear the seat belt. When driving downhill with an angle of inclination of more than 5 degrees, the accompanying person should turn the wheelchair face upwards and drive it downhill slowly and carefully.
- Do not drive the wheelchair downhill with the user upside down or the wheelchair will tip over and cause injury to the user.
- Do not drop the wheelchair by its own force of inertia on a step, it could damage the frame or wheels.
- Do not go up or down steps or stairs without the help of an assistant.
- Do not lift the wheelchair with the user seated in it, otherwise you will damage the wheelchair and cause injury to the user.
- Unless specified, do not adjust, mount or dismount the wheelchair with the user seated. Only adjust, assemble or disassemble the wheelchair when the wheelchair is not occupied.
- Do not lean your body out of the chair, to prevent the chair from tipping over.
- Check and ensure that all parts are tight and securely locked before using the wheelchair.
- Do not drive the wheelchair on escalators.
- Check and ensure that all parts are tight and securely locked before using the wheelchair.
- Do not drive the wheelchair on escalators.
- Do not step on the footrest when getting in and out of the wheelchair. Do not stand on the footrest. Otherwise it could tip over.
- Do not tilt the wheelchair or turn the wheelchair suddenly. Do not brake the wheelchair abruptly.
- The parking brake is only for parking the wheelchair. Do not use it to slow down the wheelchair or to brake the wheelchair in motion.
- Do not force the chassis or the push handle to fold the wheelchair. To fold, lift the seat of the wheelchair from the centre.

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**GUANGDONG DAYANG MEDICAL TECHNOLOGY CO., LTD**
**EC REP** **SUNGO EUROPE B.V.**
**GRUPO R. QUERALTÓ S.A.**