MANUAL DE INSTRUCCIONES

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VERONA

Manual de instrucciones

ES

Instruction manual

EN

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FR

Manuale D'Istruzioni

IT

Anweisungen

DE

Instruções

PT

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NL

Bruksanvisning

SWE

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PL

Brugsanvisning

DK

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TO KEEP FOR FUTURE REFERENCE.

NOTICE

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

1. PURPOSE OF THE PRODUCT

The purpose of the power wheelchair is to provide mobility and autonomy to people with physical limitations that prevent them from walking or moving easily on their own. The power wheelchair is designed to facilitate mobility, improve quality of life, reduce physical effort, promote independence, increase safety and comfort.

2. WARNINGS

1. **WARNING!** The power wheelchair must not be used until this manual has been read in full.
2. **WARNING!** The power wheelchair should not be used until the installation and inspection upgrades have been completed.
3. **WARNING!** It is recommended that persons with intellectual or psycho-motor developmental difficulties do not use the electric wheelchair.
4. **WARNING!** The electric wheelchair must not be disassembled or converted, nor must it be replaced with parts not manufactured by our company.
5. **WARNING!** The user shall not sit or stand up from the power wheelchair when the controller is on or the power wheelchair is operating and there is no attendant to perform a protective function.
6. **WARNING!** The electric wheelchair shall not be used when the anti-tip wheel does not deploy or is not in operation.
7. **WARNING!** The power wheelchair must not be tilted or lifted from the side during normal use.
8. **WARNING!** The user shall not stand on the footrest, in order to prevent the power wheelchair from tipping sideways.
9. **WARNING!** Do not change direction when on sloping terrain.
10. **WARNING!** The electric wheelchair cannot be used by two people at the same time, in order to avoid accidents.

11. **WARNING!** When the power wheelchair is in motion, the speed should be reduced before making turns.
12. **WARNING!** When going downhill, drive slowly in the slowest gear; when going uphill, proceed with caution and at low speed.
13. **WARNING!** It is strictly forbidden to transport objects with the wheelchair.
14. **WARNING!** It shall be verified that the wheels are correctly mounted.
15. **WARNING!** Gently pull the joystick controller and avoid moving it backwards and forwards rapidly.
16. **WARNING!** The joystick is the essential component of the wheelchair; the wheelchair should not be parked outdoors for extended periods of time. When it rains, the wheelchair should be placed indoors to avoid dampness.
17. **WARNING!** The user shall ensure that the “manual/electric” lever of the left and right motors are set to “electric” before switching on the controller.
18. **WARNING!** The “manual/electric” handle of the motor shall not be switched to the “manual” position while the wheelchair is in motion.
19. **WARNING!** If the user needs to use radio communication equipment such as mobile phones and laptop computers, the power supply to the wheelchair controller must be disconnected first because it may interfere with them.
20. **WARNING!** Power wheelchairs are suitable for flat, low-slope terrain and should avoid track surfaces with slopes of more than 6 degrees and obstacles of more than 4 cm.
21. **WARNING!** Passage through any culverts and other places with horizontal voids is strictly prohibited.
22. **WARNING!** If the power wheelchair is not to be used for an extended period of time, it must be disconnected from the power switch on the battery box.
23. **WARNING!** Our electric wheelchair must be used outdoors in communal areas.
24. **WARNING!** Driving on main roads is prohibited, as this could lead to traffic accidents. It is also prohibited by road safety regulations.

3. PRODUCT FEATURES

Scope of application

- The electric wheelchair can be used for disabled and elderly people with movement difficulties.

Product composition

- This product consists of the frame, controller, 2 motors (one on each side), a battery, a footrest, armrests, front wheels and rear wheels.


Illustration 3.1

Illustration 3.2

Structural features

- Highly efficient, energy-saving motors.
- Easy to attach folding frame.
- Smart controller: Power button, power display, universal joystick and loudspeaker.
- Safe and reliable electromagnetic braking system.
- Anti-tip device.
- Lithium battery pack (battery in its case) can be quickly disassembled. The battery can be charged independently of the buggy, i.e. outside the buggy.
- Two modes of operation: Autonomous electric driving mode and manual mode.

Technical parameters




- These parameters will change due to occupant weight, operating environment and battery usage.
- The equipment may not be used in environments where flammable anaesthetic gas mixed with oxygen or nitrous oxide is used.

Product type	Device for outdoor spaces
Maximum speed	≤6,0 km/h
Ambient temperature range	-25° , +50°
Weight capacity	≤100kg
Maximum driving distance per load	≥ 20km
Horizontal braking performance	≤1,5m
Braking safely on maximum gradient	≤ 3,6m (6°)
Battery	Lithium battery DC24V × 12 Ah
Obstacle crossing height	≥ 40mm
Maximum width of crossable recess	100 mm
Minimum turning radius	1,2 m
Normal operating conditions:	
Atmospheric pressure range: 86kPa ~ 106kPa	
Relative humidity range: 25% ~ 95%	
Internal power supply	DC24V±5V
Electrical requirements	Type B applied parts, this means that they have a basic level of protection against electric shock and are suitable for applications where there is no prolonged direct contact with the heart or internal organs.
Mode of operation	Continuous operation
Liquid protection input	IPX4
Engine power	130 W in each

4. DESCRIPTION OF SYMBOLS

Symbols relating to power wheelchair safety requirements and their meaning.

	Importer
	European Authorized Representative

	Fragile; must be handled with care.
	Manufacturer.
	Temperature range.
	Atmospheric pressure range.
	Symbol for the marking of electrical and electronic devices according to Directive 2002/96/EC.
	This device complies with (EU) 2017/745 (Medical Device Regulation).
IPX4	Splash-proof: splashes of water in any direction will not cause any harmful effect.
	Note! Review the document accompanying the wheelchair.
	The pollution control mark for electronic information products, which indicates that the period of use for environmental protection is 10 years, except for consumables.
	This sign indicates mandatory contents (must be observed). Specific obligation contents are expressed by words or pictures on or near these elements. The illustration on the left indicates the "general obligation contents".
	This sign indicates prohibited (not allowed) content. Specific prohibited contents are expressed by words or pictures on or near those elements. The illustration on the left indicates "general prohibited content".
	Application Part Type B.
	Keep out of the rain.
	Above.
	Do not roll.

	Stacking layer limit.
	Consult the instruction manual.
	Humidity range.
	Medical device.

5. FOLDING AND UNFOLDING

UNFOLDING

- Remove the cardboard and the chair will look like in picture 5.1: Place the controller in the direction shown in the illustration and lock the mechanical key in the direction shown in the illustration, as shown in illustrations 5.2 and 5.3.

Illustration 5.1

Illustration 5.2


Illustration 5.3

- Place the power wheelchair on a flat surface, hold the handle with your left hand, and press the seat tube with your right hand, with the centre of gravity biased towards the seat tube, as shown in illustration 5.4 until the wheelchair is fully extended, as shown in illustration 5.5.


Illustration 5.4


Illustration 5.5

- Fold down the footrest as shown in illustration 5.6.

- Push the retracril lever inwards by hand. Locking is complete when a “click” sound is heard, as shown in illustrations 5.7 and 5.8.


Illustration 5.6


Illustration 5.7


Illustration 5.8

- Press the knob on the anti-roll wheel and turn the anti-roll assembly anticlockwise until a “click” sound is heard. The anti-roll assembly extends, and the anti-roll chute is about 4 cm above the floor, as shown in illustrations 5.9 and 5.10.


Illustration 5.9


Illustration 5.10

FOLDING

- As shown in illustration 5.6, press the anti-roll assembly retract button and rotate the assembly clockwise until a “click” sound is heard. Once the anti-roll assembly has been pulled back, the user shall not stand on the anti-roll wheel with their foot, as shown in illustration 5.11.


Illustration 5.11

- Hold the wheel of the footrest with your left hand, and pull the retracting lever with your right hand in an upward direction, and then the wheelchair retracts as shown in illustration 5.12. After folding, place the wheelchair upright as shown in illustration 5.13.


Illustration 5.12


Illustration 5.13

6. INTRODUCTION TO THE DESCRIPTION

Description of the controller function


Illustration 6.1



Illustration 6.2


Power button


Horn button


Battery button

Speed acceleration button

Speed indicator light

Speed deceleration button

Joystick

The main function of the joystick is to control the direction and speed of the wheelchair. The direction in which the joystick is pressed is the same as the direction in which the wheelchair moves. The further the joystick is pressed away from the centre position, the faster the wheelchair moves. When the wheelchair stops, the joystick should not be pressed violently, but should be used as gently as possible. The user can control the power wheelchair to move forwards or backwards, and to turn left or right. It will automatically return to its position when the joystick is released, and the wheelchair will brake automatically. This process should be done gently.

Battery indicator

The battery indicator identifies the battery level. The fully charged status can be shown by 1 red LED, 2 yellow LEDs and 2 green LEDs. When the 2 green LED lights are on, it means that the battery is fully charged. When only the yellow or red LED lights are on, it means that the battery is low and charging should be done as soon as possible to ensure normal use of the wheelchair.

Speedometer

Displays the maximum speed setting of the wheelchair. There are five speed settings: gear 1 is the slowest gear and gear 5 is the fastest gear.

Horn button

Press this button to sound the horn.

Acceleration button

Press this button to increase the value of the speed parameter. Each time the user presses it, the gear will be increased by one gear until 5th gear is reached.

Deceleration button

Press this button to reduce the speed parameter value. Each time the user presses it, the gear is downshifted by one gear until 1st gear, which is the lowest gear, is reached.

Conversion between electric and manual driving

When the wheelchair is not to be moved by electric drive and is moved manually, turn the left and right motor handles to make them both in the "manual" position, when driven with electric power, the joysticks should all be in the "electric" position.


Illustration 6.3


Illustration 6.4

Description of the battery box

As shown in Figure 6.5, by pressing the power display button, the indicator light will turn on, which can show the current battery capacity level; once it is released, the light will turn off. When 4 green lights are on at the same time, it means the power supply is sufficient; when 3 green lights are on, it means the power supply is normal; when 2 green lights are on, it means the power supply is insufficient; when 1 green light is on, it should be charged immediately. Disassembly steps: unplug the connector from the controller, turn the knob 180 degrees counterclockwise to remove the battery box.

1. **Unlock:** turn the battery box knob counterclockwise 180 degrees to remove the battery box.


Illustration 6.5

2. Remove the battery box


Illustration 6.6

3. **Install the battery box:** insert the battery box and rotate it clockwise 180 degrees to lock it.


Illustration 6.7

Description of cargo

- It is recommended that users choose to purchase the charger recommended by our company: the lithium battery charger has a nominal output of 24 V and the charger shall comply with the requirements of IEC60601-1.

NOTE

- When charging, turn on the power switch located on the battery box of the power wheelchair, and turn off the power button on the controller.

Note: The power switch located on the battery box is in the off position when shipped from the factory. The user should ensure that the power switch is turned on when charging "I" is the on state and "O" is the off state, as shown in illustration 6.8.

- Insert the charger output connector into the bottom slot of the charger as shown in illustration 6.9.

!Note! When charging from the controller port, it is strictly forbidden to disconnect the controller power cable from the battery box; **DO NOT TOUCH** the pins inside the power connector.


Illustration 6.8


Illustration 6.9

- **Loading method 2:** Remove the battery box from the power wheelchair and insert the charger connector into the slot in the battery box as shown in illustration 6.10.


Illustration 6.10

- The circuit must not be loaded arbitrarily, in order to ensure circuit correctness.
- During the charging process, the battery circuit must not be disconnected to avoid burns or fire.

Steps for normal driving

- Activate the power switch located on the battery box and press "I" to set it to a normal operating state, as shown in illustration 6.8.
- The handles of the two engines are switched, both from the "manual" to the "electric" position.
- **Note: It is forbidden to activate the engine switch handle on slopes.**
- Press the power button on the controller to check if the e-brakes of the electric wheelchair are working. If the wheelchair cannot move, this means that the e-brake is working, otherwise you should contact the manufacturer.
- After sitting in the power wheelchair, the user should press the power button on the controller, the indicator light will turn on and the joystick should be in the middle position.
- The user needs to concentrate on controlling the wheelchair; this is especially important for novice drivers. The joystick can control direction and speed at the same time. Push the joystick slowly in the forward direction, and the electronic brake will release and then start to move. The push range of the joystick is increased to accelerate; otherwise the wheelchair will slow down.
- If the user wants to stop the power wheelchair while driving, they simply release their hand to reset the joystick position and it will come to a smooth stop; while moving forward, if the user suddenly pushes the joystick backwards or presses the power button on the controller, the wheelchair will stop immediately.
- The speed adjustment button on the controller can adjust the speed of the power wheelchair. The user should choose the maximum speed of the power wheelchair according to his/her physical condition and the road conditions on the road.
- The electric wheelchair is suitable for driving on level roads. On muddy, bumpy or uneven roads, the drive station and the control system of the power wheelchair can be damaged.

7. DAY-TO-DAY MAINTENANCE

NOTE: Before servicing, press the power switch on the battery box to the "O" position.

- It is strongly recommended that worn parts be adjusted and replaced. Professional personnel should be sought or the manufacturer should be contacted.

Inspection cycle	Daily	Weekly	Monthly	For each season of the year	Six-monthly
Battery					
Tyre pressure					
Wire					
Frame					
Controller					
Motore					
Connecting piece					
Quilting					
Change of wheels					
Electromagnetic brake					

The method of replacing worn parts is as follows (if it is difficult to replace any accessory, the manufacturer should be contacted in time and asked to replace it).

- **Method of replacing the front wheel:** Unscrew the screw with a spanner, remove the front wheel, install the front wheel, tighten the screw, adjust the screw tension and confirm that the front wheel rotates flexibly.
- **Method of replacing the rear wheel:** Personnel or professional should be sought or the manufacturer should be contacted for repair or replacement.
- **Seat (back) padding:** Use a screwdriver to unscrew the screw; remove the (back) padding from the seat, fit the (back) padding and tighten the screw with a screwdriver.
- **Method of replacing the armrest:** Use a spanner to unscrew the screw; remove the armrest; fit the new armrest and tighten the screw with a spanner.

- **Battery:** Mainly the remaining battery level should be checked. If the battery life has expired, if the battery needs to be replaced, the user can contact the supplier or the manufacturer or purchase locally a battery with the same specifications.

- **Tyre/tyre pressure:** 310 x 50 wheels are inflated to 260 kpa (maximum 325 kpa), which can be adjusted according to different weight changes and temperatures; chronic air leakage occurs in long-term storage or when not in use, which is normal; when using the wheelchair, the following methods of use should be read carefully:

A. When the air volume in the wheel is insufficient, the operating steps are as follows:

1. First press the wheels evenly by hand to make the wheels and rim fit evenly together.
2. Fill the wheels to a suitable level.

B. When the tyres are completely deflated, the operation steps are as follows:

1. Inflate the deflated tyres to a reduced level of about 30% of full capacity, and then press the tyres evenly by hand to make the wheel and rim fit evenly.

2. Fill the wheels until they reach a suitable level.

- **Wires:** Electrical parts and connecting cables, check for damage or breakage; if there is damage or breakage, contact the supplier or find a suitable person to deal with the problem; they should not attempt to carry out repairs themselves.

- **Frame:** Clean the frame covering surface with a soft step and keep it clean; lubricants should not be used to maintain the wheelchair. If you find cracks in the frame, you should contact your supplier.

- **Maintenance of the controller:** Clean the controller and joystick with a cloth dampened with a neutrally diluted detergent. Caution should be exercised when cleaning the joystick. Never use abrasive materials or alcohol-based cleaners for cleaning. When transporting the wheelchair, the controller must be protected from damage.

- **Engine:** Check for lubricant leakage or increased noise level. If this is the case, the supplier or manufacturer should be contacted.

- **Reliability and maintenance of connecting parts:** Always check whether the bolts and nuts on the wheelchair body are tight. If there is a problem, it should be addressed in time to ensure safe driving.

- **Seat cushion:** The seat cover and backrest should be cleaned with warm water and diluted soapy water, and the wheelchair should not be stored in damp places.

- **Electromagnetic brake:** Driving brake. The method of inspection is to allow the wheelchair to travel straight ahead at maximum speed on a level tarmac road and then release the joystick on the controller so that it automatically returns to the original position, and measure the distance travelled from the time of joystick release to stopping. If this distance is greater than the original value, the braking effect has been reduced. If the distance exceeds 1.5 metres, the user should contact the supplier or the manufacturer for repair of the power wheelchair.

- Use and maintenance of batteries:

- Attention should be paid to the charge level indication on the controller panel. If the light is off, the battery should be charged as soon as possible; when the red light is on, the battery is seriously insufficient and should be charged immediately to prevent the battery voltage from becoming too low and affecting the life of the battery.
- The battery is marked with obvious signs for positive and negative, and there are reliable connectors to ensure normal connection of the circuit. Non-professional personnel shall not connect the circuit randomly.
- During the charging process, the temperature of the battery will rise, but it is not allowed to exceed 45°C, charging should be stopped and continued when the temperature has dropped below 35°C. If the wheelchair is parked without being used for an extended period of time, the battery must be recharged at least once a month.

- The battery has a certain service life. After normal long-term use, if there is a normal difference between the mileage of the electric vehicle and the rated mileage, the battery should be replaced.
- The battery shall not be used when the temperature is $\geq 50^{\circ}\text{C}$ or $\leq -20^{\circ}\text{C}$.
- The battery shall be kept clean and dry. Do not hit the battery with hard objects. The battery shall be kept in good condition and out of the reach of children.
- The power switch on the battery box can cut off the power supply to the battery and reduce the natural loss of battery charge. When the wheelchair will not be in use for an extended period of time, the power switch on the battery box should be disconnected.
- "Fully charged state": The habit of keeping the battery fully charged should be developed. Depending on usage, the battery should be recharged in time to keep it in a "full state of charge" for an extended period of time.
- Waste disposal shall be managed in accordance with national regulations on environmental protection.

8. TRANSPORT AND STORAGE

1. Transport : The product shall be kept upright during transport operations, without being rained on and without rolling; care shall be taken when placing the product, and the height of the levels shall be limited.

2. Storage: The product shall be placed in a dry and ventilated place, and shall not be placed in an environment with high temperature and rapid temperature changes; the product shall be isolated from chemically corrosive substances such as acids and alkalis.

3. Environmental restrictions on transport and storage:

Ambient temperature range: $-40^{\circ}\text{C} \sim +65^{\circ}\text{C}$

Relative humidity range: 10% ~ 100%

Atmospheric pressure range: 86kPa ~ 106kPa

9. TROUBLESHOOTING AND SPECIFICATIONS

- When the power button located on the controller is pressed there is no power signal. In this case it should be confirmed that the power switch on the battery box is in the "I" position. If it is in the "O" position, the power switch on the battery box should be pressed first to set it to the "I" position. Since the lithium battery will cut off power when overcharged, the user must first press the power switch to move it to the "I" position to the "O" position, and then press it back to the "I" position.

- When an error occurs, the controller shall emit an alarm sound and flash, and the error can be assessed according to the number of alarm sounds, to act accordingly.

- Through the function of the information query device integrated in the product, the displayed diagnostic number can reflect the nature of the abnormal situation. Those abnormal conditions can still be detected without using other service tools. Sound information means: cyclic reminder alarm sound.

Sound diagnostics	Description	Recommended method of management
1	Bayeux voltaics	The battery voltage is low, and is used after charging; the battery is broken and is used after replacing the battery; or the battery cannot be charged.
2	The right engine error	Check whether the other right, the connecting device and/or the motor cable are loose.
3	The right-hand braking system	The right-hand braking system should be checked if the braking device and/or the cable is loose. And if the brake switch is damaged, or the switch makes poor contact.
4	The left engine error	If the left-hand motor, the connecting device and/or the motor cable is loose, it must be changed.
5	The left-hand braking system	The left-hand braking device shall be modified, if the device and/or the cable is loose. And if the brake switch is damaged, or the switch makes poor contact.
6	The right-hand motor has an overvoltage protection status.	The brake should be checked and the motor drive station should be checked for binding. Current is not high through the ammeter, which may be a controller problem.
7	Joystick	The joystick does not reset to its position or the joystick cable is broken, or the connector is loose.
8	Controller error	The manufacturer should be consulted for repair.

Most power wheelchair breakdowns are related to batteries, motors and controllers.

Number	Failure	Analysis of the cause and solutions
1	Charging indicator does not illuminate after pressing the power button on the controller	<ul style="list-style-type: none"> - The connection between the battery and the controller is incorrect. It must be reconnected correctly after verification - Battery voltage is too low. If it still does not light up when the battery is charged, the battery may have reached the end of its service life. The battery should be replaced - For problems with the controller, the agency, manufacturer or after-sales service should be contacted.
2	High battery voltage	The battery is charged too high. After charging, the voltage is not more than 29,4V.
3	Low battery voltage	<ul style="list-style-type: none"> - The connection between the battery and the controller connector is not reliable, it should be reconnected - The contact resistance of the wiring between the batteries is large; if the contact resistance is not large, the contact surface may oxidise or become loose. The oxide layer should be removed or the connector should be installed.
4	The engine does not work	<ul style="list-style-type: none"> - The connection between the battery and the controller is incorrect. It must be reconnected correctly after verification. - For problems with the engine, you should contact the agency, the manufacturer or the after-sales service.
5	Motor brake foul	<ul style="list-style-type: none"> - Brake connector is unreliable; must be reconnected - Electric brake coil is damaged.

6	No load indication	<ul style="list-style-type: none"> - The connection between the motor and the controller connector is not reliable, it must be reconnected. - Battery life is over or the charger has been damaged; the battery or charger must be replaced.
7	Short driving distance after each charge	<ul style="list-style-type: none"> - The battery is not completely damaged; it must be recharged. - The battery is approaching the end of its service life; the battery should be replaced.

10. DESCRIPTION OF ELECTROMAGNETIC COMPATIBILITY

NOTE:

- This product complies with the requirements of the related content of the EMC (electromagnetic compatibility) standards IEC60601-1-2 and ISO7176-21.
- The user shall assemble and operate the product in accordance with the instructions for use enclosed with the power wheelchair.
- This product may be affected by RF and portable communication equipment. To avoid strong electromagnetic interference, this product should not be used near mobile phones, ovens, microwave ovens, etc...
- The annex with guidelines and manufacturer's declaration should be consulted.

WARNING:

- This product must not be used when adjacent to or stacked with other equipment. If it is necessary to use the product when adjacent to or stacked with other equipment, you must observe and ensure that the product operates normally under the circumstances.
- Electrical or magnetic environmental guidelines in use:

Guidelines and Manufacturer's Declaration - Electromagnetic Emissions	
Emissions test	Regulatory compliance
RF Emissions CISPR 11	Group 1
RF Emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	ND
Voltage fluctuations / intermittent emissions IEC 61000-3-3	ND

Guidelines and Manufacturer's Declaration - Electromagnetic Emissions	
Immunity test	Regulatory compliance
Electrostatic Discharge (ESD) IEC 61000-4-2 ISO 7176-21	± 8 kV contact ± 15 kV Air
Electrical fast transient/bursts IEC 61000-4-4 ISO 7176-21	± 2kV for power supply lines
Pico IEC 61000-4-5 ISO 7176-21	± 1 kV line by line
Voltage drops, short interruptions and voltage variations on incoming power supply lines IEC 61000-4-11 ISO 7176-21	0% UT 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT 1 cycle 70% UT 25/30 cycles a 0° 0% UT 1 cycle 70% UT 25/30 cycles to 0°
Mains frequency (50 Hz) magnetic IEC 61000-4-8 ISO 7176-21	30A/m
Radiated EM RF fields IEC61000-4-3 ISO 7176-21	20 V/m 80MHz-2,7GHz 80% AM a 1kHz
Conducted disturbances induced by RF fields IEC 61000-4-6 ISO 7176-21	3 V/m 0,15MHz-80MHz 6V en ISM and amateur radio bands on 0,15Mhz y 80MHz 80% AM a 1 kHz
Note: UT is the a.c. mains voltage prior to the application of the test level.	

Test specifications for RECEPT Port IMMUNITY to RF wireless communications equipment.

Frequency of tests (MHz)	Banda (MHz)	Service	Service	Maximum power (w)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ±5 kHz deviation 1 kHz sinusoidal	2	0.3	28
- 710 - 745 - 780	704-787	LTE Banda 13,17	Pulse modulation 217 Hz	0.2	0.3	9
-810 -870 -930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Banda 5	Pulse modulation 18 Hz	2	0.3	28
- 1720 - 1845 - 1970	1700 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Banda 1,3,4,25; UMTS	Pulse modulation 217 Hz	2	0.3	28
2450	2400 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Banda 7	Pulse modulation 217 Hz	2	0.3	28
- 5240 - 5500 - 5785	5100 5800	Bluetooth, WLAN, 802.11 b/g/n, R- FID2450, LTE Banda 7	Pulse modulation 217 Hz	2	0.3	28

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT OR SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

11. AFTER-SALES SERVICE

- Notice of guarantee:

If the product has quality problems caused by non-human factors within one week after the date of sale, our company will be responsible for return, replacement or repair. Under the conditions of normal use and storage, if the product presents quality problems due to human factors within one year from the date of purchase, the company will provide maintenance free of charge. If the product has quality problems after one year from the date of purchase. Users may contact our company's after-sales service department, office or agency in accordance with the invoice. Our company supplies parts and components for maintenance at a reasonable charge. If the user cannot provide an invoice, the warranty period will be confirmed by the company's batch number or factory date extended by one month. If foreign users require repairs, they can send it to our company, but the shipping costs will be borne by the customer.

- **Service life** : 3 years (except for worn parts).

- The warranty period for major parts is as follows:

Núm	Party	Guarantee period
1	Frame	3 years
2	Controller	1 year
3	Engine	1 year
4	Lithium battery	1 year

- The following situations are not covered by the warranty:

1. Worn and consumable parts: seat back cushions, wheels, footrests.
2. Failures caused by unauthorised disassembly, repair or modification of the product.
3. Failure caused by accidental drops during use or handling.
4. Improper use or other accidental or human-caused damage.
5. Failure to follow the correct method in the manual in a way that causes breakdowns.
6. Damage caused by unforeseen natural disasters (such as fires, earthquakes, floods, etc).

- The device number is shown on the frame of each power wheelchair.

NOTE: If there is a need for maintenance, a circuit diagram, list of components and information necessary for such maintenance actions can be provided. If you have questions about circuit maintenance, you should contact the manufacturer.

12. ADDITIONAL REMARKS

Table of size and weight parameters of the electric wheelchair.

Model	Total length mm	Total width mm	Total height mm	Seat width mm	Seat height mm	Seat depth mm	Ar-mrest height mm	Armrest distance mm	Back height mm	Distance from the footrest to the ground mm	Folded length mm	Folded width mm	Folding height mm	Front wheel inches	Rear wheel inches	Maximum load kg	Neto weight kg	Description of features	
D130FL	980	580	900	420	460	430	230	450	410	70	470	580	900	8	12	100	21	Fixed hands and flexible feet, fold back style	
<p>Flat seat angle $\geq 2^\circ$ Angle of the backrest $\geq 3^\circ$ Distance between footrest and seat $\geq 300\text{mm}$ The angle between the leg and the supporting surface $\geq 6^\circ$ The front position of the armrest structure $\geq 250\text{mm}$ Horizontal position of the axis 400mm The quality of the dumm testy = supported weight (-2, +5)</p>																			

We reserve the right to change the technology and appearance of this product, subject to change without notice.

Packing list of the electric wheelchair

N°	Number		Quantity	Remarks
1	Electric wheelchair		1	
2	Key	Inner hexagonal mechanism	2	
	Tools	Pipe spanner	1	
3	User manual		1	

If the related parts are damaged, it is recommended to contact the manufacturer and make a return to the original factory for repair and replacement.

NOTE: Yuyue Medical will not be liable for the consequences if customers contravene the operating requirements or purchase individual access accessories.

ES - Peso máximo soportado 100 KG
EN - Maximum weight supported 100 KG
IT - Peso massimo supportato 100 KG
FR - Poids maximum supporté 100 KG
PT - Peso máximo suportado 100 KG
DE - Maximales Gewicht unterstützt 100 KG
NL - Maximaal ondersteund gewicht 100 KG
SWE - Maximal vikt stöds 100 KG
PL - Maksymalna obsługiwana waga 100 KG
DK - Maksimal vægt understøttet 100 KG

Made in P.R.C.

JIANGSU YUYUE MEDICAL EQUIPMENT &
SUPPLY CO. LTD.

No.1 Baisheng Road Development Zone, 212300
Danyang, Jiangsu, China
sales@yuyue.com.cn

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GRUPO R. QUERALTÓ S.A.