# MANUAL DE INSTRUCCIONES

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### CA-01 Plus

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## 1. MASSAGE TABLE ASSEMBLY

1. Release the buckle. Lay the table on its side and unfold.
2. Open up the legs. Leave the lower legs in the original position when you start unfolding.
3. Straighten each brace.

4. Lift the table up until all the legs are on the floor.
5. Open the headrest piece. Insert the bars of the headrest into the holes, put the pillow on the headrest.
6. Loosen the side handle, adjust the height and angle of the headrest and tighten the handle.
7. Position the armrests and check each side of the stretcher before using it.
8. To adjust the backrest, release the buckle and raise the backrest to the desired height.

## 2. STORING THE TABLE

1. Lay the table on its side and fold at an angle. Completely fold the legs under the table.
2. Close the table halves and use the latch to lock the table.

## 3. PRECAUTIONS

- The table should be folded and unfolded with the help of another person.
- Check and tighten all the screws before use.
- Always place the product on a flat and stable surface. Make sure all the legs are at the same level to avoid instability which might cause injury.
- Maintain the table periodically. Make sure all screws are secured and the cables are not deteriorated or frayed. Make sure there are no cracks.
- Do not overload the table beyond the recommended weight limit. This could result in the serious injury to the user. Do not put your fingers into the gap of the table, the buckle or the bottom of the table.
- Do not use alcohol based cleaners when cleaning the PVC cover. This could damage the material. Wipe clean with a damp cloth only.

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**HANGZHOU FIRSTWELL INDUSTRY CO., LTD.**
**GRUPO R. QUERALTÓ S.A.**