USER MANUAL SIMPLE HEAD STETHOSCOPE FD-01

NOTICE

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

The stethoscope is one of the most common used medical instruments, for clinical, hospital and household use. This Mobiclinic stethoscope is a complete set measuring instrument, mainly to listen heart, lungs, other organs and for other purposes.

DIRECTIONS

OPEN THE PRODUCT PACKAGING, THE PRODUCT SHOULD BE SET ACCORDING TO THE STRUCTURE INSTRUCTIONS. INSERT THE EARPLUGS IN THE EARS, AND PUT THE HEAD ON THE ANATOMICAL POINT YOU WANT TO AUSCULTATE.

CAUTION

- Avoid banging the head of the stethoscope to avoid damage to the diaphragm.
- After use, the stethoscope should be packed into a proper packing bag; do not squeeze the diaphragm.
- If diaphragm is damaged or the latex tube aged, they should be replaced; otherwise it will affect the device function.

CLEANING AND DISINFECTION

- Don't use oils and aggressive solvents.
- Whole stethoscope can be wiped with soap solution.
- To disinfect wipe the stethoscope with a suitable disinfecting agent. Take care of the regulations of disinfection in your respective country.
- To clean the earpieces, they can be removed from the stethoscope.
- Earpieces, rings, bells and diaphragms can be immersed similar into liquids for cleaning, but they must be dry for 100% before installing.

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